

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Quinclorac

MRID #: 410635-70

Matrix: Water

Analysis: GC/ECD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

CONFIDENTIAL

CONFIDENTIAL - SECURITY INFORMATION

CONFIDENTIAL - SECURITY INFORMATION

CONFIDENTIAL - SECURITY INFORMATION

CONFIDENTIAL - SECURITY INFORMATION

CONFIDENTIAL - SECURITY INFORMATION

The information contained in this document is classified as CONFIDENTIAL - SECURITY INFORMATION. It is intended for the use of authorized personnel only and should be protected accordingly. This information is not to be disseminated outside the organization without the express written consent of the appropriate authority.

The information contained in this document is classified as CONFIDENTIAL - SECURITY INFORMATION. It is intended for the use of authorized personnel only and should be protected accordingly. This information is not to be disseminated outside the organization without the express written consent of the appropriate authority.

The information contained in this document is classified as CONFIDENTIAL - SECURITY INFORMATION. It is intended for the use of authorized personnel only and should be protected accordingly. This information is not to be disseminated outside the organization without the express written consent of the appropriate authority.

This document is classified as CONFIDENTIAL - SECURITY INFORMATION. It is intended for the use of authorized personnel only and should be protected accordingly. This information is not to be disseminated outside the organization without the express written consent of the appropriate authority.

This document is classified as CONFIDENTIAL - SECURITY INFORMATION. It is intended for the use of authorized personnel only and should be protected accordingly. This information is not to be disseminated outside the organization without the express written consent of the appropriate authority.

BASF CORPORATION CHEMICALS DIVISION
Agricultural Chemicals
100 Cherry Hill Road
Parsippany, New Jersey 07054

410635-70

Study Title

Gas chromatographic determination of 3,7-dichloro-8-quinoline
carboxylic acid (quinclorac) in drinking water. Method 245

Data Requirements

EPA Guideline Number: 164-2

Author

Dr. F. Mayer

Study Completed On

October, 1985


Performing Laboratory

BASF Aktiengesellschaft
Agricultural Research Station
Limbergerhof, W. Germany

Registration Document No. BASF:

85/0488

This report consists of 21 pages.



85/0488 0001

Quinclorac



**MRID
NUMBER**

4106357

PR 86-5 DATA CONFIDENTIALITY CLAIM

No claim of confidentiality is made for any information contained in this study on the basis of its falling within the scope of FIFRA 10 (d) (1) (A), (B), or (C).

Company BASF CORPORATION
Agricultural Chemicals Group

Company Agent: Jack R. Graham 2/28/89
Jack R. Graham Date
Manager, Registration

85/0480

PR 86-5 GOOD LABORATORY PRACTICE CERTIFICATION

This study is not required to meet the requirements for 40 CFR 160,
Good Laboratory Practice Standards.

Company EASF CORPORATION
Agricultural Chemicals Group

Company Agent

Jack R. Graham
Manager, Registration

2/28/87

Date

558

85/6488 0003

BASF AKTIENGESELLSCHAFT
 LIMBUGERGERHOF AGRICULTURAL RESEARCH STATION

APE/RU
 Environmental Research

October 1985
 Method 245

Dr. F. Mayer


Quinclorac	Gas chromatographical determination
Drinking water	

ALL POSITIVE RESULTS OBTAINED WITH THIS METHOD HAVE IN ANY CASE TO BE VERIFIED BY MASS SPECTROMETRY, SINCE IN THIS ULTRA TRACE RANGE CONTRIBUTIONS FROM DIFFERENT INTERFERING SUBSTANCES HAVE TO BE CONSIDERED. SPECIAL CARE SHOULD BE TAKEN TO AVOID CONTAMINATION IN THE ANALYZING LABORATORY ESPECIALLY WHEN THE INVESTIGATED COMPOUNDS ARE OR WERE USED THERE.

1. INTRODUCTION

Chemical name: 3,7-dichloro-8-quinoline carboxylic acid

Structural formula:



Molecular formula: $C_{10}H_6Cl_2NO_2$


Molecular weight: 242.05

85/0488 0004

Solubility: (g substance in 100 g solvent at 20°C)

Solvent:	water	6.2×10^{-3}
	acetone	0.2
	acetonitrile	<0.1
	methylenechloride	<0.1
	ethanol	0.2
	n-hexane	<0.1

Derivatization reaction for the GC determination:



2. DESCRIPTION OF THE METHOD

The active ingredient is absorbed from the water on C₁₈-alkylated silica gel, purified on a precoated silica gel column, methylated with diazomethane and determined by gas chromatography using a capillary column and an electron capture detector (⁶³Ni-ECD).

3. EQUIPMENT

Gilson-minipuls 2 peristaltic pump (Abimed Ana-
lysentechnik GmbH, Ludwigshafener Str. 26,
4000 Düsseldorf, FRG)

PH meter (WTW-Werkstätten, D-8120 Weilheim)

Graduated flasks, 100 ml

Full pipettes, 1, 2, 5 and 10 ml

Measuring cylinders, 500 ml

Flasks, 10 ml NS 10/19 (see Annex 1 for sketch)

Erlenmeyer flasks, 500 ml

Rapid evaporator according to Rettenberger
(N-EVAP; Labotec, Wiesbaden, Gebr. Rettberg,
Rudolf-Diesel-Str. 19, 3400 Göttingen, FRG)

Accessory equipment for column clean-up

Baker-10 extraction system (No. 70 18-0) with
empty reservoir

Adapter (Order No. 7122-0, Baker, Gross-Gerau)

Collective rack for 10 ml flasks

Gas chromatograph:

e.g. Perkin Elmer F22 with ⁶³Ni-ECD

4. REAGENTS

Hexane, dist.

Acetone, dist.

Dichloromethane, dist.

Methanol, dist.

H₂SO₄ 1 M

One-way separating columns, silica gel Si

(Baker, Order No. 70 86-3)

BONDED-PHASE-Octadecyl C 18 (Baker, Order No.
70 31-0)

Empty columns, 3 ml (Baker, Order No. 71 21-3)

Instructions for preparing diazomethane

The following amounts of reagents are required for the preparation of approx. 100 ml ethereal, about 1.5% diazomethane solution in a distillation apparatus with a descending condenser and a dropping funnel: 3 g KOH, dissolved in 5 ml water, diluted with 45 ml methanol, 6 g N-methyl-N-nitroso-P-toluene sulfonamide (diactin) dissolved in 100 ml ether.

3 g KOH dissolved in 50 ml of 90% methanol in the reaction vessel are connected to the apparatus. The receiver vessel cooled with ice and water and filled with 10 ml ether is connected to the descending condenser. 6 g diactin, dissolved in 100 ml ether, are filled into the 250-ml dropping funnel and slowly introduced dropwise into the solution of KOH in water and methanol maintained at 60°C. The whole amount of diactin solution is to be dripped in during about 30 minutes. The mixture of ether and diazomethane is to be distilled off at the same rate in order to avoid an over-concentration of diazomethane in the reaction mixture. Finally, a further 10 ml ether are added via the dropping funnel and distilled off with the remaining diazomethane. Rubber or Latex gloves must always be worn when diactin and diazomethane are being handled.

Standard solutions for additional tests

Quinclorac > 99.5% (Dr. Ohnsorge, BASF AG; APE/RU)

85/0488 0007

Prepare 100 µg/ml in acetone; dilute 1 µg/ml
and 10 ng/ml with water

Standard solutions for gas chromatography

3,7-Dichloro-8-quinoline carboxylic acid methyl
ester > 99.5% (Dr. Ohnsorge, BASF AG, APE/RU)
100 µg/ml; 1 µg/ml; 10 ng/ml in acetone.

3-Ethyl-7-chloro-8-quinoline carboxylic acid
(internal standard) > 99.5% (Dr. Ohnsorge, BASF AG,
APE/RU)
500 µg/ml; 5 µg/ml; 50 ng/ml in acetone.

In order to obtain comparable signal levels,
the concentration of the internal standard must
be five times as great as that of the active
ingredient derivative.

STABILITY OF THE STANDARD SOLUTIONS.

Storage Days	At room temperature in the light	At + 4°C in the dark
	3,7-dichloro-8-quinoline carboxylic acid 200 µg/ml in acetone	
5	100%	100%
71	100%	100%
	3,7-dichloro-8-quinoline carboxylic acid methyl ester 10.2 µg/ml in acetone + hexane = 5 + 95	
1	100%	100%
7	100%	100%
20	100%	100%
272	100%	100%
	3-ethyl-7-chloro-8-quinoline carboxylic acid methyl ester in acetone	
17	96%	96%
36	103%	109%
52	110%	103%

Quinclorac is stable in water for 30 days at pH 5, pH 7 and pH 9 at 70°C.

5. ANALYTICAL PROCEDURE

The expressions in brackets used below refer to the symbols used in the calculation formula in 7.2. All pressure specifications are to be understood as the difference to atmospheric pressure. The method is described for one sample. However, 8 samples can be handled in parallel.

5.1. Solid phase extraction

5.1.1. Preparation of the column

1 g bonded phase octadecyl is introduced into an empty column with a fritted glass disc, the material is covered with the second fritted disc, and it is placed on a Luer fitting on the cover of the extraction system. A reservoir is fixed on the column by means of an adaptor. 10 ml methanol is introduced in each case in order to wet the C₁₈ material and it is sucked through the column at 20 kPa until it has reached the surface of the column packing. The methanol is displaced by 10 ml water with a pH of 2.5 (adjusted with 1 M H₂SO₄).

5.1.2. Enrichment of the active ingredient on the column

500 µ water (= G) are adjusted to a pH of 2.5 with 1 M H₂SO₄ (about 1.5 ml) in a 500-ml Erlenmeyer flask and degasified for a quarter of an hour with helium. The column prepared according to 5.1.1. is connected up to the tube of the pump and placed in the water with its opening

downward. By gently shaking the tube, the air bubble between water and C_{18} material is displaced. Then 2.5 ml water/min. is drawn through the column with the peristaltic pump.

5.1.3. Prewashing of the column

The active ingredient-laden C_{18} column is replaced on the extraction system and each of the unused openings is closed with a plastic stopper. Air is sucked through the column at 40 kPa for half a minute.

4 ml of a mixture of acetone + hexane = 3 + 97 are drawn through the column by means of the reservoir. The cover of the extraction system is removed together with the columns and the stainless steel rack with the 10-ml flasks (receivers) are placed in the basin. Before the cover is put back on, small amounts of washing liquid that may be hanging on the outlets of the Luer fittings are wiped off.

5.1.4. Elution of the active ingredient

7 ml of a mixture of acetone + hexane = 25 + 75 are applied via the reservoir and sucked through the column at 20 kPa. The eluate in the 10-ml flasks is evaporated to dryness in a stream of nitrogen on the N-EVAP at 40°C. The residue in the ultrasonic bath is dissolved in 1 ml dichloromethane.

5.2. Purification of extract on silica gel

The extract from 5.1.4. is now applied to a Baker 3-ml silica gel separating column. The flask is rinsed out with 0.5 ml DCM, and the whole solution and then, for 10 sec., air at 20 kPa are sucked through the column.

Preliminary washing is carried out by drawing 4 ml methanol + DCM = 20 + 80 through the column. Then the stainless steel rack with the 10-ml flasks as receiver is placed in the basin, residues of washing liquid are wiped off the outlets and the cover with the columns is replaced. The active ingredient is eluted with 3 ml methanol + DCM = 40 + 60 (vacuum 20 kPa) and the eluate is evaporated to dryness at 40°C in a stream of nitrogen.

5.3. Methylation

The residue from 5.2. is dissolved in 0.5 ml acetone. 2 ml of ethereal diazomethane solution are added and the whole is left to stand for one hour. The solution must then still be colored yellow; otherwise a further 2 ml diazomethane solution are added and a further hour is waited. Then the solvent is blown off in a stream of nitrogen at 30°C and the residue is dissolved in 0.5 ml of a solution of the internal standard in acetone.

6. GC conditions

Equipment: e.g. Perkin-Elmer, F 22 with ⁶³Ni-
ECD

Column: SE54, 19 m WCOT, 0.28 mm i.d., film
thickness 0.5 µm

Injection volume: 1 µl

Injector: 260°C

Oven: 200°C

Detector: 370°C

Carrier gas: He, 1000 mbar

Make-up gas: Ar + CH₄ = 90 + 10; 30 ml/min.

7. EVALUATION

7.1. Method of evaluation

The signals are evaluated via the peak height.
In the calibration curve the peak height (in
counts or mm) is plotted against the weight of
the injected standard.

A calibration curve is prepared by injecting at
least three standard amounts. For each series
of analyses two samples of water that are free
of the active ingredient (blank samples) and
two additional tests are subjected to the ana-
lytical procedure.

One of the blank samples is taken up in pure
acetone (not in solution with the internal stan-
dard) for the GC injection. If the chromatogram

of this sample is clean at the retention time of the internal standard, the series of samples is evaluated via the internal standard. For this purpose, instead of the peak heights, the quotients of the peak heights of samples (or additional tests and calibration curve) and internal standard are used. The concentration of the internal standard in the final volume must always be the same within a series of analyses!

For the additional tests water that is free of the active ingredient is mixed with a known amount of quinclorac that is in the same order of magnitude as the expected residues. The yield factor (= F) is determined from the additional tests. The evaluation can also be carried out by means of an appropriate computer program.

7.2. Calculation of the content

The content of quinclorac (R) in µg/kg is calculated by means of the following formula:

$$R = \frac{V_E \cdot W_A \cdot F}{G \cdot V_I}$$

G = sample weight in g

V_E = final volume of the extract before injection in ml

V_I = partial volume injected from V_E in µl

F = yield factor determined by investigator through additional tests

$$F = \frac{100\%}{\text{yield in \%}}$$

W_A = amount of active ingredient derived from the calibration curve in pg

85/0488 0014


8. YIELD, REPRODUCIBILITY AND LOWEST CONCENTRATION DETERMINED

Quinlorac in drinking water

Addition of quinlorac ng/kg (ppt)	Evaluated according to standard	Yield %	Mean %	Standard deviation ±	Variation coefficient ± %
200	internal	86.5; 89.0; 82.1; 84.8; 84.3	85.34	2.58	3.02
10 *	external	84.7; 92.2; 72.8; 91.1	85.20	8.90	10.45

* limit of determination

9. ANALYTICAL SCHEME



85/0488 0016

10. REMARKS

All glass apparatus must be rinsed with acetone before being used for analyses at the limit of determination.


In the event that traces of C₁₈ material go into solution in the elution of the active ingredient, they are removed again by the following silica gel column clean-up.


Technical procedure: H. SträGner

Annexes

1. 20-ml flasks for Baker high-grade steel rack
2. Gas chromatograms
3. 4th calculation

10-ml flasks for Baker stainless steel rack





- 1: standard 10 pg
 2: standard 5 pg
 3: standard 2.5 pg
 4: blank sample $V_E = 0.5$ ml, $V_I = 1$ μ l
 5: addition of ng/kg $V_E = 0.5$ ml, $V_I = 1$ μ l
 +: active ingredient derivative
 o: internal standard; on account of the following negative peak
 not suitable here for evaluation

85/0488 0019

BASF Limburgerhof --- Evaluation --- Environmental Chemistry APE/RU

Queue file name	: ST85037.LST	Report number	: 02/01
Investigator	: STRAESSNER	Date of report	: 30/10/85 07:33:47
Sample material	: WATER	Method number	: 245/01/01
Culture	:	Injection volume (µl)	: 1.0
Act. ingred./metab.	: BAS514H	Detection limit (NG/kg):	: 1.0
Internal standard	:	Blank value	: ST231 .01
Derivative determ.	: QUINCLORAC M.ESTER	Derivat. factor	: 0.945
Subst. concerned	: QUINCLORAC	Conversion factor	: 1.0
Eval. accord. to	: Peak height	Mean recovery	: 82.0
Remarks	:		

Protocol	:			
Study	:			
Study section	:	0	0	0
Sample number	:	ST231 .01	ST232 .01	ST233 .01
Descr. of sample	:	U	Z1	Z2
Sample weight (g)	:	500.0	500.0	500.0
Aliquot (%)	:	100.0	100.0	100.0
Final dilution (ml)	:	0.5	0.5	0.5
Peak area subst.	:	2818	12600	14617
Peak height subst.	:	66	402	481
Peak area i.std.	:			
Peak height i.std.	:			
Quotient sub/i.std.	:			
Residue (µg/kg)	:	< detect. limit		
Addition (NG)	:		5.0	5.0
Addition (µg/kg)	:		0.0100	0.0100
Recovery (%)	:		72.8	91.1

85/0488 0020

BASF Limburgerhof --- Evaluation --- Environmental Chemistry APE/RU

Queue file name : ST85037.LST Report number : 03/01
 Investigator : STRAESSNER Date of report : 29/10/85 13:52:04
 Sample material : WATER Method number : 245/01/01
 Culture : Injection volume (µl) : 1.0
 Act. ingred./metab. : BASS14H Detection limit (NG/kg): 1.0
 Internal standard : Blank value : ST231 .01
 Derivative determ. : QUINCLORAC M. ESTER Derivat. factor : 0.945
 Subst. concerned : QUINCLORAC Conversion factor : 1.0
 Eval. accord. to : Peak height Mean recovery : 88.5

Remarks

Values identified with * lie outside the calibration series

Protocol			
Study			
Study section			
Sample number	ST231 .01	ST236 .01	ST237 .01
Descr. of sample	U	Z3	Z4
Sample weight (g)	500.0	500.0	500.0
Aliquot (%)	100.0	100.0	100.0
Final dilution (ml)	0.5	0.5	0.5
Peak area subst.	2818	14997	18312*
Peak height subst.	66	507	546*
Peak area i.std.			
Peak height i.std.			
Quotient sub/i.std.			
Residue (µg/kg)	< detect. limit		
Addition (NG)		5.0	5.0
Addition (µg/kg)		0.0100	0.0100
Recovery (%)		84.7	92.2*

85/0488 0021