

Cover Sheet for

**ENVIRONMENTAL CHEMISTRY METHOD**

***Pesticide Name:*** Clorethoxyfos

***MRID #:*** 412906-19

***Matrix:*** Soil

***Analysis:*** GC/ECD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at [flynt.elizabeth@epa.gov](mailto:flynt.elizabeth@epa.gov).

MRID 41290619

Du Pont Report No. AMR-1194-88

Study Title

(COMMON  
CHLORETHOXY FOS)

RESIDUE METHOD FOR DETERMINATION OF FORTRESS<sup>R</sup> INSECTICIDE  
ACTIVE INGREDIENT DPX-43898 AND ITS OXON ANALOG IN-34158  
IN SOIL BY ELECTRON-CAPTURE GAS CHROMATOGRAPHY

Data Requirement

U. S. EPA Pesticide Assessment Guidelines  
Subdivision O, 171-4

Authors

Walter A. Babicki,  
George F. Barber, and  
Robert V. Slates

Study Completed On

November 10, 1988

Performing Laboratory

E. I. du Pont de Nemours & Co., Inc.  
Agricultural Products Department  
Research and Development Division  
Experimental Station  
P. O. Box 80402  
Wilmington, Delaware 19880-0402

Laboratory Project ID

AMR-1194-88

STATEMENT OF NO DATA CONFIDENTIALITY CLAIMS

No claim of confidentiality is made for any information contained in this study on the basis of its falling within the scope of FIFRA Section 10(d) (1)(A), (B), or (C).

Company: E. I. Du Pont de Nemours and Company, Inc.

Company Agent: David A. Wustner

Date: Nov. 17, 1988

Registration Specialist  
(Title)

David A. Wustner  
(Signature)

We have submitted this material to the United States Environmental Protection Agency specifically under the provisions contained in FIFRA as amended, and thereby consent to use and disclosure of this material by EPA according to FIFRA. Notwithstanding the wording of our marking "TRADE SECRET," this marking by itself conveys no supplemental claims of confidentiality under FIFRA Sections 10(a) or 10(b). In submitting this material to the EPA according to method and format requirements contained in PR Notice 86-5, we do not waive any protection or right involving this material that would have been claimed by the Company if this material had not been submitted to the EPA, nor do we waive any protection or right provided under FIFRA Section 10(g).

GOOD LABORATORY PRACTICE STATEMENT

The Good Laboratory Practice (GLP) requirements specified in 40 CFR Part 160 were not applicable to environmental fate studies at the time this study was completed. However this study was conducted in the spirit of compliance with GLP regulations.

Sponsor: E. I. du Pont de Nemours & Co., Inc.

Submitter: E. I. du Pont de Nemours & Co., Inc.

Study Director: Robert V. Slates Date: 11-10-88  
Robert V. Slates  
Senior Research Chemist

Walter A. Babicki Date: 11-10-88  
Walter A. Babicki  
Research Technician

R.V. Slates for G.F.B. Date: 11-10-88  
George F. Barber  
Senior Research Chemist

Company Representative: David A. Wustner Date: Nov. 17, 1988  
David A. Wustner  
Registration Specialist

TABLE OF CONTENTS

	<u>Page</u>
STATEMENT OF NO DATA CONFIDENTIALITY CLAIMS.....	2
GOOD LABORATORY PRACTICE STATEMENT.....	3
INTRODUCTION AND SUMMARY.....	6
Scope.....	6
SAFETY AND HANDLING CONSIDERATIONS.....	6
MATERIALS AND METHODS.....	7
Equipment.....	7
Reagents and Standards.....	8
Analytical Procedure.....	9
Sample Preparation.....	9
Determination of Sample Dry Weight.....	11
Soil Extraction and Fortification.....	11
Clean-up by Solvent Partitioning.....	12
Instrumentation.....	12
Descriptions and Operating Conditions for Gas Chromatographic Analysis.....	12
Calibration Procedures.....	14
Confirmatory Techniques.....	14
Time Required for Analysis.....	15
Methods of Calculation.....	15
RESULTS AND DISCUSSION.....	17
CONCLUSIONS.....	19
ACKNOWLEDGMENTS.....	20
CERTIFICATION.....	21

TABLE OF CONTENTS (CONTINUED)

TABLES	
I.	Validation of Extraction..... 22
II.	Recovery Data for Determination of DPX-43898 and IN-34158 in 50-g Soil Samples.... 23-25
FIGURES	
1.	Chemical Names and Structural Formulas for DPX-43898 and IN-34158..... 26
2.	Standard Curves for DPX-43898 and IN-34158..... 27
3.	Representative Chromatograms for Determination of DPX-43898 and IN-34158 in Soils..... 28-36
REFERENCES.....	37
APPENDIX	
I.	TIC Chromatogram and Mass Spectra for DPX-43898 and IN-34158..... 38-39

RESIDUE METHOD FOR DETERMINATION OF FORTRESS<sup>R</sup> INSECTICIDE  
ACTIVE INGREDIENT DPX-43898 AND ITS OXON ANALOG IN-34158  
IN SOIL BY ELECTRON-CAPTURE GAS CHROMATOGRAPHY

Walter A. Babicki, George F. Barber, and Robert V. Slates

INTRODUCTION AND SUMMARY

Scope

An analytical method is described for simultaneously determining DPX-43898 (the active ingredient in Fortress<sup>R</sup> insecticide) and its oxon analog IN-34158 in soil. Chemical names and structural formulas for DPX-43898 and IN-34158 are given in Figure 1. The method is based on extraction of the analytes from soil in a hexane/acetone mixture, removal of the acetone from the extract by two water washes, and determination of the analytes in the hexane solution by capillary gas chromatography using an electron-capture detector. The lower limit of quantitation for each analyte in soil is 0.01 ug/g. The recovery efficiencies for DPX-43898 and IN-34158 are 98% (Std. Dev. = 12%, N = 89) and 90% (Std. Dev. = 13%, N = 89), respectively.

SAFETY AND HANDLING CONSIDERATIONS

DPX-43898 and IN-34158 are highly toxic by either oral or dermal exposure. All contact of these compounds with skin, eyes, clothing, or the respiratory system (through inhalation of

the vapors) should be avoided. All work involving DPX-43898 and IN-34158 should be conducted in a hood while wearing eye protection and protective gloves. Based on permeability studies, nitrile rubber gloves should provide acceptable protection for up to 8 hours. Neither butyl nor neoprene rubber gloves provide adequate protection. Wash thoroughly with soap and water after handling these compounds and before eating or smoking. If DPX-43898 or IN-34158 does come into contact with skin, immediately remove contaminated clothing and wash the contaminated skin or hair with plenty of soap and water.

#### MATERIALS AND METHODS

The following sections include a list of suggested equipment, reagents, and recommendations for preparing stock solutions and standards.

##### Equipment

Composite soil samples using a Hobart commercial food chopper.

During soil extraction, tumble samples using a mechanical tumbler driven by a variable-speed motor. The tumbler should hold several bottles or separatory funnels, turning them end-over-end to mix their contents.

After extraction, separate soils and extracts on an IEC Centrifuge Model K (International Equipment Co., Needham Heights,



MA) with 250-mL Nalgene<sup>R</sup> wide-mouth centrifuge bottles (VWR Scientific, San Francisco, CA, Cat. No. 21020-367).

Use Millex<sup>R</sup>-SR 0.5- $\mu$ m disposable filter units (Millipore Corporation, Bedford, MA, Cat. No. SLSR025NB) to remove particulate matter from sample solutions prior to gas chromatographic analysis.

Perform analyses on a Varian Model 3700 gas chromatograph (Walnut Creek, CA) equipped with a Varian Ni-63 electron-capture detector. Perform chromatographic separations on a 25-m x 0.32-mm i.d. Ultra 1 capillary column (Hewlett Packard, Avondale, PA, Cat. No. 19091A-112) with a 0.52- $\mu$ m film thickness of crosslinked methyl silicone gum.

#### Reagents and Standards

The reference standard of DPX-43898 (compound number SD208304, Code 6-1-0-0) was synthesized and assayed at the Shell Agricultural Chemical Company, Biological Sciences Research Center in Modesto, California. Its chemical purity was 99.8%. The reference standard of IN-34158 was synthesized and assayed at the Du Pont Agricultural Products Department, Research and Development Division, Wilmington, Delaware. Its chemical purity was  $92 \pm 5\%$ .

Use HPLC-grade organic solvents (hexane and acetone) from Fisher Scientific Company (Pittsburgh, PA). Use deionized water for solvent partitioning.

Prepare separate stock solutions (100  $\mu\text{g}/\text{mL}$ ) of DPX-43898 and IN-34158 by dissolving 0.100 g of each analyte in 1000 mL of acetone.

For control sample fortifications up to and including 0.40  $\mu\text{g}/\text{g}$ , combine equal amounts of the stock solutions and dilute with hexane. Prepare solutions containing equal concentrations of DPX-43898 and IN-34158 at 0.5, 1.0, 2.0, 5.0, and 10.0  $\mu\text{g}/\text{mL}$ . For fortifications of 1.0  $\mu\text{g}/\text{g}$  and higher, separately pipet the analytes into the samples from the 100  $\mu\text{g}/\text{mL}$  stock standards.

For gas chromatographic analyses, prepare mixed standards containing equal concentrations of each analyte at 0.005, 0.010, 0.020, and 0.03  $\mu\text{g}/\text{mL}$  by serial dilution of the stock 100  $\mu\text{g}/\text{mL}$  standards in hexane. Because the analytes are relatively volatile, do not evaporate the acetone from aliquots of stock standard solutions prior to dilution. When not in use, all solutions of standards and samples should be stored in a refrigerator at 6 °C or colder.

### Analytical Procedure

#### Sample Preparation

To minimize contamination, begin compositing segments from the untreated plots and proceed to those having increasingly higher treatment rates. Composite the deepest core segments

first, because they should contain the lowest concentration of analyte residues."

To composite soil samples, lay all of the frozen soil tubes for one sampling from a single test on a clean surface in a hood. Accurately measure each tube from the soil surface, marking it at depths of 3, 6, and 12 inches. If the tube is longer than 18 inches, continue marking it in six-inch increments. In some instances, the soil may have been compressed during sampling, resulting in soil cores that are shorter than the actual sampling depth. If the deepest segments of the tubes have soil cores less than 5.5 inches long, significant compression has occurred, and the length of the deepest soil segments must be recorded in the study records.

While the soil tubes are still frozen, cut them with a sharp heavy knife, using the marks as guides to provide segments of 0-3, 3-6, 6-12, 12-18, 18-24, 24-30, and 30-36 inches. Begin compositing the soil segments by transferring the soil for the deepest segment of each core to the bowl of a Hobart chopper. Add sufficient dry ice to keep the soil frozen while operating the chopper until the soil mixture is thoroughly homogenized. Divide the homogenous mixture into three approximately equal portions and transfer each to appropriately labeled sample bottles or bags. Close the sample container loosely so dry ice can sublime without introducing moisture. After sublimation is complete, seal containers tightly and freeze the samples until analyzed.

Repeat the compositing procedure for each soil sample depth. Be sure the Hobart chopper bowl and blades are thoroughly washed with soap and water and dried before compositing each segment. Complete the compositing process with minimal delay, being sure to keep the soil frozen to minimize volatilization of analytes.

#### Determination of Sample Dry Weight

Prior to analysis of each composited sample, determine the sample dry weight as follows: Weigh a 50.0-g portion of the sample into a preweighed disposable metal foil pan and allow the sample to air-dry at ambient temperature to a constant weight. Record the final weight as variable "SDW" for use in calculating analyte concentrations.

#### Soil Extraction and Fortification

To determine DPX-43898 and IN-34158 simultaneously in composited soil samples, weigh 50.0 g of undried soil into a 250-mL Nalgene<sup>R</sup> centrifuge bottle. If the sample dry weight is greater than 47.5 g, add 10 mL of deionized water to the sample. To prepare a fortified control to determine recovery efficiency, pipet into the sample 1.0 mL of a 0.50 µg/mL (or a higher concentration standard depending on the fortification level) mixed standard solution containing equal concentrations of DPX-43898 and IN-34158 in hexane. Do not evaporate the solvent because the analytes are volatile. Add 200 mL of a

hexane/acetone mixture (1:1 by volume) to the centrifuge bottle, cap it, and tumble it on a mechanical tumbler for 2 hours at medium speed. Centrifuge the sample at 3000 rpm for 5 minutes, and decant the extract into a 500-mL glass separatory funnel.

#### Clean-up by Solvent Partitioning

To remove acetone from the extract, add 200 mL of deionized water to the extract in the separatory funnel. Shake or tumble the mixture for 2 minutes and allow the phases to separate. Drain off the lower acetone-water phase and discard it. Add another 200 mL of deionized water to the extract and repeat the water wash of the extract. Drain and discard the lower aqueous phase, and transfer the upper hexane phase to a 100-mL volumetric flask. Add enough hexane to the flask to make exactly 100 mL of sample solution. Mix the solution well and filter 3-5 mL of the sample solution through a 0.5- $\mu$ m Millex<sup>R</sup> filter for analysis. Analyze the sample solution by electron-capture gas chromatography as described in the following section.

#### Instrumentation

##### Descriptions and Operating Conditions for Gas

##### Chromatographic Analysis

Simultaneously determine DPX-43898 and IN-34158 in the sample solutions by capillary gas chromatography (GC) using a

Ni-63 electron capture detector. Analyze samples and mixed standards containing 0.005 ug/mL (or higher depending on the expected residue levels) of each analyte, using the following conditions:

Inlet: Capillary splitter system with a short section of packed precolumn for mixing.  
Split ratio: 6  
Temperature 250 °C

Column: Hewlett Packard Ultra 1  
(Crosslinked methyl silicone gum)  
Length: 25 m  
I.D.: 0.32 mm  
Film Thickness: 0.52  $\mu$ m  
Temperature: 165 °C

Detector: Ni-63 electron capture  
Temperature: 340 °C  
Make up gas: Nitrogen

Carrier Gas: Helium  
Inlet Pressure: 30 psi

Injection Volume: 1  $\mu$ L

Retention Times: IN-34158 2.2 min.  
DPX-43898 2.64 min.

### Calibration Procedures

Prepare several mixed standard solutions in hexane at concentrations that span the concentration ranges expected for the analytes in the sample solutions. Analyze the standards frequently during a series of sample analyses to provide data for standard curves. If the analyte peak height for a sample falls outside the working range of the standard curve, quantitatively dilute the sample and reinject it to keep the detector response within the working range of the standard curve.

To determine recovery efficiency, prepare and analyze at least one control sample and one control sample fortified with known quantities of the analytes before sample extraction. Follow this procedure for every 4 or 5 samples analyzed. Fortification levels should approximate the residue concentration of the samples being analyzed.

### Confirmatory Techniques

The identity of the chromatographic peaks for DPX-43898 and IN-34158 were confirmed by gas chromatographic analysis of standards with mass spectrometric detection. Mass spectra for DPX-43898 and IN-34158 are provided in Appendix I.

Time Required for Analysis

A batch of 6-8 samples can be prepared in an eight-hour day and can be analyzed overnight on an automated gas chromatograph.

Methods of Calculation

Determine residues of DPX-43898 and IN-34158 in soil samples by comparing the chromatographic peak height of each analyte with the corresponding peak heights for standards of known concentration. Calculate  $\mu\text{g/g}$  (on a dry-weight basis) for each analyte in soil samples using Equation 1. For fortified control samples, calculate the percent recovery for each analyte using Equation 2. For these calculations, separate standard curves must be prepared for DPX-43898 and IN-34158. The response factor ( $R_f$ ) for each analyte is simply the slope of its standard curve through the origin.

$$(1) \quad \begin{array}{l} \mu\text{g/g of Analyte} \\ \text{in Soils} \end{array} = \frac{[\text{PK}]_s \times \text{VS} \times \text{DF}}{R_f \times \text{VI} \times \text{SDW}}$$

$$(2) \quad \begin{array}{l} \text{Percent Recovery of} \\ \text{Analyte from Soil} \end{array} = \frac{100\% \times [\text{PK}]_R \times \text{VS} \times \text{DF}}{R_f \times \text{VI} \times \text{SP}}$$



Variables for Equations 1 and 2 are defined as follows:

$$R_f = \frac{[PK]_{STD}}{VI \times C} = \text{Instrument response factor for the analyte in } \mu\text{volts}/\mu\text{g.}$$

C = Concentration of DPX-43898 or IN-34158 in the standard solution in  $\mu\text{g/mL}$ .

DF = Dilution factor for samples requiring sample dilution. DF = 1.0 if no dilution is required.

$[PK]_R$  = Chromatographic peak height (in  $\mu\text{volts}$ ) for the analyte in a fortified recovery sample.

$[PK]_s$  = Chromatographic peak height (in  $\mu\text{volts}$ ) for the analyte in an unknown sample.

$[PK]_{STD}$  = Chromatographic peak height (in  $\mu\text{volts}$ ) for analyte in a standard solution.

SDW = Sample dry weight in g.

SP = Weight of analyte (in  $\mu\text{g}$ ) fortified into a sample.

VI = Injection volume (in mL) of the sample or standard solution.

VS = Volume (in mL) of sample solution before injection or dilution.

RESULTS AND DISCUSSION

Standard curves for electron-capture detection of DPX-43898 and IN-34158 are shown in Figure 2. Both standard curves are linear with zero intercept for 1- $\mu$ L injections of the analytes in hexane at 0.005 to 0.03  $\mu$ g/mL.

Representative chromatograms for determination of DPX-43898 and IN-34158 in soils are shown in Figure 3 for a mixed standard solution, an unfortified control, controls fortified before extraction with both analytes at 0.01, 0.20, and 6.0  $\mu$ g/g, and for soil from a corn test plot treated with Fortress<sup>R</sup> insecticide at planting. The soil for these samples was an Illinois clay loam containing approximately 5% organic matter. The analytes were equally well resolved for the other types of soils analyzed, specifically, an Iowa loam and a California loam each containing approximately 3% organic matter. Note that sample solutions for some of the chromatograms of Figure 3 were diluted prior to injection to keep the detector signals within the working range of the standard curves.

To confirm the validity of the soil extraction method, soil sample No. 17 from the aerobic soil metabolism static study (Reference 1) was obtained. This sample had been fortified with <sup>14</sup>C-labeled DPX-43898 and aged in the dark at 25 °C for 61 days. The soil was a Flanagan loam with approximately 5% organic matter. A spatula was used to homogenize the cold sample on a

tray and divide it into 4 equal portions. To determine moisture content, one portion was dried to constant weight. The other three portions of soil were combined and extracted using the procedure described in this method. Total radioactivity in the sample extracts was determined by liquid scintillation counting (LSC). Total radioactivity in the extracted soil was determined by combustion followed by LSC. Table I compares these data with the average data for replicate samples 11 and 29, which were exhaustively extracted during the aerobic soil metabolism static study (Reference 1). Based on these data, the extraction procedure of this analytical method extracted 31% of the total soil radioactivity compared with an average of 30% for the exhaustive extraction method during the metabolism study. This shows that the extraction procedure is valid for quantitatively extracting DPX-43898 from aged soil samples. Data from the aerobic soil metabolism study (Reference 1) report that the soil-bound metabolites comprise the remaining radioactivity.

Studies were conducted to determine the stability of DPX-43898 and IN-34158 at 25 °C in 13 common HPLC-grade (or equivalent grade) organic solvents. The solvents were acetone, acetonitrile, benzene, chloroform, diethyl ether, dimethylformamide, ethyl acetate, hexane, methylene chloride, methanol, 2-propanol, tetrahydrofuran, and toluene. DPX-43898 was stable (< 3% decomposition) for at least 28 days in all solvents tested with one exception; decomposition in N,N-dimethylformamide occurred with a half-life of 124 days.

IN-34158 was stable for at least 28 days in all solvents tested with the exception of methanol and N,N-dimethylformamide, where decomposition occurred with half-lives of 40 and 184 days, respectively.

Recovery efficiencies using this method for determination of DPX-43898 and IN-34158 were determined on soils fortified prior to extraction with equal concentrations of both analytes at 0.01 to 8  $\mu\text{g/g}$ . The recovery data are shown in Table II. The average recovery efficiencies were 98% (Std. Dev. = 12%, N = 89) for DPX-43898 and 90% (Std. Dev. = 13%, N = 890 for IN-34158.

#### CONCLUSIONS

The lower limit of quantitation by this method is 0.01  $\mu\text{g/g}$  for each analyte. This quantitation limit is the lowest concentration that can be reliably quantified in soils. The quantitation limits and recovery efficiencies for DPX-43898 and IN-34158 are adequate for evaluation of residues in soil samples.

ACKNOWLEDGMENTS

The authors thank J. M. Fukuto for providing Sample 17 of the aerobic soil metabolism static study to validate this method's extraction, for combusting Sample 17 after extraction, and for providing exhaustive extraction data on replicate samples 11 and 29.

The authors also thank D. L. Ryan for conducting the study to determine the stability of DPX-43898 and IN-34158 in several organic solvents.

CERTIFICATION

RESIDUE METHOD FOR DETERMINATION OF FORTRESS<sup>R</sup> INSECTICIDE  
ACTIVE INGREDIENT DPX-43898 AND ITS OXON ANALOG IN-34158  
IN SOIL BY ELECTRON-CAPTURE GAS CHROMATOGRAPHY

We, the undersigned, declare that the work described in this report was performed under our supervision, and that this report, to the best of our knowledge, provides an accurate record of the procedures and results.

Report by:

R. V. Slates                      11-10-88  
R. V. Slates                      Date  
Study Director

Walter A. Babicki                      11-10-88  
Walter A. Babicki                      Date

R. V. Slates for G.F.B.                      11-10-88  
George F. Barber                      Date

Approved by:

R. D. Collins                      11/11/88  
Ronald D. Collins  
Environmental Studies  
Supervisor

Date Study Initiated: May 22, 1987. (Protocol AMR-831-87 was signed to start Field Dissipation Study)

Date Report Issued: \_\_\_\_\_

Notebook References: NB-E49630-28 through 152  
NB-E524648-6 through 13

Storage Location of Raw Data and Final Report:  
E. I. du Pont de Nemours & Co., Inc.  
Agricultural Products Department  
Experimental Station and/or  
Du Pont Records Management Center  
P. O. Box 80402  
Wilmington, Delaware 19880-0402

Sponsor:  
E. I. du Pont de Nemours & Co., Inc.  
Agricultural Products Department  
Experimental Station  
P. O. Box 80402  
Wilmington, Delaware 19880-0402

TABLE I

VALIDATION OF EXTRACTION

(Replicate Soil Samples Aged 61 Days at 25°C)

	Activity Found, dpm		Percent of Applied Activity Found Excluding Carbon Dioxide	
	Sample 17 (This Extraction)	Samples 11 and 29 (Exhaustive Extraction)	Sample 17 (This Extraction)	Samples 11 and 29 (Exhaustive Extraction)
Total Extractable	$1.10 \times 10^6$	$1.15 \times 10^6$	31	30
Bound to Soil	$2.46 \times 10^6$	$2.72 \times 10^6$	69	70
Carbon Dioxide Trap	Not Analyzed	$5.18 \times 10^6$		
Material Balance(%)	---	70		

TABLE II

RECOVERY DATA FOR DETERMINATION OF DPX-43898  
AND IN-34158 IN 50-G SOIL SAMPLES

DPX-43898			IN-34158		
PPM Added	PPM Found	Percent Recovery	PPM Added	PPM Found	Percent Recovery
0.010	0.008	80	0.010	0.012	120
0.010	0.010	100	0.010	0.010	100
0.010	0.008	80	0.010	0.008	80
0.010	0.007	70	0.010	0.009	90
0.010	0.012	120	0.010	0.013	130
0.010	0.012	120	0.010	0.013	130
0.010	0.012	120	0.010	0.010	100
0.010	0.011	110	0.010	0.009	90
0.010	0.010	100	0.010	0.007	70
0.010	0.011	110	0.010	0.008	80
0.010	0.011	110	0.010	0.010	100
0.010	0.012	120	0.010	0.010	100
0.010	0.010	100	0.010	0.009	90
0.010	0.010	100	0.010	0.008	80
0.020	0.018	90	0.020	0.017	85
0.020	0.022	110	0.020	0.020	100
0.020	0.019	95	0.020	0.019	95
0.020	0.021	105	0.020	0.021	105
0.020	0.023	115	0.020	0.023	115
0.020	0.019	95	0.020	0.019	95
0.020	0.022	110	0.020	0.020	100
0.020	0.021	105	0.020	0.022	110
0.020	0.021	105	0.020	0.020	100
0.020	0.019	95	0.020	0.016	80
0.020	0.018	90	0.020	0.017	85
0.020	0.023	115	0.020	0.021	105
0.020	0.022	110	0.020	0.019	95
0.040	0.031	78	0.040	0.027	68
0.040	0.041	103	0.040	0.040	100
0.040	0.039	98	0.040	0.038	95
0.040	0.039	98	0.040	0.040	100
0.040	0.043	108	0.040	0.044	110
0.040	0.032	80	0.040	0.034	85
0.040	0.039	98	0.040	0.040	100
0.10	0.093	93	0.10	0.090	90
0.10	0.089	89	0.10	0.105	105
0.10	0.106	106	0.10	0.093	93
0.10	0.083	83	0.10	0.081	81
0.10	0.084	84	0.10	0.071	71
0.10	0.096	96	0.10	0.084	84



TABLE II (Continued)

RECOVERY DATA FOR DETERMINATION OF DPX-43898  
AND IN-34158 IN 50-G SOIL SAMPLES

DPX-43898			IN-34158		
PPM Added	PPM Found	Percent Recovery	PPM Added	PPM Found	Percent Recovery
0.10	0.096	96	0.10	0.079	79
0.10	0.087	87	0.10	0.090	90
0.10	0.103	103	0.10	0.110	110
0.10	0.112	112	0.10	0.102	102
0.10	0.112	112	0.10	0.096	96
0.10	0.119	119	0.10	0.103	103
0.10	0.125	125	0.10	0.109	109
0.10	0.095	95	0.10	0.090	90
0.10	0.086	86	0.10	0.075	75
0.20	0.194	97	0.20	0.179	90
0.20	0.209	105	0.20	0.191	96
0.20	0.183	92	0.20	0.158	79
0.20	0.189	95	0.20	0.159	80
0.20	0.205	103	0.20	0.161	81
0.20	0.207	104	0.20	0.186	93
0.20	0.193	97	0.20	0.170	85
0.20	0.182	91	0.20	0.180	90
0.20	0.197	99	0.20	0.181	91
0.20	0.172	86	0.20	0.154	77
0.20	0.185	93	0.20	0.171	86
0.40	0.366	92	0.40	0.357	89
0.40	0.372	93	0.40	0.365	91
0.40	0.359	90	0.40	0.308	77
0.40	0.335	84	0.40	0.375	94
0.40	0.413	103	0.40	0.360	90
0.40	0.421	105	0.40	0.299	75
1.0	1.15	115	1.0	1.15	115
1.0	0.85	85	1.0	0.75	75
1.0	0.83	83	1.0	0.72	72
1.0	1.02	102	1.0	0.90	90
1.0	0.84	84	1.0	0.69	69
2.0	1.64	82	2.0	1.53	77
2.0	1.81	91	2.0	1.54	77
2.0	2.22	111	2.0	1.86	93
2.0	1.96	98	2.0	1.53	77
2.0	1.98	99	2.0	1.68	84
2.0	1.84	92	2.0	1.61	81
4.0	3.24	81	4.0	2.86	72
4.0	3.43	86	4.0	3.12	78
4.0	3.86	97	4.0	3.41	85
4.0	3.95	99	4.0	2.97	74
4.0	3.21	80	4.0	3.00	75
4.0	3.33	83	4.0	3.25	81


TABLE II (Continued)

RECOVERY DATA FOR DETERMINATION OF DPX-43898  
AND IN-34158 IN 50-G SOIL SAMPLES

DPX-43898			IN-34158		
PPM Added	PPM Found	Percent Recovery	PPM Added	PPM Found	Percent Recover
6.0	6.69	112	6.0	5.80	97
6.0	5.05	84	6.0	4.80	80
6.0	5.49	92	6.0	4.65	78
6.0	6.14	102	6.0	5.64	94
8.0	8.00	100	8.0	6.41	80
8.0	6.84	86	8.0	6.99	87
Average = 98			Average = 90		
Std. Dev. = 12			Std. Dev. = 13		
N = 89			N = 89		


FIGURE 1

CHEMICAL NAMES AND STRUCTURAL FORMULAS  
FOR DPX-43898 AND IN-34158



DPX-43898

Phosphorothioic acid,  
0,0,diethyl-0-(1,2,2,2-tetrachloroethyl)ester



IN-34158

Phosphoric acid, diethyl(1,2,2,2-tetrachloroethyl)ester

FIGURE 2  
STANDARD CURVES FOR DPX-43898 AND IN-34158




FIGURE 3

REPRESENTATIVE CHROMATOGRAMS FOR  
DETERMINATION OF DPX-43898 AND IN-34158 IN SOILS

Analysis Name : 5 E4963087.7.1.  
.02 PPM STD (.01ug). Amount : 50.000.




FIGURE 3. CHROMATOGRAM A. MIXED STANDARD CONTAINING 0.01UG/ML OF DPX-43898 AND 0.01 UG/ML OF IN-34158. (PEAK HEIGHTS ARE APPROXIMATELY EQUIVALENT TO THOSE OF A 50-G SOIL SAMPLE CONTAINING 0.02 PPM OF EACH ANALYTE.)

Analysis Name : 5 E4963080,2,1.  
0-3" CK -1PHI. Amount : 46.500.




FIGURE 3. CHROMATOGRAM B. UNFORTIFIED CONTROL SOIL, SAMPLE DEPTH 0-3 INCHES. BOTH DPX-43898 AND IN-34158 WERE BELOW THE 0.01 PPM QUANTIFICATION LIMIT.

Analysis Name : 5 E4963087, 10.1.  
.01 PPM SPIKE. Amount : 50.000.




FIGURE 3. CHROMATOGRAM C. SOIL FORTIFIED BEFORE ANALYSIS WITH DPX-43898 AND IN-34158 AT 0.01 PPM EACH. RECOVERY WAS 100% FOR DPX-43898 AND 70% FOR IN-34158.



Analysis Name : 5 E4963090,9,1.  
0.2 PPM SPIKE. Amount : 50.000.




FIGURE 3. CHROMATOGRAM D. SOIL FORTIFIED BEFORE ANALYSIS WITH DPX-43898 AND IN-34158 AT 0.2 PPM EACH. DILUTED 1 TO 5 BEFORE INJECTION. RECOVERY WAS 103% FOR DPX-43898 AND 81% FOR IN-34158.

Analysis Name : 5 E4963091.8.1.  
6.0 PPM SPIKE. Amount : 50.000.




FIGURE 3. CHROMATOGRAM E. SOIL FORTIFIED BEFORE ANALYSIS WITH DPX-43898 AND IN-34158 AT 6.0 PPM EACH. DILUTED 1 TO 150 BEFORE INJECTION. RECOVERY WAS 102% FOR DPX-43898 AND 94% FOR IN-34158.

Analysis Name : 5 E4963087.6.1.  
12-18" 0.6oz. 7PHI. Amount : 43.300.




FIGURE 3. CHROMATOGRAM F. SOIL COLLECTED 7 DAYS AFTER TREATMENT OF THE PLOT WITH FORTRESS INSECTICIDE AT 0.6 OZ-AI/KFT. SAMPLE DEPTH 12-18 INCHES. BOTH DPX-43898 AND IN-34158 WERE BELOW THE 0.01 PPM QUANTITATION LIMIT.

Analysis Name : 5 E4963087.2.1.  
0-3" 0.6oz. 7PHI. Amount : 44.500.




FIGURE 3. CHROMATOGRAM G. SOIL COLLECTED 7 DAYS AFTER TREATMENT OF THE PLOT WITH FORTRESS INSECTICIDE AT 0.6 OZ-AI/KFT. SAMPLE DEPTH 0-3 INCHES. THIS CHROMATOGRAM SHOWS THAT IN-34158 WAS LESS THAN 0.01 PPM. DPX-43898 WAS TOO CONCENTRATED TO QUANTIFY. SEE FIGURE 3. CHROMATOGRAM H.

Analysis Name : 5 E4963087.3.1.  
0-3" 0.6oz. 7PHI. Amount : 44.500.





FIGURE 3. CHROMATOGRAM H. SOIL COLLECTED 7 DAYS AFTER TREATMENT OF THE PLOT WITH FORTRESS INSECTICIDE AT 0.6 OZ-AI/KFT. SAMPLE DEPTH 0-3 INCHES. DILUTED 1 TO 25 BEFORE INJECTION. DPX-43898 WAS 0.85 PPM.

REFERENCES


1. J. M. Fukuto, "Study Director, "Aerobic Soil Metabolism of <sup>14</sup>C-DPX-43898," Du Pont Agricultural Products Department, Document No. AMR-873-87. MRID No. 40883706.

APPENDIX I


TIC Chromatogram and Mass Spectra  
for DPX-43898 and IN-34158



TIC CHROMATOGRAM  
OF DPX-43898 AND  
IN-34158 MIXED  
STANDARD.



MASS SPECTRUM  
OF IN-34158.



MASS SPECTRUM  
OF DPX-43898.

TOTAL ION CURRENT CHROMATOGRAM AND MASS SPECTRA FOR DPX-43898 AND IN-34158. TWO MICROLITER INJECTION OF MIXED STANDARD CONTAINING 50 UG/ML EACH OF DPX-43898 AND IN-34158. HEWLETT PACKARD MODEL 5890A CHROMATOGRAPH, MODEL 5970 MASS SELECTIVE DETECTOR, HP-1 COLUMN, 25M LENGTH, 0.33 MICRON FILM THICKNESS, 0.20 MM COLUMN ID, 170 C ISOTHERMAL, COLUMN HELIUM FLOW 0.8 ML/MIN, INJECTION PORT 220 C, TRANSFER LINE 280 C.