

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Prodiamine

MRID #: 413594-05

Matrix: Soil

Analysis: GC/ECD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

1950

STATE OF TEXAS

County of _____

Know all men by these presents that _____ of the County of _____ State of Texas, for and in consideration of the sum of _____ Dollars, to _____ in hand paid by _____ the receipt of which is hereby acknowledged, have granted, sold and conveyed, and by these presents do grant, sell and convey unto the said _____ of the County of _____ State of Texas, all that certain _____

413594-65

APPENDIX IV

ANALYTICAL METHOD AM-0817

25
101

**DETERMINATION OF PRODIAMINE
AND ITS
6-AMINO-IMIDAZOLE METABOLITE
IN SOIL**

METHOD # AM-0817

**SANDOZ CROP PROTECTION CORPORATION
1300 E. TOUHY AVENUE
DES PLAINES, ILLINOIS 60018**

TABLE OF CONTENTS

| | <u>PAGE #</u> |
|--------------------------------------|---------------|
| 1. SUMMARY | 1 |
| 2. SAFETY | 1 |
| 3. MATERIALS / METHODS | |
| 3.1 Apparatus | 2 |
| 3.2 Reagents | 2 |
| 3.3 Standards | 3 |
| 3.4 Procedure | 4 |
| 3.5 Analysis | 6 |
| 3.6 Confirmatory Techniques | 8 |
| 3.7 Time Required for Analysis | 8 |
| 4. RESULTS / DISCUSSION | |
| 4.1 Precision and Accuracy | 9 |
| 4.2 Limits of Detection | 9 |
| 5. CONCLUSIONS | 9 |
| 6. CERTIFICATION | 10 |
| 7. TABLES / FIGURES | 11 |

ANALYTICAL METHOD

Report No. 15, Project 480428

SANDOZ CROP PROTECTION CORPORATION

Location: 1300 E. TOUHY AVE.
DES PLAINES, IL 60018

DEVELOPMENT QUALITY CONTROL

Method Number _____

Addendum _____

Supersedes _____

Approved ISA/S Date 4-2-88

Reviewed _____ Date _____

Reviewed _____ Date _____

Reviewed _____ Date _____

Reviewed _____ Date _____

DETERMINATION OF PRODIAMINE AND ITS 6-AMINO-IMIDAZOLE METABOLITE IN SOIL

1. SUMMARY

- 1.1 This method is applicable to the determination of prodiamine and its 6-amino-imidazole metabolite in soil. This method was developed in the residue laboratory at Sandoz Crop Protection Corporation.
- 1.2 Twenty gram aliquots of sample are extracted by shaking with methanol.
- 1.3 The samples are centrifuged and aliquots taken from the supernatant.
- 1.4 Aliquots (5 g equivalent) of the extracts are diluted with 5% NaCl solution and extracted three times with dichloromethane.
- 1.5 The samples are cleaned up by silica gel column chromatography.
- 1.6 Prodiamine and the 6-amino-imidazole are quantitated by gas chromatography using an electron capture detector (ECD). The limit of detection for prodiamine and its 6-amino-imidazole metabolite was 0.01 ppm for each compound. Confirmatory analyses are conducted by Gas Chromatography equipped with a Mass Selective Detector (GC/MSD).

2. SAFETY

- 2.1 The oral LD₅₀ of prodiamine in rats is greater than 5000 mg/kg.
- 2.2 Normal laboratory precautions are adequate for safe handling of prodiamine.
- 2.3 Pentane, ethyl ether, methanol, and toluene are flammable and should not be used near heat, sparks, or open flames.

- 2.4 All solvents should be used only in well ventilated laboratories.
- 2.5 Protective gloves should be worn during extraction and analysis.
- 2.6 Disposal of samples and standards must be done in compliance with on-site safety policies and procedures.

3. MATERIALS/METHODS

3.1 Apparatus

- 3.1.1 Water bath; 60°C.
- 3.1.2 Bottles, screw cap with polyseal[®] liner; 32 oz.
- 3.1.3 Bottles, screw cap with polyseal[®] liner; 8 oz., amber.
- 3.1.4 Centrifuge, International Equipment Company, Model CS.
- 3.1.5 Chromatographic columns; 15 mm (I.D.) x 45 cm with teflon stopcock and water jacket, Lab-Crest Scientific Division, 1531 County Line Rd., Warminster, PA 18974.
- 3.1.6 Concentrator, Kuderna Danish, 125-mL.
- 3.1.7 Condensers, Vigreux.
- 3.1.8 Dish, Pyrex, 190 x 100 mm.
- 3.1.9 Distillation receivers, 15-mL.
- 3.1.10 Evaporator, N-Evap with 40°C water bath, (Organomation Assoc.).
- 3.1.11 Funnels, separatory, 500-mL.
- 3.1.12 Drying Oven; 150°C.
- 3.1.13 Pipets, Pasteur, 9", disposable.
- 3.1.14 Platform Shaker.
- 3.1.15 Rotary Evaporator, Büchi, RotoVapor-RE
- 3.1.16 Round Bottom Flask; 250-mL.

3.2 Reagents

- 3.2.1 Dichloromethane; residue analysis grade.
- 3.2.2. Ethyl ether containing 2% ethanol preservative; residue analysis grade.

- 3.2.3 Silica gel 60 (70-200 mesh), EM Reagents, MC/B Manufacturing Chemists, Inc., 2909 Highland Ave., Cincinnati, OH 45212. 3% Deactivated; Spread a layer about 1" deep in a glass dish and activate in a drying oven for 25 hours at 250°C. Cool in a tightly capped bottle. Add 3g deionized water to 97g of silica gel and shake overnight.
- 3.2.4 Methanol; residue analysis grade.
- 3.2.5 Pentane; residue analysis grade.
- 3.2.6 Sodium chloride; reagent grade.
- 3.2.7 Sodium sulfate; anhydrous, granular, reagent grade.
- 3.2.8 Toluene; reagent grade.

3.3 Preparation of Standard Solutions

- 3.3.1 Prodiamine, (N³, N³-Dipropyl 2,4-Dinitro-6-Trifluoromethyl-m-Phenylene-diamine) - Sandoz Crop Protection Analytical Reference Standard.
- 3.3.2 6-Amino-Imidazole, (6-amino-2-ethyl-7-nitro-1-propyl-5-trifluoromethylbenzimidazole) - Sandoz Crop Protection Analytical Reference Standard.
- 3.3.3 Prodiamine and 6-amino-imidazole are light sensitive. Aliquots of 10⁻¹⁰ g/μL prodiamine standard in toluene showed 15-20% loss of prodiamine after 90 hrs. exposure to ambient laboratory light. Standard solutions and sample extracts must be stored in actinic, amber or foil wrapped glassware. Regular glassware may be used during sample extraction and clean-up, but if processing is interrupted, the extracts must be protected from light.
- 3.3.4 Accurately prepare solutions containing 100.0 mg prodiamine/100 mL toluene and 100.0 mg 6-amino-imidazole/100 mL toluene in actinic or foil wrapped 100-mL volumetric flasks. This gives stock solutions of 10⁻⁶ g/μL, (1ug/μL).
- 3.3.5 Transfer 1.0 mL of each stock solutions (10⁻⁶ g/μL) to separate 100-mL volumetric flasks and bring to the mark with toluene. These solutions (10⁻⁸ g/μL) are used for fortification of recovery samples.
- 3.3.6 Prepare a range of standards for GC/ECD and GC/MSD quantitation by diluting aliquots of the fortification standards to the mark in 100-mL volumetric flask, with toluene as described below:

| <u>Volume of 10⁻⁸ g/μl Solution</u> | <u>Concentration of Final Solution</u> |
|--|--|
| 10 mL | 10 ⁻⁹ g/μL |
| 5 mL | 5 x 10 ⁻¹⁰ g/μL |
| 1 mL | 10 ⁻¹⁰ g/μL |
| 500 μL | 5 x 10 ⁻¹¹ g/μL |
| 200 μL | 2 x 10 ⁻¹¹ g/μL |
| 100 μL | 10 ⁻¹¹ g/μL |

A set of standards for prodiamine and a set of standards for 6-amino-imidazole are produced in this way.

3.4 Procedure

3.4.1 Extraction

- 3.4.1.1 Weigh 20 g of sample in an 8 oz. screw-cap bottle.
- 3.4.1.2 If fortifying for recovery data, add an appropriate volume of both prodiamine and 6-amino-imidazole fortification standard solutions to the sample, e.g., 200 μL of 10⁻⁸ g/μL solution to 20 gm of sample = 0.1 ppm fortification. Allow 15 minutes for solvent to evaporate.
- 3.4.1.3 Add 200 mL of methanol to the 8 oz. bottle. Cap the bottle and shake for 30 minutes on a platform shaker.
- 3.4.1.4 Centrifuge for 30 minutes at 550 G or until the supernatant is clear and sediments are stable enough to allow supernatant to be decanted.
- 3.4.1.5 Store extracts in 8-oz. amber bottles in the refrigerator.

3.4.2 Partition I

- 3.4.2.1 Decant a 50-mL aliquot (5 g equivalent) of sample extract and transfer it to a 500-mL separatory funnel containing 250 mL of 5% NaCl solution and 25 mL of dichloromethane.
- 3.4.2.2 Shake for 1.0 minute. Allow phases to separate, and drain the dichloromethane through anhydrous sodium sulfate into a 250-mL round bottom flask.
- 3.4.2.3 Extract the aqueous phase two more times with 25 mL of fresh dichloromethane each time, as described in 3.4.2.2. Combine all three extracts in the 250-mL round bottom flask.

- 3.4.2.4 Wash the sodium sulfate three times with 5-10 mL of dichloromethane each time, combining all washes with the previous dichloromethane extracts.
- 3.4.2.5 Evaporate the combined dichloromethane to near dryness on the rotary evaporator using a 40°C water bath.
- 3.4.2.6 Evaporate the sample from 3.4.2.5 to dryness with a gentle stream of nitrogen. Add 5 mL of 10% ethyl ether in pentane to the 250-mL round bottom flask and thoroughly dissolve the residue. The sample is now ready for silica gel column clean-up.

3.4.3 Silica Gel Column

- 3.4.3.1 To a 250-mL separatory funnel containing 70 mL of 10% ethyl ether in pentane, slowly add 20 g of 3% water deactivated silica gel. See section 3.2.3 for deactivation procedure.
- 3.4.3.2 Shake well, and quickly drain into a chromatographic column plugged with glass wool.
- 3.4.3.3 Rinse the separatory funnel with 10 mL of 10% ethyl ether in pentane and add to the column.
- 3.4.3.4 When the silica gel is completely settled, add 1.0 cm of granular Na_2SO_4 and drain the solvent to just above the top of the Na_2SO_4 layer.
- 3.4.3.5 Apply the 5-mL 10% ethyl ether/pentane solution from 3.4.2.6, which contains the extracted residue, to the silica gel column. Allow this solution to drain to the top of the Na_2SO_4 before addition of any other solvent.
- 3.4.3.6 Rinse the 250-mL round bottom flask twice with 5-mL portions of 10% ethyl ether in pentane and transfer each rinse to the column, allowing the previous rinse to drain to the top of the Na_2SO_4 before adding the next.
- 3.4.3.7 Pass an additional 70 mL of 10% ethyl ether in pentane through the column. Discard this eluate.
- 3.4.3.8 Elute the prodiamine with 75 mL of 10% ethyl ether in pentane, collecting the eluate in a KD set-up (15 mL distillation receiver connected to a Kuderna Danish concentrator).
- 3.4.3.9 Add 1 mL of Hexane to the eluate from 3.4.3.8 (containing prodiamine). Fit a Vigreux

52

condenser to the KD set-up and concentrate the sample to about 1 mL using a 60°C water bath. Perform this concentration in an efficient hood.

- 3.4.3.10 Evaporate the last traces of solvent using a gentle stream of nitrogen. Add 5.0 mL of toluene to the 15-mL distillation receiver and thoroughly dissolve the residue.
- 3.4.3.11 Cover the distillation receiver with foil to protect the sample from light. Store the sample in the refrigerator if quantitation is not performed immediately.
- 3.4.3.12 Further elute the column with 100 mL of 50% Ethyl Ether in Pentane. Discard this eluate.
- 3.4.3.13 Elute the 6-amino-imidazole with 100 mL of ethyl ether collecting this eluate in a 250-mL round bottom flask.
- 3.4.3.14 Evaporate the 6-amino-imidazole fraction to near dryness on the rotary evaporator using a 40°C water bath.
- 3.4.3.15 Evaporate the last traces of solvent using a gentle stream of nitrogen. Add 5.0 mL of toluene to the 250-mL round bottom flask and thoroughly dissolve the residue.
- 3.4.3.16 Protect this sample from light by wrapping the tube with foil. Store the sample in the refrigerator if quantitation is not performed immediately.

3.5 Analysis

3.5.1 Gas Chromatographic Conditions

The following conditions have been shown to be suitable for analysis of Proflamine and the 6-amino-imidazole metabolite in soil. Other conditions may be acceptable provided that the analytes are separated from sample interferences and the response is linear over the range of interest. Elution time and standard linearity must be checked with a new instrument or any change in operating parameters.

- 3.5.1.1 Instrument: Hewlett-Packard 5880A with ⁶³Ni ECD
Column: 10M x 0.53 mm (I.D.) Fused silica, HP-17 (crosslinked 50% phenyl/methyl silicone), 2.0 µm film thickness (H-P #19095L-121)

Oven Temperature Profile:

Initial Value: 185°C
Initial Time: 6 minutes
Post Value: 225°C
Post Time: 5 minutes

Detector Temperature: 350°C
Injector Temperature: 250°C
Carrier Gas: He @ 6 PSI (17 mL/min)
Make-up Gas: 5% methane/argon @ 18 mL/min
Injection Mode: Splitless, 0.5 min purge
delay, split outlet-28 mL/min
Proflamine Retention Time: 2.95 min
6-Amino-Imidazole Retention Time: 4.56 min

3.5.1.2 Instrument: Hewlett-Packard 5880A with 5970 MSD
Column: 25M x 0.2 mm (I.D.) fused silica, HP-1
(Crosslinked methyl silicone), 0.11 µm
film thickness (H-P part number
190912-002)

Oven Temperature Profile:

Initial Value: 100°C
Initial Time: 0.5 min.

Level 1

Program Rate: 30°C/min.
Final Value: 190°C
Final Time: 5.5 min
Post Value: 250°C
Post Time: 5 minutes

Injector Temperature: 250°C
Interface Temperature: 250°C
Carrier: He @ 5 PSI (0.3 mL/min.)

Injection Mode: Splitless, 0.5 min. purge
delay, split outlet - 60 mL/
min.

Source vacuum: 8×10^{-6} TORR
Proflamine Retention Time: 7.25 min.
6-Amino-Imidazole Retention Time: 8.11 min.

3.5.2 Quantitation

3.5.2.1 Prepare a standard curve by injecting 2 µL
aliquots of standards of known concentration and
plotting peak height versus concentration of
injected standard on log-log paper.

3.5.2.2 Determine the concentration of analyte in a 2 µl

injected aliquot of sample from the peak height and the standard curve from 3.5.2.1..

- 3.5.2.3 Calculate the concentration of the residue in the sample using the following expression:

$$\text{PPM (ng/mg)} = \frac{\text{Ce(ng/UL)} \times \text{Vs(UL)}}{\text{Ws(mg)}}$$

Where:

PPM = Concentration of analyte in the sample, in parts per million (ng/mg).

V_s = Final volume of sample extract taking into account all dilutions, (in microliters).

C_e = Concentration of Analyte in sample extract, determined from standard curve, (in nanograms per microliter).

W_s = Weight of sample represented by sample extract, (in milligrams).

3.6 Confirmatory Techniques

The presence of proflamine and its 6-amino-imidazole metabolite can be confirmed using a G.C. equipped with a mass selective detector. The conditions described in 3.5.1.2 have been used for GC/MSD quantitation. Autotune values were used for the repeller, lenses, and electron multiplier voltages. Data was acquired in SIM mode with dwell times of 500 μ SEC for each of three ions. Quantitation of proflamine was based on ion 321 and confirmation of identity was based on the ratios of ions 279 and 333 to 321. Quantitation of the 6-amino-imidazole was based on ion 316 and confirmation of identity was based on the ratios of ions 239 and 228 to 316.

3.7 Time Required for Analysis

A single sample can be extracted and prepared for analysis by gas chromatography in a single eight hour period. A set of six samples required two eight hour days to prepare for gas chromatography. Samples were analyzed by gas chromatography using an autosampler. Using the conditions in 3.5.1.1, each analysis required about 15 minutes.

4. RESULTS AND DISCUSSION

4.1 Precision and Accuracy

Recoveries of proflamine and imidazole from fortified soil samples are shown in Table I. The average recovery was 94%, (S.D.=3.5, N=6) for proflamine and 77%, (S.D.=3.6, N=5) for 6-amino-imidazole.

4.2 Limits of Detection

The limits of detection for both proflamine and 6-amino-imidazole was 0.01 ppm by ECD. This limit of detection was determined by the sample dilution (5 g equiv./5 mL) and the lowest standard used (10^{-1} g/ μ L). Check samples were generally very clean, and the limit of detection may be lowered by concentrating the sample or extending the standard curve to lower levels with more dilute standards.

5. CONCLUSION

This method detects and quantitates proflamine and its 6-amino-imidazole metabolite in soil to a limit of detection of 0.01 ppm with excellent recoveries. This method has not been validated with a variety of soils, and specific soils, particularly those with high organic matter content, may give interferences or lower recoveries.

6. CERTIFICATION

I hereby state as author and the chemist of record that the method described herein was conducted within the framework of the GLP program at Sandoz Crop Protection Corporation. The description of this method and the supporting data (recoveries) are accurate and correct to the best of my knowledge.

KEITON SMITH
name

[Signature]
signature

SCIENTIST II
title

4/5/88
date

name

signature

title

date

I hereby state as the Study Director for the analysis and reporting portion of the above method that the contents herein are accurate and correct to the best of my knowledge.

Thomas Bode
name

[Signature]
signature

Section Leader / Residue
title


4/5/88
date

Table I. Recoveries of Proflamfene and Its 6-Amino-Imidazole Metabolite from Fortified Soil.


| Sample | Fortification Level | Proflamfene Recovery | Imidazole Recovery |
|------------------------|---------------------|----------------------|--------------------|
| B-1 CK ¹ | 0.10 | 91% | 73% |
| A-1 CK ¹ | 0.10 | 97% | 81% |
| A-2 CK ¹ | 0.10 | 94% | 74% |
| A-1 CK ¹ | 0.10 | 96% | 78% |
| B-3 CK ¹ | 0.10 | 89% | 80% |
| B02 CK ¹ | 0.10 | 98% | 58% ² |
| Average (x) | | 94.2 | 77.2 |
| Std. Deviation (s.d.): | | 3.5 | 3.6 |

¹ Recoveries from Report No. 13, Project 480428

² This value is an outlier and was not used to calculate the Average or Standard Deviation.



N,N'-dipropyl-2,4-dinitro-6-(trifluoroethyl)-
1,3-benzenediamine (Procliamina)



6-amino-2-ethyl-7-nitro-1-propyl-5-trifluoromethylbenzimidazole
(6-Amino-Imidazole)

Figure 1. Chemical Structures

SOIL ON

PROGRAMME
SOIL DISSIPATION

01-JUN-67 0144 Page 1





Figure 2. Typical Standard Curve for Proflambine
Correlation = 0.9991

—□— 1.0230020X - 0.004010

CHART SPEED 1.0 CM/MIN
ATTEN: 12E ZERO: 5% 1 MIN/TICK



TITLE: SOIL DISSEMINATION (G-10 P203)

12:25 7 MAR 68

CHANNEL NO: 4

SAMPLE: E-10

METHOD: PRO

| PKT | PKT NAME | RESULT FACTOR | TIME (MIN) | TIME OFFSET | HEIGHT COUNTS | SEP CORR | WGT |
|---------|----------|---------------|------------|-------------|---------------|----------|-------|
| 1 | PRO | 0.009367L | 2.946 | -0.000 | 105527 | 0.000 | 0.000 |
| TOTALS: | | | | -0.000 | 105527 | | |

UNIDENT AREA: 0

DETECTED PKs: 1 REJECTED PKs: 0


AMT STD: 0.10000

NOISE: 25.5 OFFSET: 629

Figure 3. Typical chromatogram of prodiamine standard, 2.0 µL of 0.1 ng/µL standard (0.2 ng) injected.

177/11

ART SPEED 1.0 CM/MIN
PEN: 128 ZERO: 5X 1 MIN/TICK



TITLE: SOIL DISSEMINATION (CO-10 FILED)

12:52 017 MAR 88

CHANNEL NO: 4

SAMPLE: 452 CM

METHOD: PFC

| PEAK NO | RET. TIME (MIN) | AREA | HEIGHT | WIDTH | OFFSET | NOISE | RES | SYN | INT |
|---------|-----------------|------|--------|-------|--------|-------|-----|-----|-----|
| 4 | 2.241 | 0 | 0 | 0 | 0.014 | 0 | 0 | 0 | 0 |

TOTALS:

0

0.014

0.0

UNIDENT AREA: 18983

DETECTED PPS: 4

REJECTED PPS: 0

DIVISOR: 5.00000


MULTIPLIER: 5.00000

NOISE: 25.5

OFFSET: 556

Figure 4. Typical chromatogram of check soil, 2.0 mg equiv. inj., <0.02 ns prodiamine detected (<0.01 ppm).

CHART SPEED 1.0 CM/MIN
ATTEN: 120 ZERO: 5% 1 MIN/TICK




TITLE: SOIL DISSEMINATION (0-10 FILE)

13:05 7 MAR 88

| CHANNEL NO: | SAMPLE: | METHOD: | | | | | |
|---------------|-----------|---------|------------|---------------|---------------|----------|------------|
| 4 | 4-2 CH | PRO | | | | | |
| PEAK NO | PEAK NAME | RESULT | TIME (MIN) | TIME OFFSET | WEIGHT COUNTS | SEC CODE | SI C (SEC) |
| 4 | PRO | 0 | 2.95c | 0.00c | 101922 | VE | 0.180 |
| TOTALS: | | 0 | | 0.00c | 101922 | | |
| UNIDENT AREA: | | 20305 | | | | | |
| DETECTED PMS: | | 4 | | REJECTED PMS: | | 0 | |
| DIVISOR: | | 5.0000 | | MULTIPLIER: | | 5.0000 | |
| NOISE: | | 25.5 | | OFFSET: | | 854 | |

Figure 5. Typical chromatogram of check soil fortified at 0.1 ppm, 2.0 mg equiv. inj., 0.188 ug procliamine detected (0.094 ppm, 94% recovery).

CHART SPEED 1.0 CM/MIN
ATTEN: 128 ZERO: 5% 1 MIN/TICK



TITLE: SOIL DISSIPATION (0-10 FEET)

CHANNELING: 4 SAMPLE: 111-2 ONLY

| PEAK NO | PEAK NAME | RESULT | TIME (MIN) | HEIGHT | SEE |
|---------|-----------|--------|------------|--------|-----|
| 2 | PRO | 0.56 | 2.953 | 60827 | VE |

TOTALS: 1 -0.007 60827

UNIDENTIFIED: 419

DETECTED PYS: 4 REJECTED PYS: 2

DIVISOR: 5.0000 MULTIPLIER: 50.0000

NOISE: 25.5 OFFSET: 640

Figure 6. Typical chromatogram of treated soil, 0.20 mg equiv. inj., 0.112 ng prodiamine detected (0.56 ppm).

44

01-MAR-68 0130 Page 1

10-20 CM

IN DISS. IN CM





Figure 7. Typical Standard Curve For 6-Amino-Imidazole
Correlation = 0.999733

43

CHART SPEED 1.0 CM/MIN
ATTEN: 128 ZERO: 5% 1 MIN/TICK



BEST AVAILABLE COPY

TITLE: SOIL DISSIPATION (10-20 R2L3)

17:47 29 FEB 68

CHANNEL NO: 4 SAMPLE: SE-11

METHOD: IMID2

| PEAK NO | PEAK NAME | RESULT FACTOR | TIME (MIN) | TIME OFFSET | HEIGHT COUNTS | SEP CODE | REJECT |
|---------|-----------|---------------|------------|-------------|---------------|----------|--------|
| 3 | IMIDAZOLE | 0.013275 | 4.576 | 0.00 | 37664 | EE | REJECT |

TOTALS: 0.00 37664

UNIDENT AREA: 3462


DETECTED PKS: 4 REJECTED PKS: 0

AMT STD: 0.05000

NOISE: 93.5 OFFSET: 1069

Figure 8. Typical Chromatogram of Imidazole Standard, 2.0 uL of a 0.05 ng/uL Standard (0.10 ng) Injected.

CHART SPEED 1.0 CM/MIN
ATTEN: 120 ZERO: 5% 1 MIN/TICK



TITLE: SOIL DISSIPATION (10-20 R2L3)

16:38 29 FEB 88

CHANNEL NO: 4

SAMPLE: A-1 CM

METHOD: IMI22

| PEAK NO | PEAK NAME | RESULT | TIME (MIN) | TIME OFFSET | HEIGHT COUNTS | CONC (PPM) | TIME (SEC) |
|---------|-----------|--------|------------|-------------|---------------|------------|------------|
| 2 | IMIDAZOLE | 0.0007 | 4.616 | 0.016 | 519 | BE | 0.10 |
| TOTALS: | | 0.0007 | | 0.016 | 519 | | |

INTEGRATED AREA: 1689

DETECTED PKS: 4 REJECTED PKS: 0

DIVISOR: 5.00000 MULTIPLIER: 5.00000

NOISE: 93.5 OFFSET: 279

Figure 9. Typical Chromatogram of Check Soil, 2.0 mg equiv. injected. <0.02 ng Imidazole Detected (<0.01 ppm).

RECEIVED
FEDERAL BUREAU OF INVESTIGATION
U. S. DEPARTMENT OF JUSTICE

RECEIVED
 FEDERAL BUREAU OF INVESTIGATION
 U. S. DEPARTMENT OF JUSTICE
 DIVISION OF IDENTIFICATION
 WASHINGTON, D. C. 20535
 DATE: 11/11/60
 TIME: 10:00 AM
 FROM: SAC, NEW YORK
 TO: DIRECTOR, FBI
 SUBJECT: [Illegible]

Special Agent in Charge, New York
100-100000-100000