

Cover Sheet for

**ENVIRONMENTAL CHEMISTRY METHOD**

***Pesticide Name:*** Endosulfan

***MRID #:*** 414686-01

***Matrix:*** Soil/Water

***Analysis:*** GC/ECD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at [flynt.elizabeth@epa.gov](mailto:flynt.elizabeth@epa.gov).

Document: Produktentwicklung, GB-C  
Ökologie  
Dr. E. Dorn  
G. Schuld

414686-01  
Date: 01.09.1989  
Report No.: CH040/87  
Page: 50 (284)

Hoechst 

Document: Analytisches Laboratorium  
Author: H.-J. Werner  
C. Klante  
H. D. Herz

Document:   
Date: 27.01.87  
Report No.: AL60/86  
Page: 0

*Alh.*

Doc 002678 

Racem determination of the active  
ingredient and Endosulfan-sulfate  
in soil, water, urine and plant material  
as well as of Endosulfan-diol  
and Endosulfan-lactone  
in soil, water and urine.

#### Outline of method

The active ingredients and its metabolites are extracted from the sample material with acetone. After dilution of the extract with sodium chloride solution, reextraction is achieved with dichloromethane. After a clean up on a Sie-heads S-X and a mini silica gel column the gaschromatographic determination is carried out using an electron capture detector.

Comment: Produktentwicklung, GB-C  
Oecologie  
Dr. E. Dorn  
G. Schulz

Date: 01.09.1989  
Report No.: CR040/87  
Page: 51 (284)

Hoechst 

Comment: Analytisches Laboratorium  
Author: E.-J. Werner  
C. Klauze  
H. D. Herz

Internal Reference Number  
Date: 27.01.87  
Page No.: AL 60/86  
Page: 1 (21)

Translation of Doc. No. A34468

Hoe 002671 (endosulfan)


Residue determination of the active  
ingredient and endosulfan-sulfate  
in soil, water, urine and plant material  
as well as of endosulfan-diol  
and endosulfan-lactone  
in soil, water and urine.

1. Introduction

1.1 Hoe 002671 (endosulfan)

Chemical name : 6,7,8,9,10,10-hexachloro-1,5,5a,6,9,9a-hexahydro-  
6,9-methano-2,4,3-benzodioxathiopyri-3-oxide  
(IUPAC)

structural formula:



molecular formula :  $C_{10}H_6Cl_6O_3S$   
molar mass : 406.9 g/mol  
solubility : low solubility in water  
readily soluble in dichloromethane,  
benzene, toluene, acetone and alcohols  
other properties : unstable in alkaline media  
endosulfan exists in two stereoisomeric  
forms, the  $\alpha$ - (Hoe 052618) and  $\beta$ -  
endosulfan (Hoe 052619).  
The active ingredient has a ratio of  
1 : 1 ( $\alpha$  :  $\beta$ ).

Document: Produktentwicklungs, GB-C  
Önologie  
Dr. E. Dorn  
G. Schulz

Date: 01.09.1989  
Report No.: CR040/87  
Page: 52 (284)




Document: Analytisches Laboratorium  
H.-J. Werner  
G. Klauke  
H. D. Marx

Date: 27.01.87  
Report No.: AL 60/86  
Page: 2 (21)

1.2 See 031327 (endosulfan-sulfate)

Chemical name : 6,7,8,9,10,10-hexachloro-1,3,5a,6,9,9a-  
hexahydro-6,9-oxathia-2,4,3-benzodioxo-chicopin-3,3-  
dioxide (IUPAC)

structural formula:



molecular formula :  $C_{12}H_6Cl_6O_6S$   
molar mass : 422.9 g/mol

solubility : low solubility in water  
readily soluble in dichloromethane,  
benzene, toluene, acetone and alcohols

Department: Produktentwicklung, GB-C  
Oktalopie  
Dr. E. Dorn  
G. Schuld

Date: 01.09.1989  
Report No.: C9040/87  
Page: 53 (284)

Hoechst 

Department: Analytisches Laboratorium  
H.-J. Werner  
G. Klante  
H. D. Marx

Date: 27.01.87  
Report No.: AL 60/86  
Page: 3 (21)

1.3 Hoe 051328 (endocyclic-lactams)

Chemical name : 4,5,6,7,8,8-hexachloro-1,3,3a,4,7,7a-  
hexahydro-4,7-methano-isobenzofuran-  
1-one (IUPAC)

structural formula:



molecular formula :  $C_9H_4Cl_6O_2$

molar mass : 356.8 g/mol

solubility : low solubility in water  
readily soluble in toluene, acetone and  
alcohols

Produktentwicklung, GB-C  
Oxologia  
Dr. E. Dorn  
G. Schulz

Date: 01.09.1989  
Report No.: CR040/87  
Page: 54 (284)

Hoechst 


Analytisches Laboratorium  
H.-J. Berner  
G. Klante  
H. D. Herz

Date: 27.01.87  
Report No.: AL 60/86  
Page: 4 (21)

1.4 Noe 031329 (endosulfan-diol)

Chemical name : 1,4,5,6,7,7-hexachloro-bicyclo-(2,2,1)-  
hept-5-ene-2,3-dimethanol (IHDAC)

structural formula:



molecular formula :  $C_9H_8Cl_6O_2$

molar mass : 360.8 g/mol

solubility : low solubility in water  
readily soluble in dichloromethane, toluene,  
acetone and alcohols



## 2. Outline of method

The active ingredients in its metabolites are extracted from the sample material with acetone. After dilution of the extract with sodium chloride solution, reextraction is made by dichloromethane. After clean up on a Bio-Beads S-X3 and a mini silica gel column gaschromatographic determination is carried out using an electron capture detector.

## 3. Apparatus

beakers 400 ml  
Ultraaturax (homogenizer)  
measuring cylinder 100 ml  
vacuum glass fritted filter funnel  
separation funnels 250, 1000 ml  
volumetric flasks 2, 25 ml  
funnel  
round bottom flasks 100, 500 ml  
volumetric pipette 5 ml  
GPC-Autozap  
GPC-column 400 mm x 25 mm filled with Bio Beads S-X3 (200-400 mesh)  
soaked with 50 % ethyl acetate/50 % cyclohexane (v/v)  
chromatographic column 25 mm x 7 mm  
pasteur pipettes  
rotary evaporator (bath temperature 40°C)  
gaschromatograph equipped with capillary and ECD  
injection bottles  
on column syringe  
disposable syringe 10 ml



Department: Analytisches Laboratorium  
 H.-J. Werner  
 G. Klauke  
 H. D. Herz

Date: 27.01.87  
 Report No.: AL 60/86  
 Page: 6 (21)

4. Reagents

- acetone Pestanal (Riedel de Haen)
- sodium chloride R.G. (Merck)
- sodium chloride solution (4 g/100 ml H<sub>2</sub>O)
- dichloromethane (coffee quality)
- Calite R.G.
- ethyl acetate, dist.
- cyclohexane R.G. (Riedel de Haen)
- 50 % ethyl acetate/50 % cyclohexane (v/v)
- n-hexane for residue analysis (Merck)
- quartz wool
- silica gel 60 (Merck 7734)
- silica gel deactivated with 1.5 % of water
- Bio-Beads S-X3 (200-400 mesh) (Bio-Rad Laboratories)
- sodium sulfate anhydrous R.G. (Merck)
- toluene, clean up by silica gel and distillation
- 95 % toluene/5 % acetone (v/v)
- MSFA N-ethyl-N-triethylsilyl-trifluoroacetamide (Macherey & Nagel)
- D-endosulfan Pestanal (Riedel de Haen)
- S-endosulfan Pestanal (Riedel de Haen)
- endosulfan-sulfate Pestanal (Riedel de Haen)
- endosulfan-lactone Pestanal (Riedel de Haen)
- endosulfan-diol Pestanal (Riedel de Haen)

5. Sampling

Aliquots from the laboratory sample reduced to small pieces are taken for analysis. Liquids are mixed well.

14-10-80/2



Hoechst 

Analitisches Laboratorium  
M.-J. Werner  
C. Klante  
H. D. Marx

Date: 27.01.87  
Page No.: AL 60/86  
Page: 7 (21)

## 6. Procedure

### 6.1 Extraction

#### 6.1.1 Plant material

The analysis of plant material was tested only with the active ingredient and the metabolite endosulfan-sulfate.

50 g (V) of sample is weighed into a beaker. 70 ml of acetone are added and the sample is macerated thoroughly with an Ultraturrax. The macerate is filtered through a glass fritted filter funnel. The extract is transferred directly into the separation funnel (1000 ml) containing 600 ml sodium chloride solution. The extraction step is repeated twice with the filter cake.

If problems arise during the filtration process, Celite might be used as filter aid.

#### 6.1.2 Soil samples

50 g (V) of sample are weighed into a separation funnel (250 ml). 70 ml of acetone are added and the sample is shaken for half an hour on a mechanical shaker. The sample is filtered through a glass fritted filter funnel. The extract is transferred directly into a separation funnel (1000 ml, containing 600 ml of sodium chloride solution. The extraction step is repeated twice with the filter cake.

Department: Produktentwicklung, GB-C  
Ökologie  
Dr. E. Dorn  
G. Schuld

Date: 01.09.1989  
Report No.: CR040/87  
Page: 58 (284)

Hoechst 

Department: Analytisches Laboratorium  
H.-J. Verner  
G. Klancs  
H. D. Marx

Date: 27.01.87  
Report No.: AL 60/86  
Page: 8 (21)

#### 6.1.3 Samples of surface waters

100 ml (V) of water are measured into a 250 ml separation funnel.

#### 6.1.4 Urine samples

30 ml (V) of urine are measured into a separation funnel (1000 ml) containing 200 ml of acetone and 600 ml of sodium chloride solution.

#### 6.2 Isolation of the active ingredient and metabolites

100 ml of dichloromethane are added to the materials resulting from step 6.1.1, 6.1.2 and 6.1.4. 100 ml of dichloromethane are added to the sample from 6.1.3. The separation funnel is shaken for 1 min. After separation of the phases the lower one is drained off over a funnel containing sodium sulfate and collected in a round bottom flask. The extraction is repeated twice with 50 ml of dichloromethane. The dichloromethane is drained off as well over the sodium sulfate and collected. The sodium sulfate is rinsed with some small portions of dichloromethane. The solvent is evaporated nearly to dryness using a rotary evaporator.

Hoechst 

6.2.1 To extracts from soil-, water- and urine-samples some n-hexane is added. The sample is once more evaporated to near dryness. The reconstitution is repeated until there is no dichloromethane left. The residue is transferred quantitatively with n-hexane into a 25 ml ( $V_1$ ) volumetric flask, filled up to the mark and shaken.

6.2.2 To extracts from plant materials (6.1.1) a mixture of 50 % ethyl acetate/50 % cyclohexane (v/v) is added. The sample is once more evaporated to near dryness. The reconstitution is repeated until there is no dichloromethane left. The residue is transferred quantitatively with a mixture of 50 % ethyl acetate/50 % cyclohexane (v/v) into a 25 ml ( $V_1$ ) volumetric flask, filled up to the mark, shaken and the suspended precipitate is allowed to settle.

A 10 ml aliquot of this solution is cleared by gel permeation as described under 6.3.1.

6.3  Clean up

6.3.1 Gel permeation

The conditions for equipment (dump and collect volumes) are to be determined prior to analysis by testing the column with a standard solution.

Suck 10 ml of the solution into a disposable syringe and inject into the 5 ml loop ( $T_1$ ) of the GPC. While the sample passes through the column the cleaned extract is collected in a 100 ml round bottom flask. The solvent is evaporated to near dryness using a rotary evaporator. Some n-hexane is added and the sample is once more evaporated to near dryness. The reconstitution is repeated until there is only n-hexane in the round bottom flask.

Document: Produktentwicklung, GB-C  
Ökologie  
Dr. E. Dorn  
G. Schuld

Date: 01.09.1989  
Report No.: CRD40/87  
Page: 60 (284)

Hoechst 

Document: Analytisches Laboratorium  
H.-J. Werner  
C. Klante  
H. D. Herz

Date: 27.01.87  
Report No.: AL 60/86  
Page: 10 (21)

### 6.3.2 Mini silica gel column

Some glass wool is plugged into the lower end of the chromatographic column. The column is filled with n-hexane. 1 g of deactivated silica gel is added into the column while tapping with a ruler. After sedimentation of the silica gel, a layer of 1 cm of sodium sulfate is added on top. The column is rinsed with 10 ml of n-hexane.

The sample from 6.3.1 is transferred in small portions onto the silica gel column each time draining the solvent until it reaches the upper part of the sodium sulfate layer.

From the n-hexane extracts (6.2.1) 3 ml ( $V_1$ ) are pipetted onto the silica gel column and the n-hexane is drained off until it reaches the upper part of the sodium sulfate layer.

The column is rinsed with n-hexane until a volume of 10 ml (transfer and liquid for rinsing) is reached.

If only  $\alpha$ -endosulfan,  $\beta$ -endosulfan and endosulfan-sulfate are to be determined the column is eluted with 10 ml toluene.

If additionally endosulfan-lactone and endosulfan-diol are to be determined in soil, water and urine all 5 compounds are eluted with 20 ml of 95 % toluene/5 % acetone (v/v). The eluates are evaporated to near dryness using a rotary evaporator and made up to a volume of 2 ml ( $V_2$ ) with toluene.

Department: Produktentwicklung, GB-C  
Ökologie  
Dr. E. Dorn  
G. Schuld

Date: 01.09.1989  
Report No.: CR040/87  
Page: 61 (284)

Hoechst 

Department: Analytisches Laboratorium  
H.-J. Warber  
C. Klante  
H. D. Herz

Date: 27.01.87  
Report No.: AL 60/86  
Page: 11 (21)

#### 6.4 Gaschromatographic determination

##### 6.4.1 Derivatisation

(This step is only necessary where endosulfan-dial has to be determined in soil, water and urine.)

0.02 ml are taken from the cleaned final extract and are diluted to 0.1 ml in a graduated injection bottle with toluene. The silylation is done after adding 0.04 ml MSTFA into the capped injection bottle and heating to 70-80°C for 15 minutes.

After cooling the sample is brought to a final volume of 1 ml. The injection bottle should not be opened anymore during this step. For sampling a short steel cannula funnel shaped by drilling at the Luer-extension is pierced through the septum of the capped injection bottle. This guide cannula allows to sample by GC-syringe. Before a sample is silylated the batch of MSTFA has to be checked for purity to make sure that the determination of endosulfan and its metabolites is not affected by contaminations. This is done by silylation of the solvent and its gaschromatographic analysis.

Experience has shown not to use the opened bottle with silylation reagent longer than a week.

Contract: Produktentwicklung, GB-C  
Osteologie  
Dr. E. Dorn  
G. Schulz

Date: 01.09.1989  
Report No.: CRD40/87  
Page: 62 (284)

Hoechst 

Contract: Analytisches Laboratorium  
H.-J. Werner  
G. Klante  
K. D. Herz

Date: 27.01.87  
Report No.: AL 60/86  
Page: 12 (21)

#### 6.4.2 GC-conditions

model : C. Erba Model 5160 (Mega)  
column: type: capillary (ICT/ISS)  
stat. phase: DB-5 (0.25  $\mu$ m film)  
length: 30 m  
internal diameter: 0.32 mm  
material: fused silica  
carrier gas: Helium 4.6 l/min  
purge gas : N<sub>2</sub>-special 35 ml/min  
injector : on column injector  
detector : electron capture detector (ECD)  
temperature: injector : 45°C  
detector : 220°C  
column : -  
temperature program: oven: standby 85°C 1 min  
ramp 1 bal. T<sub>1</sub> 200°C 0 min  
ramp 2 20 C/min T<sub>2</sub> 220°C 15 min  
ramp 3 bal. T<sub>3</sub> 250°C 3 min

#### 7. Evaluation

##### 7.1 Method of determination

The evaluation is carried out by determining the peak height ( $F_s$ ) of the sample and comparing it with the peak height ( $F_c$ ) of a control sample, spiked with standard solution. The volume of the injected sample and the volume of the control sample should be identical.

Hoechst 

### 7.2 Lower limit of the practical working range and recovery

The lower limit of the practical working range is 0.01 mg/kg. The recoveries (F) from spiking with 0.01, 0.05 and 0.1 mg/kg ranged from 70 to 100 %.

### 7.3 Calculations

The calculation is carried out in accordance with DFG directives V and XI.

$$N = \frac{V_1 \times V_2 \times V_3 \times V_4 \times N_c \times F_s}{W \times I_1 \times I_2 \times I_3 \times I_4 \times F_c} \times F$$

- N - residue in mg/kg or mg/L
- W - weight or volume of analytical sample (g or mL)
- V<sub>1</sub> - volume of extract (mL)
- I<sub>1</sub> - partial volume for clean up (mL)
- V<sub>2</sub> - final volume (mL)
- I<sub>2</sub> - partial volume for possible dilution (mL)
- I<sub>3</sub> - partial volume for possible dilution (mL)
- V<sub>3</sub> - volume for possible dilution (mL)
- V<sub>4</sub> - volume for possible dilution (mL)
- I<sub>4</sub> - volume of injected sample solution (μL)
- N<sub>c</sub> - amount of standard injected (ng)
- F<sub>s</sub> - peak height of sample (cm)
- F<sub>c</sub> - peak height of standard (cm)
- F - recovery factor

Department: Produktentwicklung, GB-C  
Ökologie  
Dr. E. Dorn  
G. Schulz

Date: 01.09.1989  
Report No.: CRD40/87  
Page: 64 (284)

Hoechst   
A34558

Author: Analytisches Laboratorium  
E.-J. Garner  
G. Klauz  
H. D. Herz

Date: 27.01.87  
Report No.: AL60/86  
Page: 16 (21)

8. Literature  
Deutsche Forschungsgemeinschaft, Rückstandsanalytik von  
Pflanzenschutzmitteln (1. - 8. Lieferung, 1983)  
VCH Verlagsgesellschaft mbH, D-6940 Weinheim, FRG
9. Comments  
No comments

*E. J. Garner*  
Author: E.-J. Garner