

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Fenpropathrin

MRID #: 416796-01

Matrix: Soil

Analysis: GC/TSD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

**CHEVRON CHEMICAL COMPANY
ORTHO AGRICULTURAL CHEMICALS DIVISION
RESIDUE CHEMISTRY LABORATORY
RICHMOND, CALIFORNIA**

**DETERMINATION OF FENPROPATHRIN AND
TWO METABOLITES IN SOIL
RESIDUE METHOD RM-225-3**

**DATE: MAY 30, 1989
FILE NO: 740.01/DANITOL**

INTRODUCTION

This method describes a procedure for the determination of fenpropathrin, [(RS)- α -cyano-3-phenoxybenzyl-2,2,3,3-tetramethylcyclopropane carboxylate], desphenylfenpropathrin, (α -cyano-3-hydroxybenzyl-2,2,3,3-tetramethylcyclopropane carboxylate) and 4'-OH-fenpropathrin, (α -cyano-3-[4-hydroxyphenoxy]benzyl-2,2,3,3-tetramethylcyclopropane carboxylate) in soil. Soil metabolism studies have shown that desphenylfenpropathrin and 4'-OH-fenpropathrin are major degradation products of fenpropathrin. This revision of method RM-225-2 allows the simultaneous analysis for all three moieties and incorporates a new cleanup procedure developed by Sumitomo Chemical Company for the analysis of the two degradates (Ref. 1 and 2). In summary the revised method involves extraction with methanol/0.1N HCl, partitioning of residues into dichloromethane, Florisil column cleanup and measurement by gas chromatography using a nitrogen/phosphorous flame ionization detector.

REAGENTS

Fenpropathrin - Reference Standard

Desphenylfenpropathrin - Reference Standard

4'-OH-fenpropathrin - Reference Standard

Acetone - Pesticide Quality

Methanol - Pesticide Quality

Dichloromethane - Pesticide Quality

Hexanes - Pesticide quality

Ethyl Acetate - Pesticide quality

Hydrochloric Acid - 0.1 N

Celite Analytical Filter-Aid

Sodium Sulfate - Anhydrous, granular, reagent grade, acetone washed and air-dried.

Sodium Chloride - 5% aqueous solution (50 g/liter of water)

Florisil - Activated at 1200°F, mix 3 parts by weight of 60-100 mesh with 2 parts 100-200 mesh. Heat overnight at 110°C prior to deactivation with 10% water (See Note 2).

EQUIPMENT

Omni-Mixer with adaptor for use with 1-pint Mason jars.

Rotary vacuum evaporators equipped with a heated water bath.

Buchner Funnels

Glass-fiber filter paper - Whatman GF/A

Separatory Funnels - 1000 ml, Teflon® stopcocks

Liquid Chromatography Column - 300 x 22 mm, Supelco # 6-4760

Gas Chromatograph (Hewlett-Packard 5890A or equivalent) equipped with both a packed column and a split/splitless capillary injection port, a nitrogen/phosphorous flame ionization detector, a HP-7673A autosampler and a HP-3392A integrator.

Column A: 15 m x 0.53 mm ID 5% phenyl-methyl silicone (1.5 µm film) megabore column (J&W Scientific; DB-5)

Column B: 15 m x 0.53 mm 50% phenyl-methyl silicone (1.0 µm film) megabore column (J&W Scientific; DB-17)

Injector: Packed column with megabore adaptor

Carrier gas flow rate: 30 ml/min (He)

Auxillary gas flow rate: 0-15 ml/min (He)

Detector gases:

Hydrogen: 3.5 ml/min

Air: 90-120 ml/min

Injector temperature: 260-280°C

Detector Temperature: 260-300°C

Column Oven Temperatures:

	<u>Column A</u>	<u>Column B</u>
Initial:	150°C, hold 1.0 min	200°C hold 2.0 min
Rate:	20°C/min	12°C 1.0 min
Final:	270°C, hold for 1 min	260°C, hold for 2 min
Desphenylfenpropathrin:	4.1 min	2.5
Fenpropathrin:	5.7 min	5.2
4'-OH-fenpropathrin:	6.9 min	7.9 min
	(Figure 1)	(Figure 2)

The above parameters serve only as a guide for chromatography of the soil samples. Changes to the operating parameters may be made to optimize the chromatography for a particular instrument or soil matrix.

ANALYTICAL METHOD

All glassware should be rinsed with pesticide-quality acetone prior to use. The acetone-rinsed glassware should then be rinsed with the solvent to be used in that particular step of the procedure.

Extraction

Weigh out a 20 g sample of soil into a pint Mason jar. Fortify a control sample for recovery purposes with 1.0 ml of a 2.0 µg/ml acetone solution of fenpropathrin, desphenylfenpropathrin and 4'-OH-fenpropathrin. Add 100 ml methanol/0.1N HCl (8/2,v/v) and blend on an Omni-mixer for five minutes. Add 1 teaspoon of Celite Analytical Filter-Aid to the sample and filter the extract through a Buchner funnel lined with a Whatman GF/A filter paper. Transfer the filter cake to the Mason jar and repeat the extraction and filtration steps once more using 100 ml of extraction solvent. Following filtration, rinse the filter cake with 50 ml of extraction solvent.

Partition

Transfer the combined extracts to a 1000 ml separatory funnel and add 250 ml of 5% aqueous NaCl solution. Partition the extract with 100 ml of dichloromethane. Filter the dichloromethane extract through sodium sulfate into a roundbottom flask. Repeat the partition and filtration steps with one additional 100 ml portion of dichloromethane. Rinse the sodium sulfate cake with 50 ml of dichloromethane. Evaporate the combined extract to dryness in a 25°C water bath using the rotary vacuum evaporator. Redissolve the sample in 3 ml of hexane/ethyl acetate (3/1, v/v).

Florisil Column Cleanup

Close the stopcock on the liquid chromatography column, plug with a wad of glass wool and add 50 ml of hexane/ethyl acetate (3/1, v/v). Weigh out 15 g of deactivated Florisil and add slowly to the column to avoid air bubbles within the bed. Allow the Florisil to settle, then cap with approximately 1-inch of sodium sulfate. Drain the solvent layer to the top of the sodium sulfate bed.

Transfer the concentrated sample extract to the column with three 3-ml rinses of mixed solvent. Allow each portion to drain to the top of the sodium sulfate layer before adding the next rinse portion. Following the last rinse, add 10 ml of mixed solvent to the column and elute. Discard all column eluates. Place a 100-ml roundbottom flask under the column and elute the sample with 55 ml of hexane/ethyl acetate (3/1, v/v). Evaporate the extract to dryness and redissolve in 2.0 ml of ethyl acetate.

MEASUREMENT

The linearity of the gas chromatographic system should be checked daily by analyzing mixed standards of fenpropathrin, desphenylfenpropathrin and 4'-OH-fenpropathrin. A concentration range of 0.1, 0.5, 1.0 and 2.0 µg/ml is recommended. The mean of the response factors (integration counts/std. conc.) should have a coefficient of variation of ±10% or less. Variance from this requirement must be approved by supervisory personnel.

Inject 1.0-4.0 µl of reference standard solution (1.0 µg/ml fenpropathrin, desphenylfenpropathrin and 4'-OH-fenpropathrin) into the gas chromatograph. Inject the same volume of sample extract. As required, dilute the sample extracts with ethyl acetate to obtain peaks within the linear range of the instrument. A suggested sequence of analysis is standard, sample, sample, standard, sample, sample, standard. . . . The coefficient of variation of the mean of the standards should be ±10% or less. Variance from this requirement must be approved by supervisory personnel.

CALCULATION

$$\text{PPM} = \frac{\text{Integration CTS. (Sample)}}{\text{Ave. Integration CTS. (Std.)}} \times 1 \mu\text{g/ml} \times \frac{2.0 \text{ ml}}{20 \text{ g}} \times \text{Dilution Factor}$$

$$\% \text{ Recovery} = \frac{\text{PPM Fortified} - \text{PPM Control}}{\text{PPM Fortification}} \times 100$$

LIMIT OF DETECTION

The limit of detection is approximately 0.01 ppm for fenpropathrin, desphenylfenpropathrin and 4'-OH-fenpropathrin for a 20 g soil sample analyzed using the described procedure.

NOTES:

1. A fortified control sample must be analyzed concurrently with each set of samples. Method recovery must be between 70-120% to be acceptable. All samples analyzed concurrently with a fortified sample that does not have an acceptable method recovery must be reanalyzed or the requirement waived by supervisory personnel.

2. Preparation of Deactivated Florisil

Weigh out 2 kg of activated Florisil and transfer to a plastic screw-capped jar. Add 200 ml deionized water and mix well. Allow the deactivated Florisil to sit overnight in the capped jar. Prior to use, verify the activity of the Florisil by analyzing a mixed 10 µg standard. Recovery must be 90% or greater for all three compounds.

3. Gas Chromatography Alternatives

Due to the complexity of multiple component analysis, two gas chromatography columns of different polarity are described in the analytical method. If matrix interference cannot be separated on either column for all three components simultaneously, split analyses may be used to analyze for one or more components on one column and the second and/or third component on the alternate column. If a matrix interference cannot be separated on either column, then narrow-bore capillary column chromatography may be used if linearity and reproducibility requirements of the method are satisfied.

REFERENCES

1. Onishi, et. al., Report No. ER-MT08911, "Residue Analytical Method for 4'-OH-fenpropathrin, One of the Metabolites of Fenpropathrin, in Soil," Mar 12, 1989.
2. Onishi, et. al., Report No. ER-MT-8912, "Residue Analytical Method for Desphenyl-fenpropathrin, One of the Metabolites of Fenpropathrin, in Soil," May 12, 1989.

A. L. Hamada

A. L. HAMADA

G. H. Fujie
G. H. FUJIE


Reviewed By: *[Signature]*

GHF:gt

cc: Technical Records

RM-22-3

Figure 1



LIST: LIST
PEAK CAPACITY: 991


ZERO = 10, 68.17
ATT 21 = 7
CHT SP = 1.0
PK WD = 0.04
THRSH = 7
AR REJ = 0

LIST: TIME 0
0.65 INTG # = 9
3.50 INTG # = -9
3.98 INTG # = 0
5.50 INTG # = 0
6.65 INTG # = 0

1.0 µg/ml Mixed Standards
Column A: 15 M x 0.53 mm DB-5 column

RM-225-3

Figure 2



LIST: LIST
PEAK CAPACITY: 991

ZERO = 10, 99.76
ATT 2† = 7
CHT SP = 1.0
PK WD = 0.04
THRSH = 6
AR REJ = 0

LIST: TIME 0
2.68 INTG # = 0
3.85 INTG # = 0
5.80 INTG # = 0
6.50 ATT 2† = 6
7.68 INTG # = 0

1.0 µg/ml Mixed Standards
Column B: 15 M x 0.53 mm DB-17 Column