

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Ethoprop

MRID #: 417124-01

Matrix: Soil

Analysis: GC/FPD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

MORSE LABORATORIES, INC.

CHEMICAL ANALYSIS AND RESEARCH

1828 FULTON AVENUE
SACRAMENTO, CALIFORNIA 95825
TELEPHONE (916) 481-3141
FAX (916) 481-2989

April 20, 1989

RHONE-POULENC, INC.

Rhone-Poulenc Analytical Method No. 172

An Analytical Residue Method for the Determination of
O-Ethyl-S, S-DI-N-Propylphosphorodithioate (Ethoprop)
in Soil by Gas Chromatography (August 1984)

Rhone-Poulenc Study No.: EC/P-89-0012

PROCEDURE:

1. Homogenize the soil sample such that a representative sub-sample may be obtained (remove stones, sticks and other extraneous materials).
 2. Weigh 50 g of soil into a 500 mL erlenmeyer.
- NOTE: See the Rhone-Poulenc Protocol Study No.: EC/P-89-0012 for Soil Moistures. We will not be following Morse Labs SOP for Soil Moistures.
3. Add 100 mL MeOH to the bottle, cover with foil and tape foil around the neck of the erlenmeyer.
 4. Blend on a roller mill for 15 minutes (Morse Labs uses a platform shaker).
 5. Allow soil to settle out. Filter erlenmeyer contents through a buchner funnel with the aid of a mild vacuum. Rinse the erlenmeyer with 30 mL and pour into the buchner funnel.
 6. Discard the filter cake.
 7. Add 5 drops of 10% decanol in acetone to filtrate.
 8. Rotovap at 30-35°C to approximately 30 mL.
 9. Transfer to a 250 mL separatory funnel. Rinse the flask with 20 mL of MeOH and pour into the separatory funnel.

OUR LETTERS AND REPORTS ARE FOR THE EXCLUSIVE USE OF THE CLIENT TO WHOM THEY ARE ADDRESSED, AND THE USE OF THE NAME OF MORSE LABORATORIES, INC. MUST RECEIVE OUR PRIOR WRITTEN APPROVAL. OUR LETTERS AND REPORTS APPLY ONLY TO THE SAMPLE TESTED AND ARE NOT NECESSARILY INDICATIVE OF THE QUALITIES OF APPARENTLY IDENTICAL OR SIMILAR PRODUCTS. SAMPLES NOT DESTROYED IN TESTING ARE RETAINED A MAXIMUM OF THIRTY DAYS. THE REPORTS AND LETTERS AND THE NAME OF MORSE LABORATORIES, INC. IS NOT TO BE USED UNDER ANY CIRCUMSTANCES IN ADVERTISING TO THE GENERAL PUBLIC.

10. Add to the separatory funnel 150 mL of the 10% aqueous sodium chloride solution.
11. Partition the aqueous solution two times with 100 mL portions of hexane.
12. Drain each 100 mL of hexane through Whatman #4 filter paper containing anhydrous sodium sulfate into a 400 mL beaker.
13. Add 5 drops of 10% decanol in acetone.
14. Place on steambath to evaporate. Be sure to place the beaker on top of the watchglass. If the beaker gets too cool, the hexane will pick up H₂O.
15. Evaporate down to 1-2 mL.
16. Transfer the 1-2 mL of hexane to a test tube. Continue the transfer using ETAC.

RESIDUE - RAW DATA SHEET

M.L. Project No. ML-89-0075-RMP
 Sponsor No. _____
 Protocol No. EE/P-89-0012

MORSE LABORATORIES, INC.
 1525 FULTON AVENUE
 SACRAMENTO, CA 95825

Laboratory No. 52791 Page No. 1 of 1
 Product Ellenprop

DATE SAMPLE RECEIVED: 8-4-89 SAMPLES STORED IN FREEZER AT -20 ± 2 °C UNLESS OTHERWISE NOTED:
 DATE SAMPLE GROUND/PREP: NA
 DATE SAMPLE FORTIFIED (STABILITIES): 8-15-89 ANALYTICAL METHOD: Rhone-Poulenc Analytical Method No. 172
 DATE SAMPLE EXTRACTED: 3-30-90 CHEMICAL SPECIES ANALYZED FOR: Ellenprop
 DATE SAMPLE ANALYZED (INSTRUMENTATION): 4-10-90 DETECTION LIMIT: 0.11 ppm/PPB NOTEBOOK: _____ WORKING SHEET: X

SAMPLE IDENTIFICATION	MATRIX	ADDED		RAW ANALYTICAL DATA										PARTS PER MILLION		% RECOVERY		
		ug	P.P.M.	gm of Sample	mL Solvent	mL Alquot	Final Vol.	Dil. Fact.	ul Injct.	PA. H. A. min.	Area. B	Corr. At/B	Net-Gr. (Sigs)	ng Found	mg Injct.		FOUND	CORR.
CK	Soil	-	-	50.0	150	150	5	1	2	0.0	-	0.0	56.2	0.000	20	0.000	0.00	86.
Method Spike		2.5	0.05	50.0	150	150	5	1	2	41.5	-	41.5	56.2	0.163	20	0.143	-	80
Stability spike Dup A		0.5	0.01	50.0	150	150	5	1	2	9	-	9	56.2	0.160	20	0.140	-	80
Stability spike Dup B		0.5	0.01	50.0	150	150	5	1	2	9	-	9	56.2	0.166	20	0.140	-	80
Stability spike Dup A		1.0	0.02	50.0	150	150	5	1	2	19	-	19	56.2	0.230	20	0.0160	-	86
Stability spike Dup B		1.0	0.02	50.0	150	150	5	1	2	19	-	19	56.2	0.211	20	0.0167	-	84
Stability spike Dup A		2.5	0.05	50.0	150	150	5	1	2	42	-	42	56.2	0.237	20	0.037	-	74
Stability spike Dup B		2.5	0.05	50.0	150	150	5	1	2	41	-	41	56.2	0.254	20	0.043	-	56

For soil results only: 1 as rec'd basis
1 as dry weight basis

Analyst: Frank Robinson
 Reviewed by: Doug Peat
Jason Clark

COMMENTS: Residual Stability Study

FILED AVAILABLE COPY