

Cover Sheet for

**ENVIRONMENTAL CHEMISTRY METHOD**

***Pesticide Name:*** Deltamethrin

***MRID #:*** 427739-03

***Matrix:*** Soil

***Analysis:*** GC/ECD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at [flynt.elizabeth@epa.gov](mailto:flynt.elizabeth@epa.gov).

ANALYTICAL METHOD FOR THE DETERMINATION  
OF DELTAMETHRIN, AND DEGRADATES  
IN SOIL BY GAS CHROMATOGRAPHY

CHEMALYSIS METHOD NO. CA004

MRID#  
427739-3

SUBMITTED TO:

HOECHST-ROUSSEL AGRI-VET COMPANY  
ROUTE 202-206  
SOMERVILLE, NEW JERSEY  
08876-1258

HRAV. PROJECT NO.

90-0091

SUBMITTED BY:

CHEMALYSIS INCORPORATED  
8510 CORRIDOR ROAD  
SAVAGE, MARYLAND 20763

STUDY DIRECTOR

ANTHONY F. GRIGOR, PH.D.

LABORATORY PROJECT NO.

900405

INITIATION DATE: July 10, 1990  
INITIATION DATE: July 26, 1991

STUDY #

0173 of 537

173

CONFIDENTIALITY STATEMENT

---

STUDY #

pg 0174 of 537

Analytical Method for the Determination of Deltamethrin, Trans-Deltamethrin and Degradates in Soil by Gas Chromatography

Good Laboratory Practice Statement

This method "Analytical Method for the Determination of Deltamethrin, Trans-Deltamethrin and Degradates in Soil by Gas Chromatography" was conducted according to EPA Good Laboratory Practice Standards, 40 CFR 160.

Anthony F. Grigor

Anthony F. Grigor, Ph.D.  
Study Director

7/26/91

Date

Mary Hagerman

Mary Hagerman  
Chemist

7/26/91

Date



Analytical Method for the Determination  
of Deltamethrin, Trans-Deltamethrin and Degradates  
in Soil by Gas Chromatography

Quality Assurance Monitoring Statement

This method has been audited, confirmed with the raw data and found to be compliance with the specifications for acceptance by Quality Assurance.

The following are dates on which the Quality Assurance unit conducted audits and the dates on which the findings of those audits were reported to management and the Study Director.

<u>Type of Audit</u>	<u>Phase Audited</u>	<u>Audit Date</u>	<u>Date Reported to Study Director &amp; Management</u>
Data	Method Development	1-21-91	1-23-91
		1-29-91	1-31-91
		2-01-91	2-01-91
		2-28-91	3-01-91
		3-01-91	3-05-91
		3-06-91	3-14-91
		3-19-91	3-20-91
		3-24-91	4-09-91
		4-08-91	4-09-91
		4-09-91	4-09-91
		5-23-91	5-29-91

-----  
Quality Assurance Coordinator

-----  
Date

STUDY #

PS 0175#537

Report PR900405



TABLE OF CONTENTS

	<u>Page</u>
CONFIDENTIALITY STATEMENT . . . . .	2
GLP COMPLIANCE WITH ALL PERSONNEL SIGN OFF . . . . .	3
QA SIGN OFF . . . . .	4
LIST OF TABLES AND FIGURES . . . . .	5
1.0 INTRODUCTION . . . . .	6
2.0 PRINCIPLE . . . . .	6
3.0 EQUIPMENT . . . . .	6
4.0 REAGENTS . . . . .	7
5.0 PROCEDURE . . . . .	7
5.1 Standard Preparation . . . . .	7
5.2 Sample Fortification . . . . .	8
5.3 Sample Extraction . . . . .	8
5.4 Deltamethrin and trans-Deltamethrin Clean-up . . . . .	8
5.5 Decamethrinic Acid Partition . . . . .	9
5.6 Analysis Conditions . . . . .	10
5.7 Calculation of Residues . . . . .	11
5.8 Time Requirements . . . . .	12
6.0 References . . . . .	12

## 1.0 INTRODUCTION

This report validates the method used for the determination of residues of deltamethrin, and its degradates trans-deltamethrin and decamethrinic Acid (Br<sub>2</sub>CA) in soil. Analysis is by GC/ECD. The method sensitivity for deltamethrin and trans-deltamethrin is 0.002 ppm. The method sensitivity for decamethrinic acid is 0.010 ppm based on the lowest concentration standard analyzed. The molecular structures of deltamethrin, trans-deltamethrin, decamethrinic acid and its pentafluorobenzyl ester derivative are presented in Figure 1.

### Chemical Information.

#### Deltamethrin, RU 22974

(1R,3R)3-(2,2-dibromovinyl)-2,2-dimethylcyclopropanecarboxylic acid  
(S) alpha-cyano-(3-phenoxyphenyl) methyl ester

#### trans-Deltamethrin, RU 25979

(1R,3S)3-(2,2-dibromovinyl)-2,2-dimethylcyclopropanecarboxylic acid  
(S)-alpha-cyano-(3-phenoxyphenyl) methyl ester

#### Decamethrinic Acid, (Br<sub>2</sub>CA), RU 23441

(1R-cis)-3-(2,2-dibromoethenyl)-2,2-dimethyl-cyclopropanecarboxylic acid

## 2.0 PRINCIPLE

Residues are extracted from soil with 1:1 hexane:acetone, filtered and aliquoted. An aliquot for the determination of deltamethrin and trans-deltamethrin is prepared by using a silica gel cartridge for clean-up. Another aliquot for the determination of decamethrinic acid is taken and derivatized using pentafluorobenzyl bromide (PFB-Br). The two final extracts are analyzed under different GC conditions.

## 3.0 EQUIPMENT

- 3.1 Round bottom flasks, 500 ml, 24/40 ground glass joints
- 3.2 Wide mouth bottles with teflon-lined caps, 500 ml
- 3.3 Volumetric flasks, 500 ml
- 3.4 Horizontal shaker
- 3.5 Buchner Funnels
- 3.6 Filtering Flasks
- 3.7 Whatman Glass Fiber Filters to fit Buchner Funnels
- 3.8 Rotary Evaporators
- 3.9 N.-Evap analytical evaporator
- 3.10 Analytichem Cat. No. 4261304 silica gel cartridges with stopcocks (or equivalent)
- 3.11 Analytichem VAC-ELUT vacuum box for use with cartridges (or equivalent)

STUDY #

pg 0173 of 537

178

#### 4.0 REAGENTS

- 4.1 Hexane, pesticide grade
- 4.2 Acetone, pesticide grade
- 4.3 Chloroform, pesticide grade
- 4.4 Ethyl ether, pesticide grade
- 4.5 Ethyl acetate, pesticide grade
- 4.6  $K_2CO_3$  (30%). Weigh 30 g  $K_2CO_3$  in 100 ml DI water.
- 4.7 Pentafluorobenzyl bromide (5% in acetone). Weigh 5 g PFB-Br in 100 ml acetone. Prepare fresh monthly. PFB-Br is a lachrymator. Prepare and use under a hood.
- 4.8 Silica gel, 5% deactivated. Dry silica overnight. Weigh out 95 g and add 5 ml water. Shake for several hours.
- 4.9 Sodium Hydroxide, 0.1 N. Weigh 4 g sodium hydroxide pellets into 1000 ml volumetric flask. Bring to volume with DI water.
- 4.10 Deltamethrin, trans-Deltamethrin, Decamethrinic Acid and the Penta-FluoroBenzyl Ester of Decamethrinic Acid analytical standards are available from Hoechst-Roussel Agri-Vet Co., Route 202-206, Somerville, NJ 08876-1283.

#### 5. PROCEDURE

##### 5.1 Standard Preparation

Deltamethrin Stock Solution (RU 22974) - Weigh 100.0 mg deltamethrin into a 100 ml volumetric flask. Bring to volume with acetone.

Transfer 1.0 ml of stock solution into a 100 ml volumetric. Bring to volume with hexane. The concentration of this solution is 10.0 ug/ml.

trans-Deltamethrin Stock Solution (RU 25979) - Weigh 100 mg trans-deltamethrin into a 100 ml volumetric flask. Bring to volume with acetone.

Transfer 1.0 ml stock solution into a 100 ml flask. Bring to volume in hexane. The concentration of this solution is 10.0 ug/ml.

Decamethrinic Acid Stock Solution (RU 22441) - Weigh 100 mg decamethrinic acid into a 100 ml volumetric flask. Bring to volume with acetone.

Transfer 1.0 ml stock solution into a 100 ml flask. Bring to volume with acetone. The concentration of this solution is 10.0 ug/ml.

Decamethrinic Acid PFB, Ester (RU 91-2688) - Weigh 100 mg decamethrinic acid PFB ester into a 100 ml volumetric flask. Bring to volume with acetone.

Transfer 1.0 ml stock solution into a 100 ml flask. Bring to volume with hexane.

STUDY #

pg 0179 of 537

170



Standards for GC Analysis - Prepare appropriate dilutions for a mixture of deltamethrin and trans-deltamethrin at 0.20 ug/ml, 0.10 ug/ml, 0.05 ug/ml, 0.02 ug/ml and 0.01 ug/ml in hexane.

Prepare appropriate dilutions of the derivatized standard at 0.20 ug/ml, 0.10 ug/ml, 0.05 ug/ml, 0.02 ug/ml and 0.01 ug/ml in hexane:ether:acetone 98:1:1 (v:v:v).

## 5.2 Sample Fortification

Control soil samples are fortified at 0.002 ppm, and 0.10 ppm with deltamethrin and trans-deltamethrin. Control samples are fortified at 0.01 ppm and 0.10 ppm with decamethrinic acid.

## 5.3 Sample Extraction Clean-Up

- 5.3.1. Weigh 40.0 grams of soil into a 500 ml bottle.
- 5.3.2. Add 300 ml of 1:1 hexane:acetone. Shake on horizontal shaker for 1 hour.
- 5.3.3. Suction filter through a glass fiber filter to remove soil. Rinse filter with 1:1 hexane:acetone. Quantitatively transfer to a 500 ml volumetric flask and bring to volume with 1:1 hexane:acetone. Split filtrate into two 250 ml aliquots into two 500 ml round bottom flasks and evaporate to near dryness by rotary evaporation at 40-45°C.

## 5.4 Deltamethrin and trans-Deltamethrin Clean-Up

- 5.4.1. Rinse and condition a silica gel SPE cartridge (Analytichem LR01304) with 10 ml chloroform followed by 10 ml of hexane.
- 5.4.2. Rinse the flask, step 5.3.3., with 10 ml of hexane and transfer quantitatively to column. Repeat with another 10 ml portion of hexane.
- 5.4.3. Using a vacuum box and stopcock, control the drip rate to 1-2 ml/min. Rinse the column with 10 ml hexane. Do not allow the column to go dry. Discard the hexane.
- 5.4.4. Add 5 ml of chloroform to the silica gel column reservoir and collect the eluate in a graduated test tube.
- 5.4.5. Concentrate the sample to near dryness using an N-Evap at a temperature of 35-40°C.
- 5.4.6. Bring sample to final volume of 2.0 ml with hexane. Transfer 1.0 ml into two GC vials for GC/ECD analysis. This fraction contains deltamethrin and trans-deltamethrin.

STUDY #

PG 0180-537

150

## 5.5 Decamethrinic Acid Partition and Clean-Up

- 5.5.1. Quantitatively transfer residue from step 5.3.3. to a 250 ml separatory funnel with 20 ml of 0.1 N NaOH.
- 5.5.2. Add 20 ml ethyl ether. Shake for 1 minute with venting. Discard ethyl ether.
- 5.5.3. Add 20 ml ethyl ether. Shake for 1 minute with venting. Discard ethyl ether.
- 5.5.4. Add 10 ml of 0.5 N HCl followed by 20 ml ethyl ether.
- 5.5.5. Shake for 1 minute with venting. Collect ethyl ether in a 500 ml round bottom flask.
- 5.5.6. Add 20 ml ethyl ether. Shake for 1 minute with venting. collect ethyl ether in the 500 ml round bottom flask.
- 5.5.7. Evaporate the sample by rotary evaporation at 40-45°C.
- 5.5.8. Add 100 ml acetone, 1 ml of PFB-Br (5% in acetone) and 300 ul of 30% K<sub>2</sub>CO<sub>3</sub>.
- 5.5.9. Reflux gently for 3 hours.
- 5.5.10. Evaporate the sample by rotary evaporation at 40-45°C.
- 5.5.11. Add 250 ml ether, 2 ml DI water and 10 ml ethyl acetate. Swirl vigorously.
- 5.5.12. Filter through 12g anhydrous sodium sulfate prewashed with ether.
- 5.5.13. Evaporate sample to near dryness by rotary evaporation at 40-45°C.
- 5.5.14. Bring sample to final volume of 10 ml in a glass scintillation vial with hexane:ether:acetone (98:1:1).
- 5.5.15. Freeze all samples until analysis.

## 5.6. Analysis Conditions

Analysis by Gas Chromatography for deltamethrin and trans-deltamethrin is accomplished using the following column and conditions:

Instrument: HP 5890 GC, with Autosampler

Column: DB-1, 30 m x 0.322 mm, film thickness 0.25 um

Injector: 200°C, splitless, purge on at 1.0 min., 20 ml/min

Detector: Electron Capture, 325°C, make-up argon/methane, 30 ml/min.

Temp Program: 100°C for 1 min.  
20°C/min.

STUDY #

130°C for 10 min.

ps 0781 of 537

181

10°C/min.  
250°C for 5 min.

Carrier Gas: Helium, 5 ml/min.

Inj. Vol: 2 uL

Under these conditions, approximate retention times are as follows:

cis-deltamethrin 19.1 min.  
trans-deltamethrin 19.5 min.

Analysis by Gas Chromatography for Br<sub>2</sub>CA is accomplished using the following column and conditions:

Instrument: HP 5890 GC, with Autosampler

Column: DB-1, 30 m x 0.322 mm, film thickness  
0.25 um

Injector: 200°C, splitless, purge on at 1 min.

Temp Program: 100°C for 1 min.  
20°C/min.  
230°C for 10 min.  
10°C/min.  
250°C for 5 min.

Carrier Gas: Helium, 5 ml/min.

Inj. Vol: 2 uL

Under these conditions, approximate retention time is as follows:

penta-fluorobenzyl ester of decamethrinic acid: 16.6 minutes

GC conditions may be changed to produce similar results.

Typical standardization data for deltamethrin, trans-deltamethrin and decamethrinic acid PFB ester are shown in Tables 1 and 2. Recovery data for fortified soil are presented in Table 3. Typical chromatograms for standards, control soil and fortified control soil are depicted in Figures 3-5.

#### Separability

The method is sensitive to 0.002 ppm for deltamethrin and trans-deltamethrin and to 0.010 ppm for decamethrinic acid PFB ester based on the lowest concentration standard analyzed.

### 5.7 Calculationn of Residues

$$\text{mg injected} = \frac{W \times V2 \times V4}{V1 \times V3 \times DF}$$

where W = 40.0 grams sample  
V1 = 500 ml initial volume  
V2 = 250 ml aliquot volume  
V3 = 2.0 mL or 10.0 final volume  
V4 = 2 uL injected  
DF = Dilution Factor

$$\text{ppm, residue} = \frac{\text{ng from curve}}{\text{mg injected}}$$

Field samples are corrected for soil moisture such that  
ppm, residue (dry weight basis)

$$= \frac{\text{ppm, residue} \times 100}{(100 - \% \text{ Moisture})}$$

### 5.8. Time Requirements

A trained analyst can complete the extraction of a set of 10-12 samples in 12 hours. A flowchart of the extraction procedures (Figure 2) indicates possible stopping points during extraction.

### 6.0 References

EN-CAS Method No.: ENC-7/89  
Maguire, R. James, et al. J. Agric. Food Chem., Vol 37, Nov. 4, 1989.

STUDY #

0183 of 537

1-8-8

LIST OF TABLES AND FIGURES

	<u>Page</u>
Table 1. Typical Standardization data for Deltamethrin, trans-Deltamethrin and Decamethrinic Acid PFB Ester . . . . .	13
Table 2. Recovery Data for Soil Fortified with Deltamethrin, trans-Deltamethrin and Decamethrinic Acid . . . . .	14
Figure 1. Molecular Structures for Deltamethrin, trans-Deltamethrin, Decamethrinic Acid and Decamethrinic Acid PFB Ester . . . . .	15
Figure 2. Flow Chart of Extraction Procedures . . . . .	17
Figure 3. Representative Chromatograms of Deltamethrin and trans-Deltamethrin . . . . .	19
Figure 4. Representative Chromatograms of Decamethrinic Acid PFB Ester . . . . .	20

Table 1. Typical Standardization Data for Deltamethrin,  
trans-Deltamethrin and Decamethrinic Acid PFB Ester.

Standard Conc.	Peak Height	ng Injected	Regression Parameters
<b>Deltamethrin</b>			
0.01 ug/ml	1592	0.02	Corr. coeff (r) 0.9920 y intercept (b) 899.7 slope (m) 42202.2
0.01 ug/ml	1457	0.02	
0.02 ug/ml	2869	0.04	
0.05 ug/ml	5415	0.10	
0.05 ug/ml	4433	0.10	
0.10 ug/ml	10818	0.20	
0.10 ug/ml	8777	0.20	
0.20 ug/ml	17415	0.40	
<b>trans-Deltametrin</b>			
0.01 ug/ml	1091	0.02	Corr. coeff (r) 0.9931 y intercept (b) 641.0 slope (m) 41756.6
0.01 ug/ml	1303	0.02	
0.02 ug/ml	2668	0.04	
0.05 ug/ml	5201	0.10	
0.05 ug/ml	4275	0.10	
0.10 ug/ml	10262	0.20	
0.10 ug/ml	8377	0.20	
0.20 ug/ml	17048	0.40	
<b>Decamethrinic Acid PFB Ester</b>			
0.01 ug/ml	1383	0.01	Corr. coeff (r) 0.9905 y intercept (b) -4205 slope (m) 408918
0.01 ug/ml	4761	0.02	
0.02 ug/ml	4633	0.02	
0.05 ug/ml	18109	0.05	
0.10 ug/ml	28297	0.10	
0.10 ug/ml	36819	0.10	
0.20 ug/ml	81020	0.20	

Table 2. Recovery Data for Soil Fortified with Deltamethrin, trans-Deltamethrin and Decamethrinic Acid.


Depth	Fortification Level	Replicate	Percent Recovery	Average
<b>Deltamethrin</b>				
0-6"	0.002 ppm	1	115	100
		2	89	
		3	96	
	0.020 ppm	1	107	89
		2	84	
		3	76	
<b>trans-Deltamethrin</b>				
0-6"	0.002 ppm	1	94	81
		2	72	
		3	76	
	0.020 ppm	1	103	102
		2	103	
		3	101	
<b>Decamethrinic Acid PFB Ester</b>				
0-12"	0.010 ppm	1	104	93
		2	84	
		3	90	
	0.100 ppm	1	96	95
		2	88	
		3	101	
	0.500 ppm	1	114	106
		2	91	
		3	113	

STUDY #

8301853531

1808


Figure 1: Chemical Structures (1 of 2)



DELTAMETHRIN

RU 22974

(1R,3R)3-(2,2-dibromovinyl)-2,2-dimethylcyclopropanecarboxylic acid (S)-alpha-cyano-(3-phenoxyphenyl) methyl ester



trans-DELTAMETHRIN

RU 26979

(1R,3S)3-(2,2-dibromovinyl)-2,2-dimethylcyclopropanecarboxylic acid (S)-alpha-cyano-(3-phenoxyphenyl) methyl ester


STUDY #

930187537

157



Figure 1: Chemical Structures (2 of 2)




DECAMETHRINIC ACID

(BR<sub>2</sub>CA)

RU 23441

(1R-cis)-3-(2,2-dibromoethenyl)-2,2-dimethylcyclopropane-carboxylic acid



DECAMETHRINIC ACID

PFB ESTER

Pentafluorobenzyl ester of (1R-cis)-3-(2,2-dibromoethenyl)-2,2-dimethylcyclopropanecarboxylic acid

Figure 2: Sample Preparation Scheme (1 of 2)

40.0 grams of soil

Add 300 ml of 1:1 hexane:acetone (pesticide grade).

Shake on a horizontal shaker for 1 hour.

Suction filter through a glass fiber filter (Whatman GF/C) to remove soil. Rinse extraction bottle and filter with 1:1 hexane:acetone.

Transfer filtrate quantitatively to a 500 ml volumetric flask and bring to volume with 1:1 hexane:acetone.\*\*

cis/trans Deltamethrin  
Clean-up

Evaporate 250 ml aliquot to near dryness by rotary evaporation at 40-45 C.

Condition a silica gel SPE cartridge with 10 ml chloroform followed by 10 ml hexane.

Transfer residue with 10 ml of hexane to column

Control the drip rate through column at 1-2 ml/min.

Rinse flask with 10 ml of hexane and quantitatively transfer to column.

Rinse the column with 10 ml hexane. Discard the hexane.

Add 5.0 ml chloroform to column reservoir and collect the eluate in a graduated test tube.\*\*

\*\*Stopping point  
Refrigerate sample.

Br<sub>2</sub>CA Partition

Evaporate 250 ml aliquot to near dryness by rotary evaporation at 40-45 C.

Transfer to a 250 ml separatory funnel with 20 ml 0.1 N NaOH.

Add 20 ml ethyl ether. Shake for 1 minute with venting. Discard ethyl ether.

Add 20 ml ethyl ether. Shake for 1 minute with venting. Discard ethyl ether.

Add 10 ml of 0.5 N HCl followed by 20 ml ethyl ether.

Shake for 1 minute with venting and collect ethyl ether.

Add 20 ml ethyl ether. Shake for one minute with venting and collect ethyl ether.

Figure 2: Sample Preparation Scheme (2 of 2)

cis/trans Deltamethrin  
Cleanup (cont.)

Concentrate to dryness  
under nitrogen at 35-40 C.

Dissolve in hexane to  
final volume of 2.0 ml.

Freeze samples.

Br<sub>2</sub>CA Partition  
(cont.)

Evaporate by rotary  
evaporation at 40-45 C.

Add 100 ml acetone. Set  
up acetone control and  
0.20 ug/ml standard.\*\*

Add 1.0 ml of PFB-Br  
(5 % acetone) and  
300 ul of 30 % K<sub>2</sub>CO<sub>3</sub>.

Reflux gently for 3 hrs

Evaporate by rotary  
evaporation at 40-45 C.

Add 250 ml ether, 2 ml  
water and 10 ml ethyl  
acetate. Swirl  
vigorously.


Filter through 12 g  
anhydrous sodium sulfate  
prewashed with ether.

Evaporate to near  
dryness by rotary evap-  
oration at 40-45 C.

Bring sample to final  
volume of 10 ml with  
hexane:ether:acetone  
(98:1:1 v:v:v).


Freeze samples.

Figure 3: Representative Chromatograms of Deltamethrin and trans-Deltamethrin (1 of 2)



Typical Chromatogram  
GC Standard: 0.01 ug/ml  
Deltamethrin and trans-Deltamethrin


0.02 ng injected, GC run #590, dated 07/23/90




Typical Chromatogram  
GC Standard: 0.20 ug/ml  
Deltamethrin and trans-Deltamethrin

0.40 ng injected, GC run #586, dated 07/23/90


Figure 3:--Representative Chromatograms of Deltamethrin and trans-Deltamethrin (2 of 2)



20 mg injected, GC run #591, dated 07/23/90



20 mg injected, GC run #592, dated 07/23/90



STUDY # 2.0 mg injected, GC run #593, dated 07/23/90

Figure 4: Representative Chromatograms of Decamethrinic Acid PFB Ester (1 of 2)








Figure 4: Representative Chromatograms of Decamethrinic Acid PFB Ester (2 of 2)



STUDY # 1.0 mg injected, GC run #5453, dated 03/04/91  
Pg 0194 of 537