

## Cover Sheet for

# ENVIRONMENTAL CHEMISTRY METHOD

***Pesticide Name:*** Fipronil

***MRID #:*** 442628-28

***Matrix:*** Soil

***Analysis:*** GC/ECD

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at [flynt.elizabeth@epa.gov](mailto:flynt.elizabeth@epa.gov).

---

**APPENDIX B**

(Analytical Method)

**FIPRONIL**

Analytical Method of Analysis for Possible Residues of Fipronil and its Nonpolar Metabolites, MB 45950, MB 46136, MB 46513 and RPA 200766, in Soil

Appendix B

Rhône-Poulenc Ag Company  
Standard Operating Procedure

SOP-90231  
Issue 1.0  
Date: August 31, 1992  
Page 1 of 20

Title: FIPRONIL  
Method of Analysis for Possible Residues of Fipronil  
and its Nonpolar Metabolites, MB 45950,  
MB 46136, MB 46513 and RPA 200766, in soil.

Objective: To provide instructions in the determination of the residues  
of Fipronil (MB 46030) and its nonpolar metabolites, MB  
45950, MB 46136, MB 46513 and RPA 200766, in soil.

Scope: This method was developed specifically for the determination  
of the residues of fipronil and its nonpolar metabolites in  
soil.

Limitations: Reliable Determination Limit = 0.005 µg/g for each analyte..

By: Ahmed S. Ibrahim Sept. 1, 1992  
Ahmed S. Ibrahim, Author Date

Approval: Patricia Beo Sept. 1, '92  
Supervisor Date  
P. Wang Sept. 1, 1992  
Management Date

Distribution: All

## Appendix B

### Rhône-Poulenc Ag Company Standard Operating Procedure

SOP - 90231  
Issue 1.0  
Date: August 31, 1992  
Page 2 of 20

#### SUMMARY

Residues of fipronil and its metabolites MB 45950, MB 46513, RPA200766 and MB46136 are extracted from the soil with acetonitrile / acetone (70/30 v/v) mixture. After shaking, the soil sample is centrifuged and a portion of the extract is dried with sodium sulfate, the analytes are adsorbed onto granular charcoal, then eluted with acetonitrile. The eluant is then concentrated to the desired volume before analysis with an Electron Capture Gas Chromatograph equipped with a DB-1701 capillary column.

1. **Safety:** All procedures must be conducted in compliance with the safety regulations of Rhône Poulenc Ag Company.

#### 2. EXPERIMENTAL

##### 2.1 Reagents

Acetone, HPLC grade, Burdick & Jackson.  
Acetonitrile, HPLC-UV grade, Burdick & Jackson.  
Activated carbon, darco 20-40 mesh, Aldrich (cat # 24,226-8).  
Sodium sulfate, anhydrous, Aldrich.  
Fipronil, MB 46030, Rhône-Poulenc.  
MB 45950, Rhône-Poulenc.  
MB 46513, Rhône-Poulenc.  
MB 46136, Rhône-Poulenc.  
RPA 200766, Rhône-Poulenc.

## Appendix B

Rhône-Poulenc Ag Company  
Standard Operating Procedure

SOP - 90231

Issue 1.0

Date: August 31, 1992

Page 3 of 20

### 2.2. Chromatographic Conditions

Instrument : Gas chromatograph (HP 5890):

Column: J&W, DB-1701, 15 M x 320  $\mu$ M I.d., 0.25  $\mu$ M film thickness (part # 123-0712).

Detector: Ni<sup>63</sup> Electron Capture.

Oven Temperature: Initial temp 50 °C., hold for 1 minute.  
Program ramp 70 °C/minute to 200 °C.  
Hold for 22 minutes. Ramp 70 °C/minute to 230 °C and hold for 17 minutes

Injection Temperature: 280 °C.

Detector Temperature: 300 °C.

Flow Rates:  
Helium (carrier) 1.9 ml/min.  
95% Methane/Argon (make-up) 50 ml/min.

Integrator: Hewlett-Packard HP3396 series 2 or equivalent.  
Due to the sharp temperature rates used in this program, the integrator must be instructed to zero to the baseline when necessary.

Waters 860 Networking Computer System, version 3.0, or equivalent...

|                          |            |              |
|--------------------------|------------|--------------|
| Typical retention Times: | MB 46513   | 10.5 minutes |
|                          | MB 45950   | 15 minutes   |
|                          | MB 46030   | 16 minutes   |
|                          | RPA 200766 | 30 minutes   |
|                          | MB 46136   | 31 minutes   |

### 2.3. Apparatus

1. Mettler PM4600 balance or equivalent.
2. Mettler AE163 balance or equivalent.
3. N-Evap, model 106, Organomation Assoc. Inc. or equivalent.
4. Beckman model TJ-6 centrifuge or equivalent.
5. Adlab Shaker, Arthur H. Thomas Company, Philadelphia, PA, USA or equivalent.

Rhône-Poulenc Project No. EC-96-339

EBT Study No. RP-01-96

Page 112

## Appendix B

### Rhône-Poulenc Ag Company Standard Operating Procedure

SOP - 90231  
Issue 1.0  
Date: August 31, 1992  
Page 4 of 20

#### 2.4. Preparation of Charcoal Columns

Pack a 10 mL Bio-Rad Poly-Prep chromatography column with approximately 2 g of granular activated carbon (20-40 mesh). Although the activated carbon can be used as purchased from the manufacturer, it is recommended that the carbon be heated at approximately 120 °C over night to prevent moisture. Top with approximately 2-3 g granular anhydrous sodium sulfate. Columns similar to the 10 mL Bio-Rad chromatography columns could be used for this procedure. Open columns should be plugged with a piece of glass wool before packing with charcoal.

#### 2.5. Preparation of standard and spiking solutions

Weigh 100 mg of each of MB 46030, MB 45950, MB 46513, RPA200766 and MB46136 individually into 100 mL volumetric flasks. Dissolve each compound in approximately 50 mL of acetonitrile and bring up the volume to the mark with acetonitrile. This is now 1000 µg/mL solution of each individual standard. Transfer 1 mL of each solution to a 100 mL volumetric flask using 1 mL pipets and bring up the volume to the mark with acetonitrile. This solution is 10 µg/mL of each compound as combined standard. Using this stock solution, make the appropriate dilutions to prepare 1 µg/mL, 0.5 µg/mL, 0.1 µg/mL, 0.02 µg/mL, 0.01 µg/mL, 0.005 µg/mL and 0.002 µg/mL solutions in acetonitrile. All dilutions must be made using grade A glassware. Use these solutions as standards and soil fortification solutions.

The 1000 µg/mL solutions could be prepared as a mixture by combining the weights of all five analytes in one volumetric flask.

#### 2.6. Extraction and Cleanup Procedure:

1. Weigh 50 g soil into a 250 mL Nalgene screw-capped bottle. Fortification of control samples for the determination of recoveries should be done at this point.
2. Add 100 mL 30% acetone in acetonitrile (UV grade) and shake vigorously by hand for one minute. Place bottle in mechanical shaker and shake for 30 minutes.

## Appendix B

### Rhône-Poulenc Ag. Company Standard Operating Procedure

SOP - 90231  
Issue 1.0  
Date: August 31, 1992  
Page 5 of 20

3. Centrifuge for 5 minutes.
4. Rinse a carbon/ $\text{Na}_2\text{SO}_4$  column above with approximately 10 mL methanol, 10 mL acetone, then 10 mL of UV grade acetonitrile.
5. Place an Erlenmeyer flask under the column. Pipet a 4 mL aliquot from the soil sample and put onto the charcoal column. Rinse the column with approximately two 5 mL portions of acetonitrile.
6. Connect a 20 mL polypropylene disposable column to the top of the charcoal column. This is used as a reservoir only. Add an additional 40-45 mL UV grade acetonitrile.
7. Place the Erlenmeyer in water bath, not to exceed 50 °C and evaporate to approximately 3 mL using a gentle stream of nitrogen.
8. Transfer to a 15 mL graduated centrifuge tube, rinsing the Erlenmeyer several times with acetonitrile.
9. Place the centrifuge tube in an N-EVAP analytical evaporator (or equivalent) and evaporate to desired volume. Any necessary dilutions should be done at this point.

#### 2.7 GC Resolution and Retention Times

Typical retention times obtained using the conditions described in section 2.2 were as follows: MB 46513 (10.5 minutes), MB 45950 (15 minutes), MB 46030 (16 minutes), RPA 200766 (30 minutes), MB 46136 (31 minutes). The DB-1701 capillary column and the chromatographic conditions used gave over 1 minute separation of MB 46030 and MB 45950. It also gave about 1 minute separation of RPA 200766 and MB 46136. These separations were not possible using a variety of other capillary columns.

#### 2.8. Quantification of Residues

1. The residues of flupyrifluoride and its metabolites in the extracts should be quantified by comparison with standard solutions of appropriate concentration (standard curve). This should be done utilizing the Waters 860 Networking Computer System, version 3.0, or equivalent. The formula and an example of the calculations performed by this system are given below. In cases when this system is not operating, the integrator output could be accepted as raw

## Appendix B

### Rhône-Poulenc Ag Company Standard Operating Procedure

SOP - 90231

Issue 1.0

Date: August 31, 1992

Page 6 of 20

data. In this case the linear regression and the residues should be calculated with a scientific calculator capable of performing this operation.

2. When necessary, the extracts should be diluted so that the peak areas or heights obtained are within the linear range of the detector (0.002-0.02 µg/ml).

#### Equation used in the calculations:

$$PPB = \frac{(Y-A)}{B} \times \frac{F}{W}$$

where:

PPB = Parts per billion

Y = Peak area

A = Intercept of standard curve

B = Slope of standard curve

F = Scale factor (final volume in milliliters)

W = Sample weight

Notes: The 4 ml aliquot of the sample extract analyzed represents 2 grams of soil since 50 grams of soil are extracted with 100 ml of solvent.

#### Example of calculations:

The following example show the calculations of the concentration of fipronil (MB 46030) in soil spiked with 10 ppb of each of the five compounds under analysis.

Peak area = 27128 area counts

Intercept of the standard curve = 1822 area counts

Slope of the standard curve = 4467 counts/ng/ml

Scale factor (final volume) = 2 ml

Sample amount = 2 grams

$$PPB \text{ of fipronil} = \frac{(27128 \text{ counts} - 1822 \text{ counts})}{4467 \text{ counts/ng/ml}} \times \frac{2 \text{ ml}}{2 \text{ g}} = 9.7 \text{ ppb}$$

Rhône-Poulenc Project No. EC-96-339

EBT Study No. RP-01-96

Page 115

## Appendix B

### Rhône-Poulenc Ag Company Standard Operating Procedure

SOP - 90231  
Issue 1.0  
Date: August 31, 1992  
Page 7 of 20

The chromatograms of the sample and standards used in these calculations are shown in the appendix provided.

#### 2.9. Recovery Determinations

Two spikes should be performed with each set, a low spike (5 or 10 ppb) and a high spike (100 ppb). In addition to these routine spikes, other spikes that exceed the highest level of residue found in the samples must also be performed.

Spike recoveries should be calculated using the following formula.

$$\% \text{ recovery} = \frac{\text{ppb found}}{\text{ppb spiked}} \times 100$$

#### Example:

From the calculations given above, the concentration of fipronil in the spiked sample as calculated is 9.7 ppb. The concentration of fipronil spiked is 10 ppb, therefore,

$$\% \text{ recovery of fipronil} = \frac{9.7}{10} \times 100 = 97\%$$

#### CRITICAL POINTS

1. Any filtration steps are not recommended because of possible interferences from the filter units. If filtration becomes necessary, it should be done using glass fiber filters.
2. Reconditioning of the capillary column becomes necessary after several injections of soil extracts. Reconditioning of the column can be done by breaking about twenty centimeters from the front of the column. To avoid having to break the column a retention gap could be used.
3. The purity of acetonitrile should be checked before used in the analysis. This purity check can be done by taking approximately 40 ml of acetonitrile and concentrating it to 2 ml using nitrogen evaporator (N-Evap) and injecting the concentrated solvent on GC.

---

Appendix B

SOP-90231  
Issue 1.0  
Date: August 31, 1992  
Page 8 of 20

APPENDIX

# Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

Version: 860/V2.3 Printed: 29-May-1992 at 07:29:31  
Sda: RTEVX4 GC Project: TWR\_EC91\_158\_EC50\_200\_P28  
Status: GC Report Using Multi-method MAY\_28\_92 Line 8

SOP-90231  
Issue 1.0  
Date: August 31, 1992  
Page 9 of 20

5 S P K 28 - M A Y - 1992 12:41:58


|         |                     |            |                    |              |
|---------|---------------------|------------|--------------------|--------------|
| Header: | Acquisition method  | FIP        | Processing method  | FIP          |
|         | Units               | PFB        | System number      | 1            |
|         | Channel             | 1          | Manual injector    |              |
|         | Injection           | 1          | Total injections   | 1            |
|         | Run time            | 35.00 min  | Sample rate        | 2.00 per sec |
|         | Injection volume    | 2 uL       | Mode               | Analysis     |
|         | Acquisition version | LAC/E/V2.3 | Processing version | V2.3         |

Description:  
Node:  
Acquired on node RTEV07 system 1 for TWR\_EC91\_158\_EC50\_200\_P28

10 ppb spike

First Plot:

GCSPK Manual Injection 1 Ch 1



| GC Results: |          |       |        |     |        |
|-------------|----------|-------|--------|-----|--------|
| Peak Name   | Ret Time | Area  | Height | Int | Amount |
| MB45513     | 10.560   | 42787 | 7346   | BB  | 9.307  |
| MB45950     | 14.994   | 31453 | 3445   | BB  | 8.130  |
| MB46030     | 16.065   | 27128 | 2422   | BB  | 9.705  |
| RPA200766   | 30.146   | 13826 | 1595   | BB  | 10.430 |
| MB46136     | 30.842   | 22266 | 2399   | BB  | 6.484  |

# Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

ersion: 860/V2.3 Printed: 29-May-1992 at 07:29:21 SOP-90231  
ode: RTFVX4 GC Project: TWR\_EC91\_168\_EC50\_200\_F28 Issue 1.0  
sters GC Report Using Multi-method MAY\_28\_92 Line 7 Date: August 31, 1992  
Page 10 of 20

T C 20 - M a y - 1 9 9 2 1 1 : 5 5 : 5 9


|                     |            |                    |              |
|---------------------|------------|--------------------|--------------|
| Header:             | FIP        | Processing method  | FIP          |
| Acquisition method  | PFB        | System number      | 1            |
| Units               | 1          | Manual injector    | 1            |
| Channel             | 1          | Total injections   | 1            |
| Injection           | 35.00 min  | Sample rate        | 2.00 per sec |
| Run time            | 2 uL       | Mode               | Analysis     |
| Injection volume    | LAC/E/V2.3 | Processing version | V2.3         |
| Acquisition version |            |                    |              |

Description:  
Node:  
Acquired on node RTPL07 system 1 for TWR\_EC91\_168\_EC50\_200\_F28

## Untreated Control UTC

First Plot:

UTC Manual Injection 1 Ch 1



| GC Results: | Ret Time | Area | Height | Int | Amount |
|-------------|----------|------|--------|-----|--------|
| MB46511     | 10.500   | -    | -      | NP  | -      |
| MB45950     | 14.900   | -    | -      | NP  | -      |
| MB46030     | 16.000   | -    | -      | NP  | -      |
| RFA200766   | 30.100   | -    | -      | NP  | -      |
| MB46136     | 30.820   | -    | -      | NP  | -      |

# Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

Version: 860/V2.3 Printed: 29-May-1992 at 07:28:47  
Code: RTPVX4 GC Project: TWR\_EC91\_168\_EC50\_200\_F28  
Users GC Report Using Multi-method MAY\_28\_92 Line 3  
SOP-90231  
Issue 1.0  
Date: August 31, 1992  
Page 11 of 20

OSTD 29 - May - 1992 0 : 13 : 06


|                     |            |                    |              |
|---------------------|------------|--------------------|--------------|
| Header:             |            |                    |              |
| Acquisition method  | FIP        | Processing method  | FIP          |
| Units               | NG/ML      | System number      | 1            |
| Channel             | 1          | Manual injector    | 1            |
| Injection           | 35.00 min  | Total injections   | 1            |
| Run time            | 2 uL       | Sample rate        | 2.00 per sec |
| Injection volume    | LAC/E/V2.3 | Mode               | Calibration  |
| Acquisition version |            | Processing version | V2.3         |

Description:  
Node:  
Acquired on node RTPLO7 system 1 for TWR\_EC91\_168\_EC50\_200\_F28

20 ppb standard

First Plot:

20STD Manual Injection 1 Ch 1



| Peak Name | Ret Time | Area  | Height | Int | Amount |
|-----------|----------|-------|--------|-----|--------|
| MB45513   | 10.548   | 90355 | 15380  | BB  | CAL    |
| MB45950   | 14.987   | 75397 | 7956   | BB  | CAL    |
| MB46030   | 16.058   | 53420 | 5201   | BB  | CAL    |
| RFA200766 | 30.163   | 25467 | 2625   | BB  | CAL    |
| MB46136   | 30.859   | 47152 | 4947   | BB  | CAL    |

# Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

Version: 860/V2.3 Printed: 29-May-1992 at 07:28:54 SOP-90231  
Node: RTFVX4 GC Project: TWR\_EC91\_168\_EC50\_200\_P28 Issue 1.0  
Waters GC Report Using Multi-method MAY\_28\_92 Line 4 Date: August 31, 1992  
Page 12 of 20

10STD 28-MAY-1992 20:22:59


Reader:  
Acquisition method FIP Processing method FIP  
Units NG/ML System number 1  
Channel 1 Manual injector  
Injection 1 Total injections 1  
Run time 35.00 min Sample rate 2.00 per sec  
Injection volume 2 uL Mode Calibration  
Acquisition version LAC/E/V2.3 Processing version V2.3

Description:  
Node:  
Acquired on node RTPL07 system 1 for TWR\_EC91\_168\_EC50\_200\_P28

10 ppb standard

First Plot:

10STD Manual Injection 1 Ch 1



| GC Results: | Peak Name | Ret Time | Area  | Height | Int | Amount |
|-------------|-----------|----------|-------|--------|-----|--------|
|             | MB45930   | 14.516   | 47809 | 8072   | BB  | CAL    |
|             | MB45950   | 14.941   | 19014 | 4196   | BB  | CAL    |
|             | MB46030   | 16.010   | 27815 | 2710   | BB  | CAL    |
|             | RPA200765 | 30.131   | 12860 | 1449   | BB  | CAL    |
|             | MB46136   | 30.823   | 28829 | 2982   | BB  | CAL    |

# Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

Revision: 860/V2.3 Printed: 29-May-1992 at 07:29:10 SOP-90231  
Code: RTFVX4 GC Project: TWR\_EC91\_168\_EC50\_200\_P28 Issue 1.0  
Users GC Report Using Multi-method NMY\_28\_92 Line 6 Date: August 31, 1992  
Page 13 of 20

STD 28-May-1992 14:13:57


|                     |            |                    |              |
|---------------------|------------|--------------------|--------------|
| Header:             |            |                    |              |
| Acquisition method  | FIP        | Processing method  | FIP          |
| Units               | NG/ML      | System number      | 1            |
| Channel             | 1          | Manual injector    |              |
| Injection           | 1          | Total injections   | 1            |
| Run time            | 35.00 min  | Sample rate        | 2.00 per sec |
| Injection volume    | 2 uL       | Mode               | Calibration  |
| Acquisition version | LAC/E/V2.3 | Processing version | V2.3         |

## Description:

Node:  
Acquired on node RTPL07 system 1, for TWR\_EC91\_168\_EC50\_200\_P28

First Plot: 8 ppb standard

STD Manual Injection 1 Ch 1



| GC Results: |          |       |        |     |        |
|-------------|----------|-------|--------|-----|--------|
| Peak Name   | Ret Time | Area  | Height | Int | Amount |
| MB45950     | 10.575   | 36369 | 6212   | BB  | CAL    |
| MB46030     | 15.007   | 31594 | 3373   | BB  | CAL    |
| MB46030     | 16.074   | 24529 | 2418   | BB  | CAL    |
| RPA200766   | 30.149   | 12073 | 1322   | BB  | CAL    |
| MB46136     | 30.843   | 20448 | 2180   | BB  | CAL    |

# Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

Version: 860/V2.3 Printed: 29-May-1992 at 07:28:38  
Mode: RTEVX4 GC Project: TWR\_EC91\_168 EC50\_200\_P28  
Waters GC Report Using Multi-methd MAY\_28\_92 Line 2

SOP-90231  
Issue 1.0  
Date: August 31, 1992  
Page 14 of 20

5 S T D

28 - M a y - 1 9 9 2

1 1 : 1 9 : 1 8

|                     |            |                    |              |
|---------------------|------------|--------------------|--------------|
| Header:             |            |                    |              |
| Acquisition method  | FIP        | Processing method  | FIP          |
| Units               | MG/ML      | System number      | 1            |
| Channel             | 1          | Manual injector    | 1            |
| Injection           | 1          | Total injections   | 1            |
| Run time            | 35.00 min  | Sample rate        | 2.00 per sec |
| Injection volume    | 2 uL       | Mode               | Calibration  |
| Acquisition version | LAC/E/V2.3 | Processing version | V2.3         |


## Description:

Node:  
Acquired on node RTPL07 system 1 for TWR\_EC91\_168 EC50\_200\_P28

5 ppb standard

## First Plot:

SSTB Manual Injection 1 Ch 1



| GC Results: | Ret Time | Area  | Height | Int | Amount |
|-------------|----------|-------|--------|-----|--------|
| Peak Name   |          |       |        |     |        |
| MB46513     | 10.744   | 22449 | 3789   | BB  | CAL    |
| MB45950     | 15.203   | 19411 | 2086   | BB  | CAL    |
| MB46030     | 16.274   | 14719 | 1457   | BB  | CAL    |
| RPA200766   | 30.243   | 7347  | 824    | BB  | CAL    |
| MB46136     | 30.942   | 13256 | 1452   | BB  | CAL    |

## Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

```

Version: 860/V1.3      Printed: 29-May-1992 at 07:29:01
Mode: RTFVX4          GC Project: TWR EC91_168 EC50_200 #28
Interns GC Report Using Multi-method MAY_28_92 Line 5

```

SOP-90231  
Issue 1.0  
Date: August 31, 1992  
Page 15 of 20

**2 3 T D**

28 - Мау - 1992

17 : 18 : 45

|                     |            |                    |              |
|---------------------|------------|--------------------|--------------|
| Header:             |            |                    |              |
| Acquisition method  | FIP        | Processing method  | FIP          |
| Units               | NG/ML      | System number      | 1            |
| Channel             | 1          | Manual injector    |              |
| Injection           | 1          | Total injections   | 1            |
| Run time            | 35.00 min  | Sample rate        | 2.00 per sec |
| Injection volume    | 2 uL       | Mode               | Calibration  |
| Acquisition version | LAC/2/V2.3 | Processing version | V2.3         |


**Description:**

Node:  
Acquired on node RTPL07 system 1 for TWR\_EC91\_168\_EC50\_200\_P28

**2 ppb standard**

**First Plot:**

25TD Manual Injection 1 Ch 1



| GC Results:<br>Peak Name | Ret Time | Area  | Height | Int | Amount |
|--------------------------|----------|-------|--------|-----|--------|
| MB46513                  | 10.511   | 10059 | 1729   | BB  | CAL    |
| MB45950                  | 14.937   | 8099  | 891    | BB  | CAL    |
| MB46030                  | 15.008   | 5971  | 617    | BB  | CAL    |
| RPA200766                | 30.126   | 2477  | 296    | BB  | CAL    |
| MB46136                  | 30.824   | 7441  | 580    | BB  | CAL    |

Rhone-Poulenc Project No. EC-96-339  
EBT Study No. RP-01-96  
Page 124

# Appendix B

Rhône-Poulenc Ag Company  
2 Alexander Drive, Research Triangle Park, NC 27709

Version: 860/V1.3 Printed: 29-May-1992 at 07:31:31  
Code: RTFVX4 GC Project: TWR\_EC91 168 EC50 200 P28  
Users GC Report Using Multi-method MAY\_28\_92 Line 20

Page 1  
User: TROBINSO

Users GC Calibration Method FIF

SOP-90231  
Issue 1.0  
Date: August 31, 1992  
Page 16 of 20

Channel 1

Channel Parameters:  
Calibrate Channel ? YES  
Retention Time Offset 0 sec  
Relative Peak Window 3%  
Absolute Peak Window 0 sec  
Update Retention Time ? NO

| Peak Name | Ret Time | Fit Type | Levels | By   | R | I | D | Group |
|-----------|----------|----------|--------|------|---|---|---|-------|
| MB46513   | 10.500   | Linear   | 5      | Area |   |   |   | 1     |
| MB45950   | 14.900   | Linear   | 5      | Area |   |   |   | 1     |
| MB46030   | 16.000   | Linear   | 5      | Area |   |   |   | 1     |
| RPA200766 | 30.100   | Linear   | 5      | Area |   |   |   | 1     |
| MB46136   | 30.820   | Linear   | 5      | Area |   |   |   | 1     |

| Peak Name | Coefficients                                                                          |
|-----------|---------------------------------------------------------------------------------------|
| MB46513   | A = 9.495617e+02<br>B = 4.495404e+03<br>R <sup>2</sup> = 0.998647                     |
| MB45950   | A = 1.075074e+03<br>B = 3.736436e+03<br>R <sup>2</sup> = 0.999539                     |
| MB46030   | A = 1.822199e+03 ← Intercept<br>B = 2.607622e+03 ← Slope<br>R <sup>2</sup> = 0.996198 |
| RPA200766 | A = 8.332415e+02<br>B = 1.245729e+03<br>R <sup>2</sup> = 0.990874                     |
| MB46136   | A = 3.225849e+03<br>B = 2.244372e+03<br>R <sup>2</sup> = 0.986522                     |

Level Table:  
Peak Name: MB46513 Retention Time: 10.500 min

| Level | Amount    | Response     | Reps | % RSD |
|-------|-----------|--------------|------|-------|
| 1     | 2.000000  | 10059.000000 | 1    | 0.00  |
| 2     | 5.000000  | 24449.000000 | 1    | 0.00  |
| 3     | 8.000000  | 26369.000000 | 1    | 0.00  |
| 4     | 10.000000 | 47803.000000 | 1    | 0.00  |
| 5     | 20.000000 | 90355.000000 | 1    | 0.00  |

# Appendix B

Version: 860/V2.5  
 Code: ATPVX4 GC Project: TWR EC91 168 EC50 200 P28  
 Users GC Report Using Multi-method MAY 28 92 Line 20

User: TROBINSO

Channel 1 (continued)

SOP-90231  
 Issue 1.0  
 Date: August 31, 1992  
 Page 17 of 20

Peak Name: MB45950 Retention Time: 14.900 min

| Level | Amount    | Response     | Reps | % RSD |
|-------|-----------|--------------|------|-------|
| 1     | 2.000000  | 8099.000000  | 1    | 0.00  |
| 2     | 5.000000  | 19411.000000 | 1    | 0.00  |
| 3     | 8.000000  | 31594.000000 | 1    | 0.00  |
| 4     | 10.000000 | 39014.000000 | 1    | 0.00  |
| 5     | 20.000000 | 75397.000000 | 1    | 0.00  |

Peak Name: MB46030 Retention Time: 15.000 min

| Level | Amount    | Response     | Reps | % RSD |
|-------|-----------|--------------|------|-------|
| 1     | 2.000000  | 5971.000000  | 1    | 0.00  |
| 2     | 5.000000  | 14719.000000 | 1    | 0.00  |
| 3     | 8.000000  | 24529.000000 | 1    | 0.00  |
| 4     | 10.000000 | 27815.000000 | 1    | 0.00  |
| 5     | 20.000000 | 53420.000000 | 1    | 0.00  |

Peak Name: RPA200766 Retention Time: 30.100 min

| Level | Amount    | Response     | Reps | % RSD |
|-------|-----------|--------------|------|-------|
| 1     | 2.000000  | 2477.000000  | 1    | 0.00  |
| 2     | 5.000000  | 7347.000000  | 1    | 0.00  |
| 3     | 8.000000  | 12073.000000 | 1    | 0.00  |
| 4     | 10.000000 | 12860.000000 | 1    | 0.00  |
| 5     | 20.000000 | 25467.000000 | 1    | 0.00  |

Peak Name: MB46136 Retention Time: 30.820 min

| Level | Amount    | Response     | Reps | % RSD |
|-------|-----------|--------------|------|-------|
| 1     | 2.000000  | 7441.000000  | 1    | 0.00  |
| 2     | 5.000000  | 13256.000000 | 1    | 0.00  |
| 3     | 8.000000  | 20448.000000 | 1    | 0.00  |
| 4     | 10.000000 | 28829.000000 | 1    | 0.00  |
| 5     | 20.000000 | 47152.000000 | 1    | 0.00  |

## Appendix B

Version: 860/V2.3  
 Model: RTFVX4  
 GC Project: TWR EC91\_168 EC50\_200\_P28  
 Users GC Report Using Multi-method MAY\_28\_92 Line 20

User: TRÖBINSO

SOP-90231

Issue 1.0


Date: August 31, 1992

Page 18 of 20


Channel 1 (continued)

Curves:

GC Calibration Method FIP Peak M046513 Channel 1



GC Calibration Method FIP Peak M043928 Channel 1




# Appendix B


ersion: 860/V2.3  
 Code: RTFV4  
 GC Project: TWR EC91 168 EC50 200 P28  
 ators GC Report Using Multi-method MAY\_28\_92 Line 20  
 Channel 1 (continued)

SOP-90231  
 Issue 1.0  
 Date: August 31, 1992  
 Page 19 of 20

GC Calibration Method FIP Peak PB46830 Channel 1



GC Calibration Method FIP Peak RP9280766 Channel 1



## Appendix B


```

Version: 860/V2.3      Printout: 23-May-1992 at 07:11:11      Page: 1
Code: RTFVX4      GC Project: TWR EGS1 168 EGSU 200 P28      User: TROBINSO
Notes: GC Report Using Multi-method MAY 28 92 Line 20

Channel 1 (continued)

GC Calibration Method FIP Peak M046136 Channel 1

```



---

**APPENDIX C**

(Control Soil)