

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Dichloropropene

MRID #: 445365-11

Matrix: Water

Analysis: GC/MS

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.


GRM: 94.11
EFFECTIVE: April 14, 1995
SUPERSEDES: New

Determination of Residues of 1,2-Dichloropropane, *cis*- and *trans*-1,3-Dichloropropene, and Trichloronitromethane in Water by Purge and Trap Extraction, Capillary Gas Chromatography and Mass Selective Detection


C. E. Kubitschek, S. C. Dolder and H. E. Dixon-White
North American Environmental Chemistry Laboratory
DowElanco
Indianapolis, Indiana 46268-1053

A. Scope


This method is applicable for the quantitative determination of residues of 1,2-dichloropropane (1,2-D), *cis*- and *trans*-1,3-dichloropropene (1,3-D), and trichloronitromethane (TCNM) in water over the concentration range of 0.05-40 ng/mL with a validated limit of quantitation of 0.05 ng/mL for each compound.




1,2-D
CAS No. 78-87-5



cis-1,3-D
CAS No. 10061-01-5



trans-1,3-D
CAS No. 10061-02-6



TCNM
CAS No. 76-06-2

B. Principle

This analytical method is based on established EPA purge and trap methodology for volatile organic analytes (VOAs) such as Method 524.2 (1). The volatile chlorinated hydrocarbons are purged from the water sample by sparging with helium, after which they are captured on a sorbent-containing trap. When purging is complete, the trap is heated and backflushed with helium, and the VOAs are desorbed and transferred to the injection port of a gas chromatograph (GC). The analytes are then separated on a capillary column and quantitated using mass selective detection (MSD).

Effective Date: April 14, 1995

GRM 94.11

C. Safety Precautions

1. Each analyst must be acquainted with the potential hazards of the reagents, products, and solvents used in this method before commencing laboratory work. SOURCES OF INFORMATION INCLUDE: MATERIAL SAFETY DATA SHEETS, LITERATURE, AND OTHER RELATED DATA. Safety information on non-DowElanco products should be obtained from the container label or from the supplier. Disposal of reagents, reactants, and solvents must be in compliance with local, state, and federal laws and regulations.
2. The analytes covered by this method, particularly TCNM, can cause severe and possibly fatal respiratory distress at air concentrations in the ppm range. All operations involving the neat analytes, or concentrated solutions of these compounds, must be carried out in a fume hood. In addition, the effluent from the GC split vent should be routed into an exhaust vent or through a carbon trap to prevent the release of the analytes into the laboratory air.

D. Equipment (Note N.1.)

1. Balance, analytical, Model AE200, Mettler Instrument Corporation, Hightstown, NJ 08520.
2. Gas chromatograph, Model 5890 Series II, Hewlett-Packard, Wilmington, DE 19808.
3. Mass selective detector, Model 5971, Hewlett-Packard, Palo Alto, CA 94304.
4. Mass selective detector data system, Model G1034B, Hewlett-Packard.
5. Purge and trap autosampler, Model 2016, Tekmar Company, Cincinnati, OH 45249.
6. Purge and trap concentrator, Model 3000, Tekmar Company.
7. Water purification system, Model Milli-Q UV Plus, Millipore Corporation, Milford, MA 01757.

E. Glassware and Materials (Note N.1.)

1. Column, capillary gas chromatography, DB-VRX, 30 m x 0.25 mm i.d., 1.4 μ m film thickness, catalog number 122-1534, J&W Scientific, Folsom, CA 95630.
2. Column inlet liner, deactivated, double gooseneck, catalog number 5181-3315, Hewlett-Packard, Kennett Square, PA 19348.
3. Filter, charcoal, catalog number 7972, Chrompack, Inc., Raritan, NJ 08869. (Note N.2.)
4. Filter, moisture, catalog number 7971, Chrompack, Inc. (Note N.2.)
5. Filter, oxygen, catalog number 7970, Chrompack, Inc. (Note N.2.)
6. Flask, 40 mL volumetric, catalog number 042029-0801, Kontes, Vineland, NJ 08360.
7. Frit sparge glassware, 25 mL, catalog number 14-3022000, Tekmar Company.
8. Gas, helium, 99.995% purity, Airco, Murray Hill, NJ 07974.
9. Syringe, gas-tight, fixed needle, 10, 100, 500, and 1000 μ L, catalog numbers 1701, 1710, 1750, and 1001, Hamilton Company, Reno, NV 89520.
10. Syringe, gas-tight, Luer Lock, 25 mL, catalog number 1025, Hamilton Company.

Effective Date: April 14, 1995

GRM 94.11

11. Syringe valve, catalog number 2-0940M, Supelco, Inc., Bellefonte, PA 16823.
12. Trap, Tenax, catalog number 12-0083-003, Tekmar Company.

F. Reagents and Chemicals (Note N.1.)

1. Reagents

- a. Hydrochloric acid, concentrated, ACS reagent grade, catalog number A144-500, Fisher Scientific, Pittsburgh, PA 15219.
- b. Internal standard, 2-bromo-1-chloropropane, 95%, compound number 23,127-4, Aldrich Chemical Company, Milwaukee, WI 53233.
- c. Methyl alcohol, purge and trap grade, catalog number 41,481-6, Sigma-Aldrich, St. Louis, MO 63178.
- d. Standards
 - (1) 1,2-dichloropropane
The 1,2-D standard used for generating the validation data contained in this method was Lot Number AGR277102, with a purity of 99.2%.
 - (2) *cis*-1,3-dichloropropene
The *cis*-1,3-D standard used for generating the validation data contained in this method was Lot Number AGR164301, with a purity of 97.1%.
 - (3) *trans*-1,3-dichloropropene
The *trans*-1,3-D standard used for generating the validation data contained in this method was Lot Number TSN100232, with a purity of 97.2%.
 - (4) trichloronitromethane
The TCNM standard used for generating the validation data contained in this method was Lot Number TSN100245, with a purity of 98.9%.
Obtain all standards from Test Substance Coordinator, DowElanco, Indianapolis, IN 46268-1053.
- e. Water, distilled/deionized, purified using a Milli-Q UV Plus purification system (Section D.7.).

2. Prepared Solutions

Hydrochloric Acid, 1/1 (v/v):

Prepare by adding 100 mL concentrated HCl to 100 mL distilled/deionized water.

G. Preparation of Standards

1. Preparation of VOA Standard Solutions

NOTE: CARRY OUT THIS PROCEDURE IN A FUME HOOD. When mixing standard solutions, do not shake them excessively, as loss of analytes may occur. Mix by capping and gently inverting the solutions approximately five times. Store all standard solutions under frozen conditions in vials with PTFE-lined lids and screw-cap closures. Do not store samples and standards together in the same freezer. Allow standard solutions to warm to room temperature prior to use.

Effective Date: April 14, 1995

GRM 94.11

- a. Tare a 100-mL volumetric flask containing approximately 80 mL of methanol. Add 0.1100-0.1300 g of 1,2-D dropwise via Pasteur pipette (approximately eight to ten drops). Make sure that the liquid falls directly into the methanol and does not run down the inside walls of the volumetric flask. Stopper the flask. Reweigh the flask and calculate the exact weight of the analyte added (Note N.3.). Bring the solution to volume with methanol to yield a stock solution of approximately 1.2 g/L 1,2-D. Calculate the exact concentration of the stock solution; if the standard is less than 97% pure, make the correction for percent purity as follows (A purity of 96% for the 1,2-D standard was assumed for the purpose of this example):

$$\frac{\text{1,2-D, g}}{\text{Total solution volume, L}} \times \frac{\% \text{ purity}}{100} = \text{corrected concentration (g/L)}$$

$$\frac{0.1100 \text{ g 1,2-D}}{0.100 \text{ L}} \times \frac{96}{100} = 1.056 \text{ g/L}$$

- b. Repeat step G.1.a. for *cis*-1,3-D, *trans*-1,3-D, and TCNM, preparing a stock solution of approximately 1.2 g/L for each analyte.
- c. Using gas-tight syringes, transfer 1/C mL of each stock solution (where C = the concentration of the stock solution in g/L) into a single 100-mL volumetric flask containing approximately 80 mL of methanol. Dilute to volume with methanol to obtain a standard solution containing all four analytes, each at a concentration of 10.0 mg/L (or 10000 ng/mL). For example, a stock solution of 1.056 g/L 1,2-D would be diluted as follows:

$$1/1.056 = 0.947 \text{ mL} = 947 \mu\text{L}$$

$$(947 \mu\text{L})(1 \times 10^{-6} \text{ L}/\mu\text{L})(1.056 \text{ g/L})(1/100 \text{ mL}) = 0.00001 \text{ g/mL} = 10000 \text{ ng/mL}$$

- d. Transfer 10.0 mL of the 10000 ng/mL standard solution into a 100-mL volumetric flask. Dilute to volume with methanol to obtain a standard solution containing all four analytes each at a concentration of 1000 ng/mL.

2. Preparation of the Internal Standard

- a. Following the procedure outlined above in Sections G.1.a. and c., prepare a 10 mg/L solution of the internal standard, 2-bromo-1-chloropropane. (It is not necessary to correct for percent purity of the internal standard.)
- b. Transfer 10.0 mL of the 10 mg/L internal standard solution to a 40-mL volumetric flask and dilute to volume with methanol to obtain a 2.50 mg/L internal standard solution.

3. Preparation of the Calibration Standard Solutions

Prepare aqueous calibration standards of the VOAs by transferring aliquots of the standard solutions (Sections G.1.-2.) into the appropriate volumetric flasks containing distilled/deionized water, as outlined below. Use an appropriately-sized gas-tight syringe to transfer the standard solution and inject it into the water below the narrow neck of the flask. Because aqueous solutions of VOAs are not stable in any container with headspace, each calibration standard should be prepared just prior to purge and trap analysis.

Effective Date: April 14, 1995

GRM 94.11

a. Low range calibration standard solutions (Section J.1.):

Aliquot of 1000 ng/mL Standard Soln. ^a	Final Calibration Std. Volume	Calibration Std. Final Conc.	Equivalent Sample Conc. ^b
μL	mL	ng/mL	ng/mL
5.00	200	0.0250	0.0250
5.00	100	0.0500	0.0500
10.0	100	0.100	0.100
25.0	100	0.250	0.250
50.0	100	0.500	0.500
100.0	100	1.00	1.00

^a Section G.1.d.

^b The concentration in an aqueous calibration standard directly corresponds to the same concentration in a field sample, unless the sample has been diluted prior to analysis.

Transfer a 5.0- μL aliquot of the 10.0 mg/L internal standard solution (Section G.2.a.) to each calibration standard (with the exception of the 0.025 ng/mL standard, which requires 10.0 μL of internal standard as a result of its final volume) to obtain an internal standard concentration of 0.500 ng/mL.

b. High range calibration standard solutions (Section J.1.):

Aliquot of 10000 ng/mL Standard Soln. ^a	Final Calibration Std. Volume	Calibration Std. Final Conc.	Equivalent Sample Conc. ^b
μL	mL	ng/mL	ng/mL
10.0	100	1.00	1.00
20.0	100	2.00	2.00
40.0	100	4.00	4.00
100.0	100	10.0	10.0
200.0	100	20.0	20.0
400.0	100	40.0	40.0

^a Section G.1.c.

^b The concentration in an aqueous calibration standard directly corresponds to the same concentration in a field sample unless the sample has been diluted prior to analysis.

Transfer a 50.0- μL aliquot of the 10.0 mg/L internal standard solution (Section G.2.a.) to each calibration standard to obtain an internal standard concentration of 5.00 ng/mL.

H. Purge and Trap

1. Purge and Trap Concentrator

Install the Tenax trap (Section E.12.) and sparge glassware (Section E.7.) on the Tekmar purge and trap concentrator (Section D.6) and autosampler (Section D.5.) following the manufacturer's recommended procedure.

Effective Date: April 14, 1995

GRM 94.11

2. Typical Operating Conditions

Instrumentation:	Tekmar 3000 Concentrator Tekmar 2016 Autosampler
Purge/Carrier Gas:	helium
Purge Flow	40 mL/min
Desorb Flow	20 mL/min
Trap Pressure	5 psi
Tekmar Controller Menu Settings:	
3000 Transfer Line	140 °C
3000 Valve	140 °C
2016 Transfer Line	130 °C
2016 Valve	130 °C
Moisture control system temp.	130 °C
Trap purge ready temp.	30 °C
Trap purge temp. setting	20 °C (i.e., not heated above ambient temp. during purge)
Sample heater	Off
Prepurge	0.00 min
Preheat time	0.00 min
Purge time	11.00 min
Dry purge time	2.00 min
Moisture control system desorb temp.	50 °C
GC Start	Desorb start
Cryo Focuser	Not applicable
GC cycle time	23 min
Desorb preheat	170 °C
Desorb time	2.00 min
Desorb temp.	180 °C
Sample drain	off
Bake time	10 min
Bake temp.	185 °C
Bake gas bypass	On
Bake gas bypass delay	2.00 min
Moisture control system bake	180 °C

I. Gas Chromatography/Mass Spectrometry**1. Column (Note N.4.)**

Install the column inlet liner (Section E.2.) and the capillary column (Section E.1.) in the split/splitless injection port of the gas chromatograph/mass spectrometer (GC/MSD) following the manufacturer's recommended procedure.

Effective Date: April 14, 1995

GRM 94.11

2. Typical Operating Conditions

Instrumentation: Hewlett-Packard Model 5890 Series II GC
 Hewlett-Packard Model 5971 Mass Selective Detector
 Hewlett-Packard Model G1034B Data System Software

Column: J&W Scientific fused silica capillary
 DB-VRX liquid phase
 30 m x 0.25 mm i.d.
 1.4 µm film thickness

Temperatures:
 Column 35 °C for 1.0 min
 35 °C to 140 °C at 9 °C/min
 140 °C for 0.10 min
 140 °C to 210 °C at 20 °C/min
 210 °C for 2.0 min

Injector 200 °C
 Interface 230 °C

Carrier Gas: helium

Head Pressure 6 psi
 Linear Velocity approximately 30 cm/sec @ 35 °C

Injection Mode: Split (Note N.4.)

Splitter Flow 20 mL/min
 Septum Purge Off (capped)

Detector: electron impact selected ion monitoring

Calibration Program midmass autotune (Note N.5.)
 Electron Multiplier 1600 volts

Ions Monitored:

Compound	<i>m/z</i> , Quantitation	<i>m/z</i> , Confirmation
1,2-D	63	76
<i>cis</i> -1,3-D	75	112
<i>trans</i> -1,3-D	75	112
TCNM	119	82
2-bromo-1-chloropropane	77	—

Dwell Time: 75 msec

Full scan mass spectra of the above analytes and the internal standard are shown in Figures 1-5.

3. Typical Chromatograms

Typical chromatograms of a standard, control sample, and a 0.0500 ng/mL recovery sample for water are illustrated in Figures 6-17.

Effective Date: April 14, 1995

GRM 94.11

J. Calibration

1. General Approach

Because of the wide range of standard concentrations (0.0250-40.0 ng/mL) used in this method, it is unlikely that any single calibration curve will provide accurate quantitation over the entire range. In general, deviations from the calculated curve will be most severe at the low end of the curve, which affects measurements near the LOQ. To improve quantitation, it is useful to divide the calibration range into two subranges and to produce a separate calibration curve for each. Thus, the low range will encompass VOA concentrations of 0.0250-1.00 ng/mL, and the high range 1.00-40.0 ng/mL.

A calibration check (Section J.4.) must be carried out at the beginning of each twelve-hour period during which samples are analyzed in order to confirm calibration of the instrumentation.

2. Initial Calibration

- a. Prior to analyzing any standards, bake out the trap at approximately 185 °C and set the GC oven temperature to 220 °C for 10 minutes.
- b. Select either the low or high range for calibration. After the system has been returned to its starting conditions, load the first aqueous calibration standard (Section G.3.) in the selected range by removing the plunger from the 25-mL gas-tight Luer Lock syringe (Section E.10.), attaching the closed syringe valve (Section E.11.) to the Luer lock, and pouring the standard into the open end of the syringe until it is nearly full. (Do not draw the standard up into the syringe.) Replace the syringe plunger in the barrel, open the valve, expel any air from the syringe, and adjust the volume to exactly 25.0 mL. Open the loading valve on the purge and trap device and load the standard into the sparge tube. Analyze the standard by purge and trap using GC/MSD as described in Sections H. and I.
- c. Repeat J.2.b. for each standard in the calibration range. After each standard has been analyzed, rinse the sparge tube with two approximately 25-mL aliquots of distilled/deionized water before loading the next standard. Standards should always be analyzed from lowest to highest concentration to minimize carryover. (Note N.6.)

3. Calibration Curve

- a. Following analysis of the range of calibration standards described in G.3.a. (low range) or G.3.b. (high range), determine the peak areas for 1,2-D (*m/z* 63), *cis*-1,3-D (*m/z* 75), *trans*-1,3-D (*m/z* 75), TCNM (*m/z* 119), and 2-bromo-1-chloropropane (IS) (*m/z* 77).
- b. For each analyte, prepare a standard curve by plotting the concentration (ng/mL) on the abscissa (x-axis) and the standard/IS peak area ratio (Quantitation Ratio) on the ordinate (y-axis) as shown in Figures 18-21. Using regression analysis, determine the equation for the curve with respect to the abscissa.

For example, using power regression (2) with the 1,2-D data from Figure 18:

$$Y = \text{constant} \times X^{\text{(exponent)}}$$

Effective Date: April 14, 1995

GRM 94.11

Solving for X:

$$X = \left(\frac{Y}{\text{constant}} \right)^{1/\text{exponent}}$$

where: X = concentration (ng/mL), Y = Quantitation Ratio

$$\text{1,2-D Conc.} = \left(\frac{\text{Quantitation Ratio}}{\text{constant}} \right)^{(1/\text{exponent})}$$

ng/mL

$$\text{1,2-D Conc.} = \left(\frac{\text{Quantitation Ratio}}{1.37594} \right)^{(1/1.04500)}$$

ng/mL

- c. Typical calibration curves for the determination of each of the analytes in water are shown in Figures 18-21.

4. Calibration Check

- a. For the low calibration range, prepare a 0.100 ng/mL calibration standard (Section G.3.a.); for the high range, prepare a 2.00 ng/mL calibration standard (Section G.3.b.).
- b. Analyze the calibration check standard and calculate the amount of each analyte present, using the equation for the current calibration curve for each analyte. (Section J.3.)
- c. Calculate the absolute value of the percent difference in the calculated and the theoretical standard concentration as follows:

$$|(1 - (\text{Calculated}/\text{Theoretical}))| \times 100$$

If the difference in the calculated and theoretical value for each analyte in the calibration check standard is less than 10%, the existing curves are considered to be valid. Sample analysis may proceed.

- d. If the difference in the calculated and theoretical value for any analyte exceeds 10%, the original calibration is no longer considered valid, and a new initial calibration (Section J.2.) must be carried out prior to continuing with sample analysis. A new initial calibration must always be carried out for all analytes prior to sample analysis whenever major maintenance (column change, source cleaning, filament or multiplier replacement, trap replacement, etc.) is performed on the purge and trap GC/MSD system. Minor maintenance (replacing GC septum or inlet liner, removing a portion of the upper end of the column, etc.) does not automatically necessitate recalibration, but must be examined on a case-by-case basis.

Effective Date: April 14, 1995

GRM 94.11

K. Determination of Recovery of 1,2-D, cis-1,3-D, trans-1,3-D, and TCNM in Water

1. Preparation of Recovery Samples

- a. Following analysis of the last calibration standard or check standard, rinse the sample sparge tube(s) twice with distilled/deionized water and analyze a reagent blank consisting of 25.0 mL distilled/deionized water. Pour the water into the open barrel of the Luer Lock syringe with attached syringe valve, and adjust the volume to 25.0 mL (Section J.2.b.). Remove the syringe valve, pull the plunger back slightly, and inject 5.0 μ L of the 2.5 mg/L internal standard (Section G.2.b.) through the Luer Lock tip into the sample, if the analyses are being conducted in the low-level calibration range. If the analyses are being conducted in the high-level calibration range, add 50.0 μ L of the 2.5 mg/L internal standard to the sample in the same manner. Load the sample into the sparge tube and begin the analysis. Calculate the levels of analytes detected (Section J.3.). If the concentration of any analyte exceeds 30% of the targeted limit of quantitation for the analysis (Section M.1.c.), additional blanks should be run until a "clean" blank is obtained.
- b. Obtain a control sample of ground or surface water from the field sampling location. Analyze an aliquot of the control (Section K.1.a.) to demonstrate that none of the analytes of interest are detectable in the control sample.
- c. Prepare fortified samples at appropriate concentrations over the range of 0.050-40.0 ng/mL as described in Section G.3. for the standards, using the natural ground or surface water instead of distilled/deionized water.

It may be necessary to preserve natural water samples with 1/1 HCL (Section F.2.) (two drops per mL) if bacteria are present that degrade the specific analytes of interest. (Note N.7.)
- d. Analyze the fortified recovery samples by purge and trap using GC/MSD as described in Sections H. and I.

2. Calculation of Percent Recovery

- a. Determine the peak areas for 1,2-D (m/z 63, 76), cis-1,3-D (m/z 75, 112), trans-1,3-D (m/z 75, 112), TCNM (m/z 119, 82), and IS (m/z 77) for each calibration standard analyzed as part of the current calibration curve (Section J.2.).
- b. For each standard, calculate the confirmation ratio for each of the four analytes. The average confirmation ratio for each analyte will be used to confirm the presence of that analyte in the water samples.

For example, using the data for 1,2-D from Figure 6:

$$\text{Confirmation Ratio} = \frac{\text{peak area of quantitation ion}}{\text{peak area confirmation ion}}$$

$$\text{Confirmation Ratio} = \frac{\text{peak area } m/z \text{ 63}}{\text{peak area } m/z \text{ 76}}$$

$$\text{Confirmation Ratio} = \frac{1415}{511}$$

$$\text{Confirmation Ratio} = 2.77$$

Effective Date: April 14, 1995

GRM 94.11

Positive confirmation of the presence of each analyte is indicated when the confirmation ratio for the samples is in the range of $\pm 15\%$ of the average found for the standards in the current calibration curve.

- c. Prepare a standard curve for each of the analytes and determine the equations for the calibration curves using regression analysis as described in J.3.
- d. Determine the net concentration of 1,2-D in each recovery sample by first subtracting the average quantitation ratio of the control samples from that of the recovery sample. Substitute the net quantitation ratio obtained into the calibration curve equation and solve for the concentration.

For example, using power regression and the 1,2-D data from Figures 10, 14, and 18:

$$\text{1,2-D Conc. (ng/mL)} = \left(\frac{(\text{net quantitation ratio})}{1.37594} \right)^{1/1.0450}$$

$$\text{1,2-D Conc. (ng/mL)} = \left(\frac{(0.0673)}{1.37594} \right)^{1/1.0450}$$

$$\text{1,2-D Conc.} = 0.0557 \text{ ng/mL}$$

- e. Determine the percent recovery by dividing the concentration found for each recovery sample by the theoretical concentration added.

$$\text{Recovery} = \frac{\text{Concentration Found}}{\text{Concentration Added}} \times 100\%$$

$$\text{Recovery} = \frac{0.0557 \text{ ng/mL}}{0.0500 \text{ ng/mL}} \times 100$$

$$\text{Recovery} = 111\%$$

- f. For each analyte, determine the concentration found and corresponding percent recovery as described for 1,2-D in Section K.2.d.-e.

The average recovery for each analyte in a given sample set will be used to correct for daily method efficiency.

L. Determination of VOAs in Water

1. Calibrate the instrument/check calibration over the appropriate calibration range for the sample set as described in Section J.
2. Prepare and analyze reagent blanks, controls, and recovery samples as described in Section K.1. Rinse the sample sparge tube(s) between samples with distilled/deionized water. If standards or recovery samples spiked at 1.00 ng/mL or higher were analyzed prior to initiating the analysis of a set of water samples, analyze a blank in that sparge tube as in section K.1.a. in order to ensure that carryover is not a problem.

Effective Date: April 14, 1995

GRM 94.11

3. Allow the field samples to warm to room temperature. Examine the sample vials for headspace; any sample contained in a vial which is not completely full (free of air bubbles) is an invalid sample and should be discarded.
4. If it is anticipated that the concentration is likely to exceed the upper limit of the calibration range, the water sample should be diluted with distilled/deionized water by an appropriate factor prior to analysis. If the concentration of analytes in a group of samples is unknown, a 1:100 dilution of at least one of the samples should be analyzed first to estimate the sample concentration and prevent contaminating the instrument with a very high-level sample. (Note N.8.)
5. Load a sample into a gas-tight syringe as in Section K.1.a. After the sample volume has been adjusted to exactly 25.0 mL, transfer 5 µL (for the low calibration range) or 50.0 µL (for the high calibration range) of the 2.5 mg/L internal standard solution into the sample using a gas-tight syringe. Add the internal standard solution by removing the syringe valve, pulling the plunger back slightly, and injecting the internal standard through the Luer-Lock tip.
6. Open the loading valve on the purge and trap device and transfer the sample into the sparge tube. Analyze the sample by purge and trap using GC/MSD as described in Sections H. and I. (The purge and trap autosampler permits the loading and subsequent automated analysis of a set of 16 samples.)
7. Determine the concentration of each volatile analyte of interest by substituting the quantitation ratio into the equation for the corresponding standard calibration curve (Section J.3.), calculating the uncorrected (gross) result.

For example, using the 1,2-D data from Figures 14 and 18, the uncorrected concentration is calculated as follows:

$$\text{1,2-D Conc. (ng/mL)} = \left(\frac{(\text{1,2-D quantitation ratio})}{1.37594} \right)^{1/1.0450} \times \text{any dilution factor}$$

$$\text{1,2-D Conc. (ng/mL)} = \left(\frac{0.0673}{1.37594} \right)^{1/1.0450} \times \text{any dilution factor}$$

$$\text{1,2-D Conc. (ng/mL)} = 0.0557 \times 1$$

$$\text{1,2-D Conc.} = 0.0557 \text{ ng/mL}$$

8. The uncorrected results as determined in section L.7. must be reported. In addition, the results can be corrected for method recovery using the following procedure:
 - a. Calculate the mean % recovery of each analyte for the recovery samples analyzed with the treated samples on the same day.
 - b. Determine the corrected analyte concentration in the water samples as follows (Assume a mean % recovery for 1,2-D of 105 for the purpose of this example.):

$$\text{1,2-D Conc. (corrected ng/mL)} = \frac{\text{uncorrected 1,2-D (ng/mL)} \times 100}{\% \text{ Recovery}}$$

Effective Date: April 14, 1995

GRM 94.11

$$\text{1,2-D Conc. (corrected ng/mL)} = \frac{0.0557 \times 100}{105} = 0.0530$$

9. When analyzing samples using a high range calibration, any sample giving a peak area ratio response for one or more of the analytes which is greater than 10 percent below that of the lowest standard in the current calibration curve must be reanalyzed using the low level calibration range (Section G.3.a.). Any sample giving a peak area ratio response greater than 10 percent above that of the highest standard in the current calibration curve must be reanalyzed, following appropriate dilution of the sample with distilled/deionized water to bring it within the calibration range. (Note N.8.) The additional dilution must then be accounted for in the calculation of the result (Section L.7.).

M. Results and Discussion

1. Method Validation

a. Recovery Levels and Precision

A method validation study was conducted to determine the recovery levels and the precision of the method for 1,2-D, *cis*-1,3-D, *trans*-1,3-D, and TCNM in water; the results are summarized in Tables I-IV, respectively.

Recovery values of 1,2-D from water samples fortified over the concentration range 0.0500 to 40.0 ng/mL averaged 103% with one standard deviation equal to 5% (Table I).

Recovery values of *cis*-1,3-D from water samples fortified over the concentration range 0.0500 to 40.0 ng/mL averaged 102% with one standard deviation equal to 4% (Table II).

Recovery values of *trans*-1,3-D from water samples fortified over the concentration range 0.0500 to 40.0 ng/mL averaged 100% with one standard deviation equal to 6% (Table III).

Recovery values of TCNM from water samples fortified over the concentration range 0.0500 to 40.0 ng/mL averaged 104% with one standard deviation equal to 6% (Table IV).

b. Standard Curve Linearity

The average coefficient of determination (r^2) for the power regression equations describing the detector response as a function of the standard calibration curve concentration were greater than 0.99 for all four analytes.

c. Calculated Limits of Quantitation and Detection

Following established guidelines (3), the limits of quantitation (LOQ) and detection (LOD) were calculated using the standard deviation from the 0.0500 ng/mL recovery results for each of the four analytes. The LOQ was calculated as ten times the standard deviation ($10s$), and the LOD was calculated as three times the standard deviation ($3s$) of the results of the analysis of eight samples. The results are summarized in Tables V-VIII.

The calculated statistics support an LOQ between 0.0049 and 0.024 ng/mL for the four analytes, lower than the targeted method LOQ of 0.0500 ng/mL; however, results should not be quantified at levels below which no recovery samples have been analyzed.

Effective Date: April 14, 1995

GRM 94.11

2. Confirmation of Residue Identity

Confirmation of the presence of residues is described in Section K.2.b. For each of the four analytes, confirmation is by comparison of the retention time (GC) as well as the confirmation ratios resulting from selected ion monitoring (mass spectrometry). When the peak retention time in a sample matches that of the standards, positive confirmation of the presence of the analyte is indicated when the confirmation ratio for the sample is in the range of $\pm 15\%$ of the average found for the corresponding standards. If additional confirmation is required beyond that discussed in this method, the mass spectrum of each of the four analytes contain additional ions that may be used for confirmation.

3. Assay Time

A typical analytical run would consist of a minimum of six calibration standards encompassing the expected range of sample concentrations, a reagent blank, a control (a non-fortified sample), a minimum of two fortified controls (one of which must be at the target LOQ), and sixteen samples. This typical analytical run could be prepared in approximately 8 hours, with the purge and trap/chromatographic analysis continuing into the same evening.

As indicated in Section J., once the instrument has been calibrated, a calibration check can be done every twelve hours by analyzing a single standard to verify that the existing calibration curve is still valid. This reduces the amount of time required for daily calibration, and allows the analysis of a greater number of samples per day.

N. Notes

1. Equipment, glassware, materials, reagents, and chemicals considered to be equivalent to those specified may be substituted with the understanding that their performance must be confirmed by appropriate tests. Common laboratory supplies are assumed to be readily available and are, therefore, not listed.
2. The filters are used in the carrier gas supply lines to purify the helium entering the gas chromatograph and sparging apparatus.
3. An alternate method for transferring the neat standards to the methanol contained in the volumetric flask is to tare a gas-tight syringe, draw an amount of the neat standard into the syringe, reweigh the syringe determining the weight of the standard, and transfer the standard from the syringe to the flask by injecting it below the surface of the methanol contained in the flask. Rinse the syringe into the flask to be certain that the analyte contained in the syringe is completely transferred to the flask.
4. The purge and trap concentrator is interfaced with the gas chromatograph by cutting the GC helium inlet line approximately 4 cm from its entry into the injection port. The cut end of the helium line is routed to the concentrator, and the concentrator's transfer line is connected to the remaining short section of tubing leading to the injector. The flow from the purge and trap concentrator is 20 mL/minute and the flow through the column is approximately 1 mL/minute; therefore, a 19:1 split will occur. The injector should be run in the split mode at all times. While this technique appears to sacrifice sensitivity by splitting off 95% of the analytes, it permits the use of normal-dimension capillary columns (0.18 or 0.25 mm i.d.) rather than megabore columns which necessitate the use of jet separators or cryogenic interfaces at the GC/MSD. By virtue of the superior resolution possible with the 0.25 mm i.d. capillary column, the desired limits of quantitation (0.0500 ng/mL) are readily attained for all analytes.

Effective Date: April 14, 1995

GRM 94.11

5. The mass spectrometer should be tuned prior to analysis of a set of calibration standards and generation of a calibration curve. Once the instrument has been calibrated it should not be tuned again until just prior to analysis of a fresh set of calibration standards for generation of a new calibration curve.
6. It is possible to load all standards in a calibration curve onto the Tekmar 2016 autosampler prior to beginning the analysis; however, this then requires rinsing of all sparge tubes used and confirming the absence of carryover prior to analyzing samples. Using a single sparge tube to analyze consecutively loaded standards confines the potential for contamination and need for cleaning to a single tube.
7. Over the course of developing this method, it was discovered that levels of TCNM in distilled/deionized water were stable at room temperature for several hours; however, equal levels in natural surface water decreased rapidly. This was attributed to bacterial degradation and preservation by acidification of the surface water samples to approximately pH 2 stopped the loss of the analyte. (No degradation was observed for the other three analytes in surface water, either acidified or unacidified.) Surface water from any test site should be evaluated to determine the need for preservation prior to collection of the field samples. If necessary, acidification must be done at the time of sample collection for field samples, and at the time of fortification for calibration standards and recovery samples.
8. Following removal of an aliquot for analysis from a sample vial, the remaining sample is no longer valid due to the presence of headspace in the vial. As a result, any time a sample must be reanalyzed, the aliquot for reanalysis must be taken from a sealed replicate of the initial sample.

Effective Date: April 14, 1995

GRM 94.11

O. References

1. Eichelberger, J. W.; Munch, J. W.; Bellar, T. A., "Measurement of Purgeable Organic Compounds in Water by Capillary Column Gas Chromatography/Mass Spectrometry", Method 524.2, 1992, Environmental Monitoring Systems Laboratory, Office of Research and Development, U.S. Environmental Protection Agency.
2. Stolz, W. L., DowElanco Laboratory Notebook, B065, p 23, 1993, unpublished data of DowElanco.
3. Stolz, W. L., DowElanco Laboratory Notebook, B154, p 14, 1993, unpublished data of DowElanco.
4. Stolz, W. L., DowElanco Laboratory Notebook, B154, p 30, 1993, unpublished data of DowElanco.
5. Stolz, W. L., Report FA&PC935247, 1993, unpublished data of DowElanco.
6. *HP-41C/41CV Standard Applications Handbook*, Hewlett-Packard Publication No. 00041-90402, 1982, pp 42-48.
7. Keith, L. H.; Crummett, W. B.; Deegan, J.; Libby, R. A.; Taylor, J. T.; Wentler, G., "Principles of Environmental Analysis", *Anal. Chem.*, 1983, 55, pp 2210-2218.

The information herein is presented in good faith, but no warranty, expressed or implied, is given nor is freedom from any patent owned by DowElanco or by others to be inferred. In the hands of qualified personnel, the procedures are expected to yield results of sufficient accuracy for their intended purposes, but recipients are cautioned to confirm the reliability of their techniques, equipment, and standards by appropriate tests. Anyone wishing to reproduce or publish the material in whole or in part should request written permission from DowElanco.

Effective Date: April 14, 1995

GRM 94.11

Table I. Recovery of 1,2-Dichloropropane from Water

Sample Number	Date of Analysis	1,2-D, ng/mL		Percent Recovery	
		Added	Found		
14422301	06-May-1994	0.000	0.000	--	
14422301	11-May-1994	0.000	0.000	--	
14422301	13-May-1994	0.000	0.000	--	
14422301	06-May-1994	0.0500	0.0521	104	
14422301	06-May-1994	0.0500	0.0535	107	
14422301	06-May-1994	0.0500	0.0499	100	
14422301	06-May-1994	0.0500	0.0557	111	
14422301	06-May-1994	0.0500	0.0554	111	
14422301	06-May-1994	0.0500	0.0516	103	
14422301	06-May-1994	0.0500	0.0542	108	
14422301	06-May-1994	0.0500	0.0494	99	
14422301	13-May-1994	0.250	0.252	101	
14422301	13-May-1994	0.250	0.249	99	
14422301	13-May-1994	1.00	1.08	108	
14422301	13-May-1994	1.00	1.07	107	
14422301	11-May-1994	5.00	4.81	96	
14422301	11-May-1994	5.00	5.14	103	
14422301	11-May-1994	10.0	10.1	101	
14422301	11-May-1994	10.0	10.2	102	
14422301	11-May-1994	40.0	38.4	96	
14422301	11-May-1994	40.0	39.6	99	
				\bar{x} =	103
				s =	5
				n =	18

Effective Date: April 14, 1995

GRM 94.11

Table II. Recovery of *cis*-1,3-Dichloropropene from Water

Sample Number	Date of Analysis	<i>cis</i> -1,3-D, ng/mL		Percent Recovery	
		Added	Found		
14422301	06-May-1994	0.000	0.000	--	
14422301	11-May-1994	0.000	0.000	--	
14422301	13-May-1994	0.000	0.000	--	
14422301	06-May-1994	0.0500	0.0502	100	
14422301	06-May-1994	0.0500	0.0518	104	
14422301	06-May-1994	0.0500	0.0514	103	
14422301	06-May-1994	0.0500	0.0513	103	
14422301	06-May-1994	0.0500	0.0512	102	
14422301	06-May-1994	0.0500	0.0514	103	
14422301	06-May-1994	0.0500	0.0515	103	
14422301	06-May-1994	0.0500	0.0508	102	
14422301	13-May-1994	0.250	0.251	100	
14422301	13-May-1994	0.250	0.251	100	
14422301	13-May-1994	1.00	1.11	111	
14422301	13-May-1994	1.00	1.10	110	
14422301	11-May-1994	5.00	4.65	93	
14422301	11-May-1994	5.00	4.92	98	
14422301	11-May-1994	10.0	9.90	99	
14422301	11-May-1994	10.0	10.1	101	
14422301	11-May-1994	40.0	39.1	98	
14422301	11-May-1994	40.0	39.4	98	
				\bar{x} =	102
				<i>s</i> =	4
				<i>n</i> =	18

Effective Date: April 14, 1995

GRM 94.11

Table III. Recovery of *trans*-1,3-Dichloropropene from Water

Sample Number	Date of Analysis	<i>trans</i> -1,3-D, ng/mL		Percent Recovery	
		Added	Found		
14422301	06-May-1994	0.000	0.000	—	
14422301	11-May-1994	0.000	0.000	—	
14422301	13-May-1994	0.000	0.000	—	
14422301	06-May-1994	0.0500	0.0495	99	
14422301	06-May-1994	0.0500	0.0474	95	
14422301	06-May-1994	0.0500	0.0490	98	
14422301	06-May-1994	0.0500	0.0504	101	
14422301	06-May-1994	0.0500	0.0470	94	
14422301	06-May-1994	0.0500	0.0472	94	
14422301	06-May-1994	0.0500	0.0483	97	
14422301	06-May-1994	0.0500	0.0485	97	
14422301	13-May-1994	0.250	0.255	102	
14422301	13-May-1994	0.250	0.257	103	
14422301	13-May-1994	1.00	1.12	112	
14422301	13-May-1994	1.00	1.13	113	
14422301	11-May-1994	5.00	4.60	92	
14422301	11-May-1994	5.00	4.86	97	
14422301	11-May-1994	10.0	9.92	99	
14422301	11-May-1994	10.0	10.4	104	
14422301	11-May-1994	40.0	39.5	99	
14422301	11-May-1994	40.0	40.0	100	
				\bar{x} =	100
				s =	6
				n =	18

Effective Date: April 14, 1995

GRM 94.11

Table IV. Recovery of Trichloronitromethane from Water

Sample Number	Date of Analysis	TCNM, ng/mL		Percent Recovery	
		Added	Found		
14422301	06-May-1994	0.000	0.000	--	
14422301	11-May-1994	0.000	0.000	--	
14422301	13-May-1994	0.000	0.000	--	
14422301	06-May-1994	0.0500	0.0527	105	
14422301	06-May-1994	0.0500	0.0505	101	
14422301	06-May-1994	0.0500	0.0564	113	
14422301	06-May-1994	0.0500	0.0561	112	
14422301	06-May-1994	0.0500	0.0531	106	
14422301	06-May-1994	0.0500	0.0532	106	
14422301	06-May-1994	0.0500	0.0547	109	
14422301	06-May-1994	0.0500	0.0553	111	
14422301	13-May-1994	0.250	0.240	96	
14422301	13-May-1994	0.250	0.247	99	
14422301	13-May-1994	1.00	1.08	108	
14422301	13-May-1994	1.00	1.08	108	
14422301	11-May-1994	5.00	4.76	95	
14422301	11-May-1994	5.00	4.97	99	
14422301	11-May-1994	10.0	10.2	102	
14422301	11-May-1994	10.0	10.3	103	
14422301	11-May-1994	40.0	38.2	96	
14422301	11-May-1994	40.0	38.7	97	
				\bar{x} =	104
				<i>s</i> =	6
				<i>n</i> =	18

Effective Date: April 14, 1995

GRM 94.11

Table V. Calculated Limits of Detection and Quantitation for the Determination of 1,2-Dichloropropane in Water

Sample Number	Date of Analysis	1,2-D, ng/mL	
		Added	Found
14422301	06-May-1994	0.0500	0.0521
14422301	06-May-1994	0.0500	0.0535
14422301	06-May-1994	0.0500	0.0499
14422301	06-May-1994	0.0500	0.0557
14422301	06-May-1994	0.0500	0.0554
14422301	06-May-1994	0.0500	0.0516
14422301	06-May-1994	0.0500	0.0542
14422301	06-May-1994	0.0500	0.0494
		\bar{x} =	0.053
		s =	0.0024
		LOD ^a (3s) =	0.0072
		LOQ ^b (10s) =	0.024

^a LOD = Limit of Detection.

^b LOQ = Limit of Quantitation.

Effective Date: April 14, 1995

GRM 94.11

Table VI. Calculated Limits of Detection and Quantitation for the Determination of *cis*-1,3-Dichloropropene in Water

Sample Number	Date of Analysis	<i>cis</i> -1,3-D, ng/mL	
		Added	Found
14422301	06-May-1994	0.0500	0.0502
14422301	06-May-1994	0.0500	0.0518
14422301	06-May-1994	0.0500	0.0514
14422301	06-May-1994	0.0500	0.0513
14422301	06-May-1994	0.0500	0.0512
14422301	06-May-1994	0.0500	0.0514
14422301	06-May-1994	0.0500	0.0515
14422301	06-May-1994	0.0500	0.0508
		\bar{x} =	0.051
		s =	0.00049
		LOD ^a (3s) =	0.0015
		LOQ ^b (10s) =	0.0049

^a LOD = Limit of Detection.

^b LOQ = Limit of Quantitation.

Effective Date: April 14, 1995

GRM 94.11

Table VII. Calculated Limits of Detection and Quantitation for the Determination of *trans*-1,3-Dichloropropene in Water

Sample Number	Date of Analysis	<i>trans</i> -1,3-D, ng/mL	
		Added	Found
14422301	06-May-1994	0.0500	0.0495
14422301	06-May-1994	0.0500	0.0474
14422301	06-May-1994	0.0500	0.0490
14422301	06-May-1994	0.0500	0.0504
14422301	06-May-1994	0.0500	0.0470
14422301	06-May-1994	0.0500	0.0472
14422301	06-May-1994	0.0500	0.0483
14422301	06-May-1994	0.0500	0.0485
		\bar{x} =	0.048
		s =	0.0012
		LOD ^a (3s) =	0.0036
		LOQ ^b (10s) =	0.012

^a LOD = Limit of Detection.

^b LOQ = Limit of Quantitation.

Effective Date: April 14, 1995

GRM 94.11

Table VIII. Calculated Limits of Detection and Quantitation for the Determination of Trichloronitromethane in Water



Sample Number	Date of Analysis	TCNM, ng/mL	
		Added	Found
14422301	06-May-1994	0.0500	0.0527
14422301	06-May-1994	0.0500	0.0505
14422301	06-May-1994	0.0500	0.0564
14422301	06-May-1994	0.0500	0.0561
14422301	06-May-1994	0.0500	0.0531
14422301	06-May-1994	0.0500	0.0532
14422301	06-May-1994	0.0500	0.0547
14422301	06-May-1994	0.0500	0.0553
		\bar{x} =	0.054
		s =	0.0020
		LOD ^a (3s) =	0.0060
		LOQ ^b (10s) =	0.020

^a LOD = Limit of Detection.

^b LOQ = Limit of Quantitation.

Effective Date: April 14, 1995

GRM 94.11





1,2-D
Formula: C₃Cl₂H₆
Molecular Weight: 112

Figure 1. Mass Spectrum of 1,2-Dichloropropane

Effective Date: April 14, 1995

GRM 94.11





cis-1,3-D
Formula: C₃Cl₂H₄
Molecular Weight: 110

Figure 2. Mass Spectrum of *cis*-1,3-Dichloropropene

Effective Date: April 14, 1995

GRM 94.11





trans-1,3-D
Formula: C₃Cl₂H₄
Molecular Weight: 110

Figure 3. Mass Spectrum of *trans*-1,3-Dichloropropene

Effective Date: April 14, 1995

GRM 94.11





TCNM
Formula: CCl_3NO_2
Molecular Weight: 163

Figure 4. Mass Spectrum of Trichloronitromethane

Effective Date: April 14, 1995

GRM 94.11




2-Bromo-1-chloropropane (internal standard)
Formula: C_3ClBrH_6
Molecular Weight: 156

Figure 5. Mass Spectrum of 2-Bromo-1-chloropropane

Effective Date: April 14, 1995

GRM 94.11



Data File : STD050.D
ALS Bottle : 1
Date : 2 May 94 9:29 am
Integration: 11/30/94
Data Path : C:\CHEMPC\DATA\CK0502\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: TELONE VOLATILES, 50 NG/L, & ISTD, 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHKEK


Internal Standard Retention Time	:	9.35
Peak Area (2-Bromo-1-Chloropropane)	:	23118
1,2-D Retention Time	:	7.42
Peak Area (M/Z 63)	:	1415
Peak Area (M/Z 76)	:	511
1,2-D Confirmation	:	
Ratio of M/Z 63/76	:	2.7691
1,2-D Quantitation	:	
Ratio of M/Z 63/ISTD	:	0.0612

1,2-D Concentration: 0.0500 ng/mL
Average 1,2-D Standard Confirmation Ratio: 2.76

Figure 6. Typical Chromatogram of a 0.0500 ng/mL 1,2-D Standard Equivalent to 0.0500 ng/mL 1,2-D in Water

Effective Date: April 14, 1995

GRM 94.11



Data File : STD050.D
ALS Bottle : 1
Date : 2 May 94 9:29 am
Integration: 11/30/94
Data Path : C:\CHEMPC\DATA\CK0502\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: TELONE VOLATILES, 50 NG/L, & ISTD, 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHK


Internal Standard Retention Time	:	9.35
Peak Area (2-Bromo-1-Chloropropane)	:	23118
<i>cis</i> -1,3-D Retention Time	:	8.33
Peak Area (M/Z 75)	:	1972
Peak Area (M/Z 112)	:	424
<i>cis</i> -1,3-D Confirmation	:	
Ratio of M/Z 75/112	:	4.6509
<i>cis</i> -1,3-D Quantitation	:	
Ratio of M/Z 75/ISTD	:	0.0853

cis-1,3-D Concentration: 0.0500 ng/mL
Average *cis*-1,3-D Standard Confirmation Ratio: 5.16

Figure 7. Typical Chromatogram of a 0.0500 ng/mL *cis*-1,3-D Standard Equivalent to 0.0500 ng/mL *cis*-1,3-D in Water

Effective Date: April 14, 1995

GRM 94.11



Data File : STD050.D
ALS Bottle : 1
Date : 2 May 94 9:29 am
Integration: 11/30/94
Data Path : C:\CHEMPC\DATA\CK0502\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: TELONE VOLATILES, 50 NG/L, & ISTD, 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHKEK


Internal Standard Retention Time	:	9.35
Peak Area (2-Bromo-1-Chloropropane)	:	23118
trans-1,3-D Retention Time	:	8.89
Peak Area (M/Z 75)	:	1200
Peak Area (M/Z 112)	:	301
trans-1,3-D Confirmation	:	
Ratio of M/Z 75/112	:	3.9867
trans-1,3-D Quantitation	:	
Ratio of M/Z 75/ISTD	:	0.0519

trans-1,3-D Concentration: 0.0500 ng/mL
Average trans-1,3-D Standard Confirmation Ratio: 4.33

Figure 8. Typical Chromatogram of a 0.0500 ng/mL trans-1,3-D Standard Equivalent to 0.0500 ng/mL trans-1,3-D in Water

Effective Date: April 14, 1995

GRM 94.11



Data File : STD050.D
ALS Bottle : 1
Date : 2 May 94 9:29 am
Integration: 11/30/94
Data Path : C:\CHEMPC\DATA\CK0502\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: TELONE VOLATILES, 50 NG/L, & ISTD, 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHKEK


Internal Standard Retention Time	:	9.35
Peak Area (2-Bromo-1-Chloropropane):	:	23118
Chloropicrin Retention Time	:	9.56
Peak Area (M/Z 119):	:	699
Peak Area (M/Z 82):	:	187
Chloropicrin Confirmation	:	
Ratio of M/Z 119/82:	:	3.7380
Chloropicrin Quantitation	:	
Ratio of M/Z 119/ISTD:	:	0.0302

TCNM Concentration: 0.0500 ng/mL
Average TCNM Standard Confirmation Ratio: 4.40

Figure 9. Typical Chromatogram of a 0.0500 ng/mL TCNM Standard Equivalent to 0.0500 ng/mL TCNM in Water

Effective Date: April 14, 1995

GRM 94.11



Data File : CTRL0506.D
ALS Bottle : 1
Date : 6 May 94 4:01 pm
Integration: 11/30/94 (llx)
Data Path : C:\CHEMPC\DATA\CK0506\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: IMMOKALEE 14422301 WATER W/ ISTD AT 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHK

Internal Standard Retention Time : 9.37
Peak Area (2-Bromo-1-Chloropropane): 19909


No 1,2-D Found

1,2-D Concentration: 0.0 ng/mL
Average 1,2-D Standard Confirmation Ratio: 2.76

Figure 10. Typical Chromatogram of a Control Water Sample for the Determination of 1,2-D

Effective Date: April 14, 1995

GRM 94.11



Data File : CTRL0506.D
ALS Bottle : 1
Date : 6 May 94 4:01 pm
Integration: 11/30/94 (HLL)
Data Path : C:\CHEMPC\DATA\CK0506\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: IMMOKALEE 14422301 WATER W/ ISTD AT 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHK

Internal Standard Retention Time : 9.37
Peak Area (2-Bromo-1-Chloropropane): 19909


No cis-1,3-D Found

cis-1,3-D Concentration: 0.0 ng/mL
Average cis-1,3-D Standard Confirmation Ratio: 5.16

Figure 11. Typical Chromatogram of a Control Water Sample for the Determination of cis-1,3-D

Effective Date: April 14, 1995

GRM 94.11



Data File : CTRL0506.D
ALS Bottle : 1
Date : 6 May 94 4:01 pm
Integration: 11/30/94 *ckk*
Data Path : C:\CHEMPC\DATA\CK0506\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: IMMOKALEE 14422301 WATER W/ ISTD AT 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHK

Internal Standard Retention Time : 9.37
Peak Area (2-Bromo-1-Chloropropane): 19909


No trans-1,3-D Found

trans-1,3-D Concentration: 0.0 ng/mL
Average *trans*-1,3-D Standard Confirmation Ratio: 4.33

Figure 12. Typical Chromatogram of a Control Water Sample for the Determination of *trans*-1,3-D

Effective Date: April 14, 1995

GRM 94.11



Data File : CTRL0506.D
ALS Bottle : 1
Date : 6 May 94 4:01 pm
Integration: 11/30/94 *allie*
Data Path : C:\CHEMPC\DATA\CK0506\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: IMMOKALEE 14422301 WATER W/ ISTD AT 500 NG/L
Sample Info: 25.0 ML PURGED
Operator : CRAIG KUBITSCHKEK


Internal Standard Retention Time : 9.37
Peak Area (2-Bromo-1-Chloropropane): 19909
No Chloropicrin Found

TCNM Concentration: 0.0 ng/mL
Average TCNM Standard Confirmation Ratio: 4.40

Figure 13. Typical Chromatogram of a Control Water Sample for the Determination of TCNM

Effective Date: April 14, 1995

GRM 94.11



Data File : 0101004.D
 ALS Bottle : 1
 Date : 6 May 94 6:11 pm
 Integration: 11/30/94
 Data Path : C:\CHEMPC\DATA\CK0506\
 Instrument : GC/MSD - GC S/N 3118A35302
 Sample Name: IMMOKALEE RECOVERY SAMPLE 4, 50 NG/L VOAS
 Sample Info: SAMPLE 14422301, 25.0 ML PURGED
 Operator : CRAIG KUBITSCHK


Internal Standard Retention Time	:	9.37
Peak Area (2-Bromo-1-Chloropropane):	:	22728
1,2-D Retention Time	:	7.44
Peak Area (M/Z 63):	:	1530
Peak Area (M/Z 76):	:	498
1,2-D Confirmation		
Ratio of M/Z 63/76:	:	3.0723
1,2-D Quantitation		
Ratio of M/Z 63/ISTD:	:	0.0673

1,2-D Concentration: 0.0557 ng/mL
 Average 1,2-D Standard Confirmation Ratio: 2.76
 1,2-D Recovery: 111%

Figure 14. Typical Chromatogram of a Control Water Sample Fortified with 0.0500 ng/mL 1,2-D

Effective Date: April 14, 1995

GRM 94.11



Data File : 0101004.D
ALS Bottle : 1
Date : 6 May 94 6:11 pm
Integration: 11/30/94
Data Path : C:\CHEMPC\DATA\CK0506\
Instrument : GC/MSD - GC S/N 3118A35302
Sample Name: IMMOKALEE RECOVERY SAMPLE 4, 50 NG/L VOAS
Sample Info: SAMPLE 14422301, 25.0 ML PURGED
Operator : CRAIG KUBITSCHK


Internal Standard Retention Time	:	9.37
Peak Area (2-Bromo-1-Chloropropane):	:	22728
cis-1,3-D Retention Time	:	8.35
Peak Area (M/Z 75):	:	2038
Peak Area (M/Z 112):	:	419
cis-1,3-D Confirmation	:	
Ratio of M/Z 75/112:	:	4.8640
cis-1,3-D Quantitation	:	
Ratio of M/Z 75/ISTD:	:	0.0897

cis-1,3-D Concentration: 0.0513 ng/mL
Average cis-1,3-D Standard Confirmation Ratio: 5.16
cis-1,3-D Recovery: 103%

Figure 15. Typical Chromatogram of a Control Water Sample Fortified with 0.0500 ng/mL cis-1,3-D

Effective Date: April 14, 1995

GRM 94.11



Data File : 0101004.D
 ALS Bottle : 1
 Date : 6 May 94 6/11 pm
 Integration: 11/30/94
 Data Path : C:\CHEMPC\DATA\CK0506\
 Instrument : GC/MSD - GC S/N 3118A35302
 Sample Name: IMMOKALEE RECOVERY SAMPLE 4, 50 NG/L VOAS
 Sample Info: SAMPLE 14422301, 25.0 ML PURGED
 Operator : CRAIG KUBITSCHK


Internal Standard Retention Time	:	9.37
Peak Area (2-Bromo-1-Chloropropane)	:	22728
trans-1,3-D Retention Time	:	8.92
Peak Area (M/Z 75)	:	1252
Peak Area (M/Z 112)	:	309
trans-1,3-D Confirmation		
Ratio of M/Z 75/112	:	4.0518
trans-1,3-D Quantitation		
Ratio of M/Z 75/ISTD	:	0.0551

trans-1,3-D Concentration: 0.0504 ng/mL
 Average trans-1,3-D Standard Confirmation Ratio: 4.33
 trans-1,3-D Recovery: 101%

Figure 16. Typical Chromatogram of a Control Water Sample Fortified with 0.0500 ng/mL trans-1,3-D

Effective Date: April 14, 1995

GRM 94.11



Data File : 0101004.D
 ALS Bottle : 1
 Date : 6 May 94 6:11 pm.
 Integration: 11/30/94 *HOW*
 Data Path : C:\CHEMPC\DATA\CK0506\
 Instrument : GC/MSD - GC S/N 3118A35302
 Sample Name: IMMOKALEE RECOVERY SAMPLE 4, 50 NG/L VOAS
 Sample Info: SAMPLE 14422301, 25.0 ML PURGED
 Operator : CRAIG KUBITSCHK


Internal Standard Retention Time	:	9.37
Peak Area (2-Bromo-1-Chloropropane):	:	22728
Chloropicrin Retention Time	:	9.58
Peak Area (M/Z 119):	:	799
Peak Area (M/Z 82):	:	181
Chloropicrin Confirmation	:	
Ratio of M/Z 119/82:	:	4.4144
Chloropicrin Quantitation	:	
Ratio of M/Z 119/ISTD:	:	0.0352

TCNM Concentration: 0.0561 ng/mL
 Average TCNM Standard Confirmation Ratio: 4.40
 TCNM Recovery: 112%

Figure 17. Typical Chromatogram of a Control Water Sample Fortified with 0.0500 ng/mL TCNM

Effective Date: April 14, 1995

GRM 94.11




1,2-D Conc. ng/mL	1,2-D Quantitation Ratio
0.0250	0.0285
0.0500	0.0612
0.100	0.124
0.250	0.328
0.500	0.676
1.00	1.34

Power Regression Equation: $X = (Y/1.37594)^{1/1.0450}$
X = ng/mL, Y = Quantitation Ratio
Coefficient of Determination (r²): 0.9998

Figure 18. Typical Calibration Curve (Low Range) for the Determination of 1,2-Dichloropropane in Water Samples

Effective Date: April 14, 1995

GRM 94.11




<i>cis</i> -1,3-D Conc. ng/mL	<i>cis</i> -1,3-D Quantitation Ratio
0.0250	0.0446
0.0500	0.0853
0.100	0.171
0.250	0.451
0.500	0.949
1.00	1.88

Power Regression Equation: $X = (Y/1.87793)^{1/1.02433}$
 X = ng/mL, Y = Quantitation Ratio
 Coefficient of Determination (r^2): 0.9996

Figure 19. Typical Calibration Curve (Low Range) for the Determination of *cis*-1,3-D in Water Samples

Effective Date: April 14, 1995

GRM 94.11



<i>trans</i> -1,3-D Conc. ng/mL	<i>trans</i> -1,3-D Quantitation Ratio
0.0250	0.0302
0.0500	0.0519
0.100	0.101
0.250	0.264
0.500	0.573
1.00	1.13


Power Regression Equation: $X = (Y/1.09966)^{1/1.00183}$
X = ng/mL, Y = Quantitation Ratio

Coefficient of Determination (r^2): 0.9976

Figure 20. Typical Calibration Curve (Low Range) for the Determination of *trans*-1,3-Dichloropropene in Water Samples

Effective Date: April 14, 1995

GRM 94.11



TCNM Conc. ng/mL	TCNM Quantitation Ratio
0.0250	0.0145
0.0500	0.0302
0.100	0.0674
0.250	0.182
0.500	0.415
1.00	0.884

Power Regression Equation: $X = (Y/0.88292)^{(1/1.11905)}$
 X = ng/mL, Y = Quantitation Ratio
 Coefficient of Determination (r^2): 0.9998

Figure 21. Typical Calibration Curve (Low Range) for the Determination of Trichloronitromethane in Water Samples