

Cover Sheet for

## **ENVIRONMENTAL CHEMISTRY METHOD**

***Pesticide Name:*** Pymetrozine Degrad (GS23199 & CGA294849)

***MRID #:*** 444113-35

***Matrix:*** Soil

***Analysis:*** HPLC/UV

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at [flynt.elizabeth@epa.gov](mailto:flynt.elizabeth@epa.gov).

VOLUME 43 OF 58 OF SUBMISSION

ANALYTICAL METHOD

CGA-215944

METHOD TITLE

Analytical Method for the Determination of GS-23199 and  
CGA-294849, Metabolites of CGA-215944, in Soil by High Performance  
Liquid Chromatography with UV Detection Including Validation Data

DATA REQUIREMENT

40 CFR 158, Subdivision N, 164-1

AUTHOR

John D. Vargo, Ph.D.

METHOD COMPLETION DATE

March 7, 1997

PERFORMING LABORATORY

Novartis Crop Protection, Inc.  
(formerly Ciba Crop Protection)  
Post Office Box 18300  
Greensboro, North Carolina 27419

PROJECT ID


Ciba Method No. AG-666  
Ciba Study No. 376-96

SUBMITTER/SPONSOR

Novartis Crop Protection, Inc.  
(formerly Ciba Crop Protection)  
Post Office Box 18300  
Greensboro, NC 27419

VOLUME 1 OF 1 OF STUDY

Page 1 of 84



STATEMENT OF NO DATA CONFIDENTIALITY CLAIM

No claim of confidentiality is made for any information contained in this study on the basis of its falling within the scope of FIFRA Section 10 (d) (1) (A), (B) or (C).

COMPANY: Novartis Crop Protection, Inc.

COMPANY AGENT: Richard Pence

TITLE: Senior Regulatory Manager


SIGNATURE:

Richard Pence

DATE:

10/10/97

These data are the property of Novartis Crop Protection, Inc. and, as such, are considered to be confidential for all purposes other than compliance with FIFRA Section 10. Submission of these data in compliance with FIFRA does not constitute a waiver of any right to confidentiality that may exist under any other statute or in any other country.



SPONSOR CERTIFICATION OF GOOD LABORATORY PRACTICE

The Good Laboratory Practice Compliance Statement as defined by 40 CFR Part 160, found on page 33 of this volume is truthful and accurate.


W.T. Beidler

W. T. Beidler, Ph.D.  
Manager Environmental Residue Studies and  
Agent of Submitter/Sponsor

10-10-97


Date

Submitter/Sponsor: Novartis Crop Protection, Inc.  
PO Box 18300  
410 Swing Road  
Greensboro, North Carolina 27419



CERTIFICATION OF AUTHENTICITY

This report contains an unaltered copy of Ciba Analytical Method No. AG-666 (except for changes required to comply with PR Notice 86-5).

  
\_\_\_\_\_  
John D. Vargo, Ph.D., Scientist III  
Study Director  
Environmental Residue Studies  
Novartis Crop Protection, Inc.

10-10-97  
Date

## TABLE OF CONTENTS

Title Page.....	1
Statement of <u>No</u> Data Confidentiality Claim .....	2
Sponsor Certification of Good Laboratory Practice.....	3
Certification of Authenticity.....	4
Ciba Method No. AG-666, "Analytical Method for the Determination of GS-23199 and CGA-294849, Metabolites of CGA-215944 , in Soil by High Performance Liquid Chromatography with UV Detection Including Validation Data" .....	7

Ciba Analytical Method No. AG-666

Analytical Method for the Determination of GS-23199 and  
CGA-294849, Metabolites of CGA-215944 , in Soil by High Performance  
Liquid Chromatography with UV Detection Including Validation Data

ANALYTICAL METHOD FOR THE DETERMINATION OF  
GS-23199 AND CGA-294849, METABOLITES OF CGA-215944,  
IN SOIL BY HIGH PERFORMANCE LIQUID CHROMATOGRAPHY  
WITH UV DETECTION INCLUDING VALIDATION DATA

METHOD NO. AG-666

SPONSOR AND TESTING FACILITY:

Ciba-Geigy Corporation  
Ciba Crop Protection  
Environmental Fate and Effects Department  
Environmental Residue Studies  
410 Swing Road  
P. O. Box 18300  
Greensboro, NC 27419-8300

Project Number: 344001

Protocol Number: 376-96  
and Amendment 1

Study Initiation Date: December 2, 1996

Study Director:  
J. D. Vargo, Ph.D.

Approved By:  
W. T. Beidler

Title: Scientist III,  
Environmental Residue  
Studies

Title: Manager,  
Environmental Residue  
Studies

Signature: *John D. Vargo*

Signature: *W.T. Beidler*

Completion Date: 3-7-97

Date: 3-7-97



TABLE OF CONTENTS

	<u>PAGE NO.</u>
I. INTRODUCTION/SUMMARY .....	6
A. SCOPE .....	6
B. PRINCIPLE .....	6
II. MATERIALS AND METHODS .....	7
A. APPARATUS .....	7
B. REAGENTS AND ANALYTICAL STANDARDS .....	9
C. SAFETY AND HEALTH .....	11
D. ANALYTICAL PROCEDURE .....	11
1.0 SOIL MOISTURE DETERMINATION .....	11
2.0 SOIL EXTRACTION/CLEANUP .....	12
E. INSTRUMENTATION .....	17
1.0 DESCRIPTION AND OPERATING CONDITIONS: HPLC .....	17
2.0 CALIBRATION AND STANDARDIZATION .....	17
F. INTERFERENCES .....	17
G. CONFIRMATORY TECHNIQUES .....	17
H. TIME REQUIRED .....	18
I. MODIFICATIONS AND POTENTIAL PROBLEMS .....	18
J. PREPARATION OF STANDARD SOLUTIONS .....	19

8

TABLE OF CONTENTS  
(Continued)

	<u>PAGE NO.</u>
K. METHODS OF CALCULATION .....	21
1.0 DETERMINATION OF RESIDUES IN SAMPLES.....	21
2.0 DETERMINATION OF RESIDUES IN FORTIFIED SAMPLES.....	21
3.0 CALCULATIONS.....	22
III. RESULTS AND DISCUSSION .....	24
IV. CONCLUSION .....	26
V. CERTIFICATION .....	27
VI. CERTIFICATION OF GOOD LABORATORY PRACTICES .....	27
VII. QUALITY ASSURANCE STATEMENT .....	28
VIII. TABLES AND FIGURES .....	29
TABLE I. SOIL CHARACTERIZATION.....	29
TABLE II. HPLC SYSTEM AND OPERATING CONDITIONS: C18 ANALYTICAL COLUMN.....	30
TABLE III. HPLC SYSTEM AND OPERATING CONDITIONS: CN CONFIRMATION COLUMN.....	31
TABLE IV. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: C18 COLUMN.....	32
TABLE V. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: CN COLUMN.....	36

TABLE OF CONTENTS  
(Continued)

	<u>PAGE</u> <u>NO.</u>
TABLE VI. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: C18 COLUMN .....	38
TABLE VII. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: CN COLUMN .....	42
TABLE VIII. CALIBRATION AND RECOVERY DATA FOR FORTIFIED NEW YORK SOIL: C18 COLUMN .....	44
TABLE IX. CALIBRATION AND RECOVERY DATA FOR FORTIFIED NEW YORK SOIL: CN COLUMN .....	46
TABLE X. SUMMARY DATA FOR FORTIFIED SOIL: C18 COLUMN .....	48
TABLE XI. SUMMARY DATA FOR FORTIFIED SOIL: CN COLUMN .....	50
FIGURE 1. CHEMICAL NAMES AND STRUCTURES .....	52
FIGURE 2. AG-666 FLOW DIAGRAM FOR SOIL .....	53
FIGURE 3. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS: C18 ANALYTICAL COLUMN .....	54
FIGURE 4. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: C18 ANALYTICAL COLUMN .....	57
FIGURE 5. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: C18 ANALYTICAL COLUMN .....	60
FIGURE 6. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: C18 ANALYTICAL COLUMN .....	63

TABLE OF CONTENTS  
(Continued)

	<u>PAGE</u> <u>NO.</u>
FIGURE 7. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS: CN CONFIRMATION COLUMN .....	66
FIGURE 8. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: CN CONFIRMATION COLUMN .....	69
FIGURE 9. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: CN CONFIRMATION COLUMN .....	72
FIGURE 10. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: CN CONFIRMATION COLUMN .....	75
IX. REFERENCES .....	78

I. INTRODUCTION/SUMMARY

A. Scope

This method is used for the determination of GS-23199 and CGA-294849, metabolites of CGA-215944, in soil. The compounds are separated by high performance liquid chromatography (HPLC) and detected by UV absorption detection. The structures, chemical names, and Chemical Abstracts Registry numbers of the analytes are presented in Figure 1.

The limit of detection (smallest standard amount injected during the chromatographic run) is 2 ng for both analytes for the C18 analytical analyses and 1 ng for both analytes for the CN confirmation analyses. The limit of determination (the lowest fortification specified by the method which gives adequate recovery according to EPA guidelines) is 10 ppb in soil.

This is the second method issued for the analysis of the metabolites GS-23199 and CGA-294849 in soil (the first method issued was AG-653). After issuance of method AG-653, several problematic soils were encountered which contained significant interferences for the analytes. In this method, additional cleanup steps have been added that remove more of the potential interferences, in addition to a new reversed phase HPLC analysis system which permits the resolution of those interferences which still are present after the additional cleanup steps. While method AG-653 performed well with the soils used for method validation, this new method will do a much better job of removing and/or resolving interferences that may occur when different soils are encountered.

B. Principle

Soil samples (20 g) are reflux extracted with 20% (v/v) water/methanol. The samples are centrifuged and filtered. Methanol is removed via rotary evaporation until approximately 20 mL of extract remains. The extract is made basic with ammonium

hydroxide and then passed through a SAX solid phase extraction (SPE) column (the analytes are not retained). The non-retained fraction is placed on a rotary evaporator and methanol is removed until only aqueous remains. The aqueous is acidified and then passed through a C18 SPE attached piggy back style to a SCX SPE (the analytes are not retained.) The non-retained eluate is placed on a rotary evaporator and solvent is removed until only 2-4 mL of water remains. Methanol is added and the residue is transferred to a concentration tube where the methanol is removed until only aqueous remains. Water is added to adjust the final sample volume to a precalibrated mark. The sample is injected onto a reversed phase HPLC system with the analytes detected by UV absorbance. A flow diagram for the method is presented in Figure 2.

## II. MATERIALS AND METHODS

### A. Apparatus

- 1.0 Balance, analytical (Sartorius R160P) or equivalent.
- 2.0 Beaker, glass, 150-ml (Fisher cat. #02-540J) or equivalent.
- 3.0 Bottle, amber Boston round, with Polyseal-lined cap (Fisher cat. #05-563-2E) or equivalent.
- 4.0 Bottle, polypropylene, (Fisher cat. #05-562-23) or equivalent with cap. Appropriate size for soil extractions. Must be centrifugable.
- 5.0 Centrifuge, Sorvall Superspeed RC5-B (DuPont Instruments cat. #55228-9) or equivalent, with 6-place GSA rotor head (DuPont, Sorvall GSA cat. #08136) or equivalent.
- 6.0 Concentration tube, 50-mL (Fisher cat. #05-538-40B) or equivalent.

- 7.0 Cylinder, graduated, 50-ml, 100-mL, and 1000-mL (Fisher cat. #08-556C, #08-556D, #08-556G), or equivalent.
- ✓8.0 Filter, paper, for filtering soil extracts prior to rotary evaporation, 24-cm prepleated circles, Reeve Angel 802 (Fisher cat. #09-832D) or equivalent.
- ✓9.0 Filter, sample, for filtering final sample prior to analysis, Whatman Anotop 25 Inorganic Membrane Filter, 0.2  $\mu$ m pore, 25 mm diameter (Whatman cat. #6809-2022).
- 10.0 Flasks, round bottom, 250-ml (Fisher cat. #10-067E) and 100-mL (Fisher cat. #10-067D), or equivalent.
- needed* 11.0 Funnel, filter, 147-mm (Fisher cat. #10-373B) or equivalent.
- 12.0 Mixer, vortex (Fisher cat. #12-810-10) or equivalent.
- 13.0 Pasteur pipet, disposable (Fisher cat. #13-678-7C) or equivalent.
- 14.0 Pipets, glass, class A certified, assorted volumes. These pipets are used when an exact addition of liquid is required (i.e., final addition of solvent to samples).
- 15.0 Pipetters, Oxford BenchMate adjustable, 40-200  $\mu$ l volume range (Fisher cat. #21-231), 200-1000  $\mu$ l volume range (Fisher cat. #21-229) or equivalent. (Note: These adjustable pipetters may only be used for addition of liquid where an exact volume added is not critical, i.e., addition of acid.)
- 16.0 Rotary evaporator, Buchi (Fisher cat. #09-548-105F) or equivalent, with rotary evaporator traps (Fisher cat. #K570210-0124) or equivalent.

- 17.0 Ultrasonic bath, (Fisher cat. #15-336-6) or equivalent.
- 18.0 Vials, 1.5-ml (Sun Brokers, Inc. cat. #200-002) or equivalent, with Teflon-lined, crimp-top seals (Sun Brokers, Inc. cat. #200-152) or equivalent.

B. Reagents and Analytical Standards

All reagents are stored at room temperature. Solid analytical standards are stored in a freezer (temperature <-10°C).

- 1.0 Acetonitrile, HPLC grade (Fisher cat. #A998-4) or equivalent.
- 2.0 Ammonium hydroxide, certified ACS plus grade (Fisher cat. #A669S-500) or equivalent.
- 3.0 Extraction solvent: 20% (v/v) water in methanol. Add 800 mL of methanol to 200 mL of purified water.
- 4.0 C18 SPE extraction column, 1 gram size (Varian cat. #1225-6001) or equivalent.
- 5.0 SAX SPE extraction column, 1 gram size (Varian cat. #1225-6013) or equivalent.
- 6.0 SCX SPE extraction column, 1 gram size (Varian cat. #1225-6011) or equivalent.
- 7.0 Formic acid, 90%, laboratory grade (Fisher cat. #A119P-500) or equivalent.
- 8.0 Formic acid, 0.1%: mix 1.0 mL of formic acid with 999 mL of purified water.
- 9.0 Hexane, HPLC grade (Fisher cat. #H302-4) or equivalent.
- 10.0 Methanol, HPLC grade (Fisher cat. #A452-4) or equivalent.



- 11.0 Mobile phase (A) for C18 column: 3% methanol/water. Mix 30 mL of methanol with 970 mL of purified water.
- 12.0 Mobile phase (B) for C18 column: 30% water/acetonitrile. Mix 700 mL of acetonitrile with 300 mL of water.
- 13.0 Mobile phase (A) for CN column: 5/10/85% methanol/1-propanol/hexane. Mix 50 mL of methanol with 100 mL of 1-propanol. Add 850 mL of hexane to this mixture and mix the contents.
- 14.0 Mobile phase (B) for CN column: 20/20/60% methanol/1-propanol/hexane. Mix 200 mL of methanol with 200 mL of 1-propanol. Add 600 mL of hexane to this mixture and mix the contents.
- 15.0 1-Propanol, certified grade (Fisher cat. #A414-4), or equivalent.
- 16.0 Rinse solution for C18 SPE: 0.1/5/94.9% formic acid/methanol/water. Mix 1.0 mL of formic acid with 50 mL of methanol and 950 mL of water.
- 17.0 Rinse solution for SAX SPE: 0.5/20/79.5% ammonium hydroxide/water/methanol. Mix 5.0 mL of ammonium hydroxide with 995 mL of the soil extraction solvent (20/80% water/methanol).
- 18.0 Sample diluent for C18 analytical analysis: 5% methanol/water.
- 19.0 Sample diluent for CN confirmation analysis: 5/10/85% methanol/1-propanol/hexane. Combine 10 mL of 1-propanol with 5 mL of methanol. Add 85 mL of hexane. Mix the contents.
- 20.0 Water, HPLC grade, purified in-house with a HYDRO™ purification system or equivalent.

21.0 GS-23199 and CGA-294849, Ciba-Geigy Corp.,  
P. O. Box 18300, Greensboro, NC 27419-8300.

C. Safety and Health

Whereas most of the chemicals used and analyzed for in this method have not been completely characterized, general laboratory safety is advised (e.g., safety glasses, gloves, etc. should be used). Formic acid and ammonium hydroxide are irritants and should be used in a well-ventilated area (i.e., a fume hood).

D. Analytical Procedure

1.0 Soil Moisture Determination

Soil characterization data for the soils used in this validation study are presented in Table I.

- 1.1 Label and record the actual weight of an appropriate-sized glass beaker or aluminum weighing pan that will be used to determine the soil moisture content.
- 1.2 Add approximately 10-20 g of soil sample to the beaker or pan. Record the weight of the container plus wet soil.
- 1.3 Place the sample in an oven set at 100-120°C and let it dry overnight, or 12-16 hours.
- 1.4 Remove the sample and allow it to cool to room temperature.
- 1.5 Record the weight of the container plus dry soil.
- 1.6 Calculate the moisture content using the equation:

$$m = \frac{W_{1.2} - W_{1.5}}{W_{1.2} - W_{1.1}}$$

where  $m$  is the moisture content expressed in decimal form (i.e., 0.1 = 10%),  $W_{1.1}$  is the weight of the container (from Step 1.1),  $W_{1.2}$  is the weight of wet soil plus container (from Step 1.2), and  $W_{1.5}$  is the weight of the dry soil plus container (from Step 1.5).

## 2.0 Soil Extraction/Cleanup

Soil samples must be homogenized prior to analysis using suitable sample preparation techniques.

- 2.1 Weigh and record  $20 \pm 0.1$  g of soil sample and place in a 250-mL round bottom flask.
- 2.2 Sample fortification, if required for this particular sample, is to be done at this time (refer to Section II.K.2.0).
- 2.3 Add 100 mL of the soil extraction solvent. Swirl the contents briefly. Attach a reflux condenser to the flask and heat under reflux for one hour. Permit the extract to cool prior to centrifugation and filtering.
- 2.4 Transfer the sample to an appropriate size polypropylene centrifugeable bottle. Centrifuge the sample at approximately 9,000 RPM for 10 minutes, or at an alternate speed and time if the results are considered satisfactory.
- 2.5 Decant the sample extract through filter paper into a 100-mL graduated cylinder. Record the volume of extract.
- 2.6 Transfer the sample to a 250-mL round bottom flask. Rinse the graduated cylinder with approximately 5 mL of methanol and add to the sample. Add approximately 5 mL of water to each

sample to help prevent it from going dry during the rotary evaporation step.

- 2.7 Place the sample on a rotary evaporator with a water bath temperature of approximately 40 to 45°C. Use a solvent trap to minimize losses due to bumping. (Note: Periodic venting of the sample is required to prevent losses due to bumping.) Remove the methanol until approximately 20 mL of extract remains.
- 2.8 Remove the sample from the rotary evaporator. Add 100 µL of ammonium hydroxide.
- 2.9 Pass the sample through a preconditioned SAX SPE column, collecting the non-retained eluate in a 100-mL round bottom flask. A vacuum SPE reservoir may be used to improve flow through the SPE column. The sample loading speed should not exceed a fast drip rate. (Preconditioning of the SAX SPE is done by passing one column volume of methanol and then one column volume of the basic SAX rinse solution through the column. Do not permit the column to go dry before passing the sample through the SPE column.) Add approximately 5 mL of the basic SAX SPE rinse solution to the round bottom flask previously containing the sample. Swirl to rinse and dissolve any residues remaining on the glass. Transfer the rinsate via disposable pipette to the SAX SPE column and pass through, collecting in the 100-mL round bottom flask with the rest of the sample.
- 2.10 Place the sample on a rotary evaporator, water bath temperature of approximately 40 to 45°C, and remove methanol until only aqueous remains. Do not permit the extract to dry. A final volume of approximately 10 mL is desired.

(Purified water may be added to ensure the extract does not dry during this step.)

- 2.11 Remove the sample from the rotary evaporator. Add 50  $\mu$ L of formic acid to each sample. Sonicate and vortex samples to ensure all residue is in solution.
- 2.12 Load the sample onto a preconditioned C18 SPE extraction column attached piggy back style to a SCX SPE column. Collect the non-retained eluate, which contains both analytes, in a 100-mL round bottom flask. (Note: The SPE columns are preconditioned by passing one column volume each of methanol and 0.1% formic acid through the columns. Discard the rinse solutions. Add approximately 2 mL of 0.1% formic acid to the lower SCX column to ensure it will not dry while the sample is loaded and eluted.) The sample loading speed should not exceed a fast drip rate.
- 2.13 Add approximately 5 mL of 0.1/5/94.9% <sup>(C18 & WFE)</sup> formic acid/methanol/water to the 100-mL round bottom flask in which the rotary evaporation step was done. Vortex the solvent along the sides of the flask to dissolve any residues. Load this rinse onto the SPE columns and collect the eluate in the flask containing the sample from Step 2.12. Disconnect the C18 SPE column from the SCX column and then rinse the SCX column with one column volume of methanol, collecting the eluate with the sample in the 100-mL round bottom flask.
- 2.14 Place the sample on a rotary evaporator with a water bath temperature of approximately 40 to 45°C and remove solvent until approximately 3-5 mL of solvent remains. Do not permit the

sample to go dry. Use methanol to azeotrope the water, if needed.

- 2.15 Add 5 mL of methanol to the sample. Sonicate and vortex the flask well to ensure all of the residue is dissolved. Transfer the sample to a precalibrated 50-mL concentration tube using a disposable Pasteur pipet. Repeat with a second 5-mL rinse of the flask with methanol and add the rinsate to the sample. (Precalibration of the concentration tube is done by marking the meniscus line on each concentration tube after adding 4.0 mL of the sample diluent.)
- 2.16 Remove methanol from the sample until 2-3 mL of water remains. Water may be added as needed to prevent the sample from drying. (Note: It is important that all of the methanol be removed from the sample to ensure good chromatographic peak shape on the C18 analysis column.)
- 2.17 Add 0.2 mL of methanol to the sample. Dilute the sample with water to the 4.0-mL calibration mark. Additional dilution of the sample may be done using 5% methanol/water, if necessary.
- 2.18 Analyze using the reversed phase C18 HPLC system with UV detection. Refrigerate the sample if it will not be analyzed the same day the sample was processed.
- 2.19 For confirmation analyses using the CN normal phase HPLC system, first measure an aliquot of sample (from Step 2.17) and place in a 100-mL round bottom flask (a minimum volume of 3 mL is recommended). Add 20 mL of acetonitrile to the sample. (This is done to azeotrope the water from the sample.)

Remove the solvent from the sample via rotary evaporation with a water bath temperature of approximately 35°C until approximately 2-3 mL of sample remain. **(Warning: Do not let the sample go dry or significant losses of CGA-294849 will occur.)** Add approximately 20 mL of acetonitrile and evaporate solvent until approximately 2-3 mL remains. Repeat a third time with another 20 mL of acetonitrile, removing solvent until 2-3 mL remains. Add 5 mL of methanol to the sample, swirl the contents, and transfer to a 50-mL concentration tube via disposable pipette that has been precalibrated for the volume of sample being processed. Rinse the round bottom with 5 mL of methanol and transfer to the sample tube. Add approximately 1 mL of 1-propanol to the sample. Place the sample on a rotary evaporator with a water bath temperature of approximately 35°C and remove solvent until a volume of 1-propanol remains that is approximately 15% of the total sample volume that was processed. (Example: For a 3.0 mL sample, evaporate solvent until approximately 0.45 mL remains.) Add several drops (3-5) of methanol to the sample. Dilute to the calibration volume mark using hexane. Sonicate and vortex mix the sample thoroughly. Additional dilutions may be done using the CN sample diluent, if desired. Filter the sample with an Anotop sample filter, if necessary. (Note: A white suspended residue is frequently observed in soil samples prior to filtration. This residue will precipitate to the bottom of the vial in a short period of time and may not require filtration.) Analyze the sample by normal phase HPLC with the CN column and UV detection. Refrigerate the sample if it will not be analyzed the same day the sample was processed.

E. Instrumentation

1.0 Description and Operating Conditions: HPLC

See Tables II and III for descriptions of the HPLC systems for the C18 analytical and CN confirmation systems.

2.0 Calibration and Standardization

2.1 Determine the retention time of the analytes by injecting a standard solution into the HPLC. During a series of analyses, the analyte retention time should vary no more than 2% from its mean value, on a daily basis.

2.2 Calibrate the instrument by constructing a calibration curve from detector response (chromatographic peak height or area) and the amount of analyte injected. The response curve can be constructed manually or, preferably, by generation of a linear regression equation by use of a computer or appropriate calculator. Typical standard calibrations will be presented in Tables IV - IX for the analytical and confirmation systems. Typical standard chromatograms will be presented in Figure 3 for the C18 column and in Figure 7 for the CN column.

F. Interferences

1.0 There are no known interferences originating from the sample cleanup procedure. However, interferences can originate from impure chemicals, solvents, contaminated glassware, and the HPLC water supply.

G. Confirmatory Techniques

1.0 Confirmation of residues can be obtained by analysis of the samples on a second HPLC



system which utilizes a cyano (CN) column.  
See Table III for a description.

H. Time Required

- 1.0 The sample extraction and cleanup procedure can be completed for a set of eight samples in an eight-hour working day.
- 2.0 Each HPLC analysis requires approximately 20 minutes for the C18 analysis system and 19 minutes for the CN confirmation system.

I. Modifications and Potential Problems

- 1.0 Analytical Method AG-666 was validated only for the soil types listed in the final method. Other soil types, or soil samples from different locations, may exhibit binding or interference problems which were not observed with these samples.
- 2.0 "Bumping" is observed for soil samples during the solvent removal steps via rotary evaporation. Periodic venting of the vacuum and the use of solvent traps helps minimize inadvertent losses during these steps. A 500-mL round bottom flask may be necessary for the initial evaporation of solvent from the raw extract (Step 2.7) if severe bumping or foaming occurs in a 250-mL round bottom flask.
- 3.0 No analyte stability or solubility problems have been observed when solutions have been prepared and stored as detailed in Section II.J.
- 4.0 No analytes have been observed binding to the Whatman Anotop 25 sample filters during the final sample filtration step. It is unknown whether the analytes will bind to other brands/types of sample filters.
- 5.0 The compositions of the mobile phases were optimized for the columns used for the

analytical and confirmation analysis systems. The compositions may need to be altered from the conditions used in this method if columns of different manufacture are used.

- 6.0 During method development trials it was noted that interferences were encountered with certain soils which did not permit quantitation with the CN confirmation analysis system. The interferences were resolved from the analytes, however, with the reversed phase C18 analysis system.
- 7.0 Both analytes are very polar and have minimal retention on reversed phase columns. The YMC ODS-AQ column offers superior retention and peak shape over all other brand columns that have been tried.
- 8.0 Significant losses of CGA-294849 may occur if the samples are permitted to go to complete dryness during rotary evaporation steps.

J. Preparation of Standard Solutions

Stock solutions are stored in amber bottles in a freezer (<-10°C) when not in use. Fortification and HPLC standards are stored refrigerated in amber bottles when not in use. No analyte stability or solubility problems have been observed in the standard solutions used in this study.

- 1.0 Prepare individual 100 ng/µL stock solutions for each analyte. Weigh approximately 10.0 mg of analyte. Determine the appropriate volume of methanol to add using the equation presented below. The concentration of the analytical standard is corrected for its chemical purity.

$$V \text{ (mL)} = \frac{w(\text{mg}) \times P}{C \text{ (ng/uL)}} \times 10^3$$

Where V is the volume of methanol needed; W is the weight, in mg, of the solid analytical standard; P is the purity, in decimal form,

of the analytical standard; C is the desired concentration of the final solution, in ng/ul; and  $10^3$  is a conversion factor.

For example:

The volume of methanol required to dilute 9.9 mg of an analyte, of 98.0% purity, to a final concentration of 100 ng/uL is:

$$V(\text{mL}) = \frac{9.9 \text{ mg} \times 0.98}{100 \text{ ng/uL}} \times 10^3 = 97.02 \text{ mL}$$

- 2.0 Prepare a 20 ng/ $\mu$ L mixed standard solution in methanol by pipetting 10.0 mL of each analyte (from its 100 ng/ $\mu$ L stock solution in Step 1.0) into a 50-mL volumetric flask and diluting to the mark with methanol. Store the solution in an appropriate size amber bottle. This solution is used to prepare all subsequent dilutions.
- 3.0 Fortification standards are prepared by dilution of the 20 ng/ul mixed standard with methanol. The concentrations of the solutions to be prepared will depend upon the desired fortification level(s). Fortification standards should be prepared such that no more than 1.0 mL of the fortification solution is added to a sample. (Example: For a 20 g soil sample, the addition of 1.0 mL of a 0.2 ng/ $\mu$ L fortification solution will result in a fortification level of 10 ppb.)
- 4.0 C18 analytical standard.  
A 1.0 ng/ $\mu$ L analytical standard for the C18 HPLC system is prepared by pipetting 1.0 mL of each 100 ng/ $\mu$ L stock solution along with 3.0 mL of methanol into a 100-mL volumetric flask and then diluting to the mark with purified water. Subsequent dilutions are prepared by dilution of this solution with 5% methanol/water.

5.0 CN confirmation standard.

A 1.0 ng/ $\mu$ L analytical standard for HPLC calibration use is prepared by pipetting 0.5 mL from each 100 ng/ $\mu$ L stock solution into a 50-mL volumetric flask, adding 1.5 mL of methanol, 5.0 mL of 1-propanol, and diluting to the 50-mL mark with hexane. Subsequent serial dilutions are made with 5/10/85% methanol/1-propanol/hexane to prepare additional calibration standards.

K. Methods of Calculation

1.0 Determination of Residues in Samples

1.1 Inject the sample solution from Step II.D.2.17 or II.D.2.19 into the analysis system. The sample solution may be diluted if the analyte response exceeds the range of the calibration curve. The amount of analyte injected (ng) is determined by entering the value of the chromatographic peak height, or area, in the calibration response curve (Step II.E.4.2) and calculating (by computer, calculator, or manual means) the corresponding value of nanograms injected. Typical chromatograms for control and fortified soil samples will be presented in Figures 4-6 and 8-10 for the C18 and CN columns, respectively.

2.0 Determination of Residues in Fortified Samples

Validate the method for each set of samples analyzed by including a control sample and one or more control samples fortified prior to the extraction procedure with 10 ppb or more of each analyte in soil.

2.1 Add an appropriate volume of a fortification solution (from Step II.J.3.0) to the sample prior to any of the cleanup steps. The total volume of

the added fortification solution should not exceed 1.0 mL.

- 2.2 Proceed with the sample cleanup procedure (Step II.D.2.3).

### 3.0 Calculations

Calculations may be performed by computer program or manually as follows (soil concentrations are based on their wet weight):

- 3.1 Calculate the analyte concentration (in ppb) for field samples from equation (1):

$$(1) \text{ ppb analyte} = \frac{\text{ng analyte found}}{\text{g sample injected}} \times \frac{1}{R}$$

where R is the recovery factor expressed in decimal form (i.e., 0.8 = 80%) and is calculated from equation (4), and the chemical purity of the analytical standard has been accounted for in the preparation of the standard solutions. The use of the "1/R" recovery correction factor is optional and left to the discretion of the study director.

The grams of sample injected for soil is calculated from equation (2).

$$(2) \text{ g sample injected} = \frac{g}{V_e + V_{std} + (m \times g)} \times \frac{V_a \cdot V_i}{V_f}$$

where, g is the grams of soil (wet weight) used,  $V_a$  is the aliquot volume of extracted sample used for analysis,  $V_e$  is the volume of extract solvent used,  $V_{std}$  is the volume (mL) of fortification standard added (if any),  $V_i$  is the volume (mL) injected onto the HPLC column, m is the percent moisture in the sample, expressed in decimal form

(ex. 0.1 = 10%), and  $V_f$  is the final volume (mL) of the cleaned-up sample (from Step II.D.2.17 or II.D.2.19). (Note: the term "(m x g)" is a dilution correction factor due to the moisture in the soil, where 1.0 g = 1.0 mL.

The recovery factor, expressed as a percentage (R%), is calculated from fortification experiments and is presented in equation (3).

$$(3) R\% = \frac{\text{ppb analyte found} - \text{ppb analyte (control)}}{\text{ppb analyte added}} \times 100\%$$

The amount (ppb) of analyte found is calculated from equation (4).

$$(4) \text{ppb analyte found} = \frac{\text{ng analyte found}}{\text{g sample injected}}$$

Residues of metabolites found in test samples may also be expressed as parent equivalents by multiplying the amount found by the ratio of the molecular weight of CGA-215944 to that of the metabolite (equation (5)).

$$(5) \text{ppb CGA-215944 equiv.} = \text{ppb metabolite} \times \frac{\text{MW (p)}}{\text{MW (m)}}$$

where MW(p) is the average molecular weight of CGA-215944 (217.2) and MW(m) is the average molecular weight of the metabolite, 127.1 for GS-23199, and 142.1 for CGA-294849.

- 3.2 The accuracy of the method is determined by the average recovery of the analytes fortified into the test substrate. The precision is estimated by the relative standard deviation of the determined concentration.

III. RESULTS AND DISCUSSION

This method was validated under Protocol 376-96<sup>1</sup> and Amendment 1 for the analysis of control and fortified control soil. The objective of Protocol 376-96 was to validate "Draft" Analytical Method AG-666 for the determination of GS-23199 and CGA-294849 in soil with a limit of determination 10 ppb.

Recovery data for fortified soil samples are presented in Tables IV, VI, and VIII for the C18 analysis column and in Tables V, VII, and IX for the CN confirmation column. These tables contain raw data for both the samples and calibration standards which permits the manual calculation of recovery values. (Attempts to duplicate the calculations will be subject to round-off errors.) Typical chromatograms for analytical standards and for the soil samples are presented in Figures 3-10.

Limits of detection (defined as the lowest standard used in the calibration curves) of 2 ng for the C18 column analyses and 1 ng for the CN column analyses were achieved. The limit of determination (defined as the lowest fortification which provides acceptable recovery values) was 10 ppb for both analysis systems. The accuracy of the method is measured by the mean recovery values obtained at each fortification level. The precision of the method is estimated by the relative standard deviations from the mean recovery values obtained at each fortification level.

Excellent mean recovery values were obtained at the three fortification levels for the C18 column analysis (Table X). At the 10 ppb fortification level, the average recoveries and percent relative standard deviations for data combined for the three soil types were  $91 \pm 16.1\%$  for GS-23199 (n=17) and  $84 \pm 11.0\%$  for CGA-294849 (n=20). At the 100 ppb fortification level, the average recoveries and percent relative standard deviations for data combined for the three soil types were  $84 \pm 21.9\%$  for GS-23199 (n=10) and  $83 \pm 6.8\%$  for CGA-294849 (n=10). At the 1000 ppb fortification level, the average recoveries and percent relative standard deviations for data combined for the three soil types were  $90 \pm 6.1\%$  for GS-23199 (n=10) and  $89 \pm$

4.1% for CGA-294849 (n=10). The relative standard deviation for GS-23199 at the 100 ppb fortification level improves from 21.9% to 16.8% if one uncharacteristically low recovery of 50% is omitted from the statistical calculations.

Small quantities of interfering residues were observed for GS-23199 in the soil control samples for California (2.7 ppb) and New York (1.6 ppb). A small interference was observed for CGA-294849 in the New York soil (1.9 ppb). Recovery values were corrected for the small amounts of interfering residues.

For the CN confirmation analyses, recovery values for GS-23199 were very good while recovery values for CGA-294849 were low, although acceptable, due to volatility losses during rotary evaporation steps as the samples are converted from aqueous to organic solvent (Table XI). At the 10 ppb fortification level, the average recoveries and percent relative standard deviations for data combined for the three soil types were  $99 \pm 26.3\%$  for GS-23199 (n=11) and  $70 \pm 9.8\%$  for CGA-294849 (n=9). The %RSD for GS-23199 improves from 26.3% to 13.3% if the recovery values of samples VAL23 (154%) and VAL24 (50%) are not used. These recovery values are outliers that are uncharacteristic. At the 100 ppb fortification level, the average recoveries and percent relative standard deviations for data combined for the three soil types were  $86 \pm 12.8\%$  for GS-23199 (n=6) and  $81 \pm 9.3\%$  for CGA-294849 (n=6). At the 1000 ppb fortification level, the average recoveries and percent relative standard deviations for data combined for the three soil types were  $93 \pm 17.8\%$  for GS-23199 (n=6) and  $86 \pm 13.7\%$  for CGA-294849 (n=6).

Small quantities of interfering residues were observed during CN analyses for GS-23199 in the soil control samples for Georgia (0.7 ppb) and New York (3.3 ppb). Recovery values were corrected for these residues. No interferences were observed in any soil for CGA-294849.

Method AG-666 proved to be accurate, rugged, and reliable for the C18 analyses of both analytes in soil. The CN confirmation analysis procedure was found to be not as reliable as the C18 analysis, but adequate as a



confirmation procedure. The highly polar nature of these two analytes makes the analysis technique very difficult and challenging.

An independent laboratory validation of this method will be conducted at a future date. This data will be included with the data package that will be submitted when the active ingredient is submitted for registration.

Reference substance ID, test system ID, protocol amendments, protocol deviations, and circumstances affecting the quality and integrity of data will be reported in the Residue Test Report<sup>2</sup>, RI-MV-016-96, Report Number 1. All raw data associated with this study and the original final report and protocol will be archived in the Agricultural Group Archive Facility at Ciba-Geigy Corporation, Greensboro, NC. All non-study specific data (i.e., instrument logbooks, etc.) will be stored in the previously mentioned archives when all entry pages are filled or when the logbook is replaced. Soil samples will be archived in the Biochemistry Group Sample Storage Facility, Greensboro, NC, until the registration studies have been accepted by the EPA and QA verification has been performed.

#### IV. CONCLUSION

Method AG-666 was found to be an accurate and reliable method for the determination of GS-23199 and CGA-294849 in soil.

V. CERTIFICATION

This report and experimental results included in this study, laboratory Project I.D. AG-666, are certified to be authentic accounts of the experiments.

3-7-97  
Date

John D. Vargo  
John D. Vargo, Ph.D.  
Scientist III  
Environmental Residue Studies  
Environmental Fate and  
Effects Department  
910-632-7525

Ciba Crop Protection  
Ciba-Geigy Corporation  
Post Office Box 18300  
Greensboro, NC 27419-8300

VI. CERTIFICATION OF GOOD LABORATORY PRACTICES

The analytical work reported in AG-666 was performed in accordance with Good Laboratory Practice Standards, 40 CFR Part 160.

3-7-97  
Date

John D. Vargo  
John D. Vargo, Ph.D.  
Study Director

3-7-97  
Date

W. T. Beidler  
W. T. Beidler  
Manager  
Environmental Residue Studies  
Environmental Fate and  
Effects Department  
910-632-2976

Ciba Crop Protection  
Ciba-Geigy Corporation  
Post Office Box 18300  
Greensboro, NC 27419-8300

VII.

**QUALITY ASSURANCE STATEMENT**

**Report Title:** Analytical Method for the Determination of GS-23199 and CGA-294849, Metabolites of CGA-215944, in Soil by High Performance Liquid Chromatography with UV Detection Including Validation Data

**Study Director:** J. D. Vargo

**Protocol Number:** 376-96

**Method Number:** AG-666

**Project Number:** 344001

Pursuant to Good Laboratory Practice Regulations, this statement verifies that the aforementioned study was inspected and/or audited and the findings reported to Management and to the Study Director by the Quality Assurance Unit on the dates listed below.

<b><u>INSPECTION/AUDIT TYPE</u></b>	<b><u>INSPECTION/AUDIT DATE(S)</u></b>	<b><u>REPORTING DATE</u></b>
Protocol Audit	11/11/96	11/11/96
In-Progress Inspection	12/18/96	12/18/96
Final Report Audit	2/5,6,10,11/97	2/11/97

Prepared by: Teresa S. Cox Date: 2/11/97

Teresa S. Cox  
Senior Quality Assurance Auditor  
Quality Assurance Unit  
Ciba Crop Protection  
Ciba-Geigy Corporation

VIII. TABLES AND FIGURES

TABLE I. SOIL CHARACTERIZATION

County	<u>California</u>	<u>Georgia</u>	<u>New York</u>
	Madera	Mitchell	Columbia
Soil Depth	0-6"	0-6"	0-6"
pH	5.6	6.8	5.8
Cation Exchange Capacity (meq/100 g)	7.2	6.8	7.77
% Organic Matter	0.2	0.9	2.35
% Water Holding Capacity @ 1/3 Bar	7.5	8.5	21.39
% Sand	71	80	39.2
% Silt	21	10	42.4
% Clay	8	10	18.4
Soil Classification	Sandy Loam	Loamy Sand	Loam
Bulk Density (g/cc)	1.43	1.38	1.26
Soil Moisture Percent	6.0	9.7	11.7

The soil samples were collected and the soil characterized under Protocol Numbers: 189-95 for New York, 190-95 for California, and 191-95 for Georgia. Soil moisture percentages will be obtained under Protocol 376-96.

TABLE II. HPLC SYSTEM AND OPERATING CONDITIONS:  
C18 ANALYTICAL COLUMN

Instrumentation:

Perkin-Elmer Model Series 410 Gradient Pump  
Perkin-Elmer Model ISS 200 Autosampler  
Eppendorf Model CH-30 Column Heater  
Perkin Elmer Model LC-95 UV Absorbance Detector

Operating Conditions:

Column Heater: 30°C  
Detection Wavelength: 265 nm  
Injection Volume: 100 µl  
Mobile Phase Flow Rate: 1.5 ml/min  
Column: YMC ODS-AQ (YMC, Inc.),  
25 cm x 4.6 mm, dp = 5 µm, equipped with an  
Upchurch (#A-318) pre-column filter (0.5 µm),  
or equivalent, and a YMC ODS-AQ guard column.  
Mobile Phase A: 3% methanol/water  
Mobile Phase B: 70% acetonitrile/water

Gradient Program:

<u>Time</u>	<u>%A</u>	<u>%B</u>
0	100	0
8	100	0
8.5	0	100
11.5	0	100
12	100	0
21.5	100	0

Total Run Time: 21.5 min.  
Analyte Retention Times:

CGA-294849 6.3 min  
GS-23199 9.9 min

TABLE III. HPLC SYSTEM AND OPERATING CONDITIONS:  
CN CONFIRMATION COLUMN

Instrumentation:

Perkin-Elmer Model Series 4 Gradient Pump  
Perkin-Elmer Model ISS 200 Autosampler  
Eppendorf Model CH-30 Column Heater  
Perkin Elmer Model LC-95 UV Absorbance Detector

Operating Conditions:

Column Heater: 30°C  
Detection Wavelength: 265 nm  
Injection Volume: 50 µl  
Mobile Phase Flow Rate: 1.5 ml/min  
Column: Spherisorb CN (Phase Separations, Inc.),  
25 cm x 4.6 mm, dp = 5 µm, equipped with an  
Upchurch (#A-318) pre-column filter (0.5 µm),  
or equivalent  
Mobile Phase A: 5/10/85% Methanol/1-Propanol/Hexane  
Mobile Phase B: 20/20/60% Methanol/1-Propanol/Hexane

Gradient Program:

<u>Time</u>	<u>%A</u>	<u>%B</u>
0	100	0
6.5	100	0
7	0	100
11	0	100
11.5	100	0
19	100	0

Total Run Time: 19 min.

Analyte Retention Times:

GS-23199 4.5 min  
CGA-294849 7.5 min

TABLE IV. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: C18 COLUMN

**GS-23199**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL1	0 (blank)	20.00	95.0	4	0.475	-	0	< 2	< 10	-
VAL2	0 (control)	20.00	90.0	4	0.445	-	0	< 2	< 10	-
VAL3	10	20.00	90.0	4	0.440	9.51	6334	3.43	7.78	78
VAL4	10	20.00	90.0	4	0.440	9.58	7171	3.96	8.98	90
VAL5	10	20.00	90.0	4	0.440	*	*	*	*	*
VAL6	10	20.00	90.0	4	0.440	*	*	*	*	*
VAL7	100	20.00	92.0	10	0.180	9.48	29363	18.0	100	100
VAL8	100	20.00	91.0	10	0.178	9.56	31358	19.3	108	108
VAL9	1000	20.00	89.0	50	0.035	9.62	51669	32.1	922	92
VAL10	1000	20.00	92.0	50	0.036	9.56	52299	32.5	903	90

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 6.0%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.1 mL

\* Analyte could not be quantitated due to interference.

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	9.57	4086	2.00
0.05	100	5	9.64	9092	5.17
0.1	100	10	9.62	16865	10.1
0.5	100	50	9.61	79103	49.5
1.0	100	100	9.50	159214	100

slope = 1579.25  
 Y-intercept = 925.214  
 corr. coeff. = 0.99998

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

TABLE IV. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: C18 COLUMN

**CGA-294849**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL1	0 (blank)	20.00	95.0	4	0.475	-	0	< 2	< 10	-
VAL2	0 (control)	20.00	90.0	4	0.445	-	0	< 2	< 10	-
VAL3	10	20.00	90.0	4	0.440	6.18	7080	4.10	9.30	93
VAL4	10	20.00	90.0	4	0.440	6.24	6306	3.61	8.20	82
VAL5	10	20.00	90.0	4	0.440	6.24	6981	4.03	9.16	92
VAL6	10	20.00	90.0	4	0.440	6.24	6289	3.60	8.17	82
VAL7	100	20.00	92.0	10	0.180	6.18	26305	16.2	89.9	90
VAL8	100	20.00	91.0	10	0.178	6.22	22993	14.1	79.2	79
VAL9	1000	20.00	89.0	50	0.035	6.27	52210	32.5	933	93
VAL10	1000	20.00	92.0	50	0.036	6.25	51371	32.0	887	89

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 6.0%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.1 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	6.27	3241	1.68
0.05	100	5	6.28	8850	5.21
0.1	100	10	6.26	16731	10.2
0.5	100	50	6.27	79944	49.9
1.0	100	100	6.22	159577	100

slope = 1589.84  
 Y-intercept = 567.829  
 corr. coeff. = 0.99999

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.



TABLE IV. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: C18 COLUMN

**GS-23199**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Sample Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL1C	0 (blank)	20.00	95.0	4	0.475	10.08	757	< 2 (0.69)	< 10 (1.46)	-
VAL2C	0 (control)	20.00	89.0	4	0.440	10.06	1505	< 2 (1.17)	< 10 (2.66)	-
VAL3C	10	20.00	92.0	4	0.450	9.96	7659	5.10	11.3	87
VAL4C	10	20.00	90.0	4	0.440	9.97	7725	5.15	11.7	90
VAL5C	10	20.00	90.5	4	0.443	9.92	6978	4.67	10.5	79
VAL6C	10	20.00	89.0	4	0.435	9.93	7061	4.72	10.8	82
VAL7C	100	20.00	89.5	10	0.175	9.97	29220	18.9	108	105
VAL8C	100	20.00	90.5	10	0.177	9.98	25968	16.8	94.9	92
VAL9C	1000	20.00	88.5	50	0.035	9.94	54098	34.8	1000	100
VAL10C	1000	20.00	92.0	50	0.036	9.96	53680	34.5	959	96

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 6.0%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.1 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	9.93	3060	2.16
0.05	100	5	9.94	7622	5.08
0.1	100	10	10.13	15029	9.81
0.5	100	50	10.06	77693	49.9
1.0	100	100	10.08	156255	100

slope = 1564.56  
 Y-intercept = -324.5281  
 corr. coeff. = 0.99999

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

TABLE IV. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: C18 COLUMN

**CGA-294849**

Sample Code	Amount Added (ppbl)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppbl)	% Recovery
VAL1C	0 (blank)	20.00	95.0	4	0.475	-	0	< 2	< 10	-
VAL2C	0 (control)	20.00	89.0	4	0.440	-	0	< 2	< 10	-
VAL3C	10	20.00	92.0	4	0.450	6.36	5607	3.47	7.72	77
VAL4C	10	20.00	90.0	4	0.440	6.34	5310	3.28	7.46	75
VAL5C	10	20.00	90.5	4	0.443	6.33	5361	3.32	7.49	75
VAL6C	10	20.00	89.0	4	0.435	6.33	6217	3.87	8.88	89
VAL7C	100	20.00	89.5	10	0.175	6.36	23229	14.8	84.4	84
VAL8C	100	20.00	90.5	10	0.177	6.38	23262	14.8	83.5	84
VAL9C	1000	20.00	88.5	50	0.035	6.37	51805	33.1	956	96
VAL10C	1000	20.00	92.0	50	0.036	6.37	51723	33.0	918	92

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
m = % moisture = 6.0%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.1 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	6.35	3894	2.38
0.05	100	5	6.35	7890	4.94
0.1	100	10	6.35	15535	9.84
0.5	100	50	6.36	77640	49.7
1.0	100	100	6.38	156409	100

slope = 1559.42  
Y-intercept = 189.023  
corr. coeff. = 0.99998

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

41

TABLE V. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: CN COLUMN

**GS-23199**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL1C	0 (blank)	20.00	95.0	4	0.237	-	0	< 1	< 10	-
VAL2C	0 (control)	20.00	89.0	4	0.220	-	0	< 1	< 10	-
VAL3C	10	20.00	92.0	4	0.225	4.38	4305	2.52	11.2	112
VAL4C	10	20.00	90.0	4	0.220	4.38	3755	2.18	9.88	99
VAL5C	*	*	*	*	*	*	*	*	*	*
VAL6C	10	20.00	89.0	4	0.218	4.39	3315	1.90	8.74	87
VAL7C	100	20.00	89.5	10	0.088	4.39	14549	8.88	101	101
VAL8C	100	20.00	90.5	10	0.089	4.39	13877	8.47	95.6	96
VAL9C	1000	20.00	88.5	50	0.017	4.39	32719	20.2	1170	117
VAL10C	1000	20.00	92.0	50	0.018	4.39	32312	19.9	1110	111

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 6.0%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.05 mL

\* Sample not analyzed due to lab prep accident destroying sample.

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	50	1	4.44	1612	0.844
0.05	50	2.5	4.40	3990	2.32
0.1	50	5	4.39	8677	5.24
0.5	50	25	4.39	40834	25.2
1.0	50	50	4.40	80510	49.9

slope = 1609.04  
 Y-intercept = 253.723  
 corr. coeff. = 0.99995

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

42

TABLE V. CALIBRATION AND RECOVERY DATA FOR FORTIFIED CALIFORNIA SOIL: CN COLUMN

**CGA-294849**

Sample Code	Amount Added (ppbl)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppbl)	% Recovery
VAL1C	0 (blank)	20.00	95.0	4	0.237	-	0	< 1	< 10	-
VAL2C	0 (control)	20.00	89.0	4	0.220	-	0	< 1	< 10	-
VAL3C	10	20.00	92.0	4	0.225	7.40	2539	1.40	6.22	62
VAL4C	10	20.00	90.0	4	0.220	7.38	2773	1.54	7.01	70
VAL5C	*	*	*	*	*	*	*	*	*	*
VAL6C	10	20.00	89.0	4	0.218	7.37	2441	1.34	6.15	61
VAL7C	100	20.00	89.5	10	0.088	7.38	11387	6.88	78.5	79
VAL8C	100	20.00	90.5	10	0.089	7.39	13585	8.24	93.0	93
VAL9C	1000	20.00	88.5	50	0.017	7.39	28232	17.3	999	100
VAL10C	1000	20.00	92.0	50	0.018	7.38	28439	17.4	968	97

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 $m = \% \text{ moisture} = 6.0\%$

$V_e = \text{volume of extraction solvent used} = 100 \text{ mL}$

$V_{std} = \text{volume of fortification solution} = 1 \text{ mL (for fortified samples only)}$

$V_i = \text{HPLC injection volume} = 0.05 \text{ mL}$

\* Sample not analyzed due to lab prep accident destroying sample.

**Calibration Standards**

Concentration (ng/ $\mu$ L)	Injection Volume ( $\mu$ L)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	50	1	7.46	1611	0.825
0.05	50	2.5	7.35	4134	2.39
0.1	50	5	7.37	8678	5.20
0.5	50	25	7.38	40994	25.2
1.0	50	50	7.38	80874	49.9

slope = 1615.50  
 Y-intercept = 279.287  
 corr. coeff. = 0.99996

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

43

TABLE VI. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: C18 COLUMN

**GS-23199**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL11	0 (blank)	20.00	95.0	4	0.475	-	0	< 2	< 10	-
VAL12	0 (control)	20.03	92.0	4	0.452	-	0	< 2	< 10	-
VAL13	10	20.03	93.0	4	0.452	9.43	5077	3.32	7.35	73
VAL14	10	20.10	90.0	4	0.439	9.57	5980	3.90	8.87	89
VAL15	10	20.08	92.0	4	0.449	9.51	8347	5.39	12.0	120
VAL16	10	20.02	93.0	4	0.452	9.48	6216	4.04	8.94	89
VAL17	100	20.01	91.0	10	0.177	9.48	18177	11.6	65.6	66
VAL18	100	20.01	92.0	10	0.179	9.41	14039	8.99	50.3	50
VAL19	1000	20.01	92.0	50	0.036	9.50	47769	30.3	848	85
VAL20	1000	20.00	91.0	50	0.035	9.53	48362	30.7	868	87

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 $m = \% \text{ moisture} = 9.7\%$

$V_e = \text{volume of extraction solvent used} = 100 \text{ mL}$

$V_{std} = \text{volume of fortification solution} = 1 \text{ mL (for fortified samples only)}$

$V_i = \text{HPLC injection volume} = 0.1 \text{ mL}$

**Calibration Standards**

Concentration (ng/mL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	9.60	2805	1.89
0.05	100	5	9.60	7680	4.97
0.1	100	10	9.60	15512	9.92
0.5	100	50	9.53	79552	50.4
1.0	100	100	9.50	157661	99.8

slope = 1581.52  
y-intercept = -180.5007  
corr. coeff. = 0.99998

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

44

TABLE VI. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: C18 COLUMN

**CGA-294849**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL11	0 (blank)	20.00	95.0	4	0.475	-	0	< 2	< 10	-
VAL12	0 (control)	20.03	92.0	4	0.452	-	0	< 2	< 10	-
VAL13	10	20.03	93.0	4	0.452	6.10	7085	4.64	10.3	103
VAL14	10	20.10	90.0	4	0.439	6.19	5891	3.89	8.86	89
VAL15	10	20.08	92.0	4	0.449	6.23	6570	4.32	9.63	96
VAL16	10	20.02	93.0	4	0.452	6.18	6088	4.02	8.88	89
VAL17	100	20.01	91.0	10	0.177	6.19	24961	15.9	89.7	90
VAL18	100	20.01	92.0	10	0.179	6.16	20587	13.1	73.3	73
VAL19	1000	20.01	92.0	50	0.036	6.19	51108	32.3	903	90
VAL20	1000	20.00	91.0	50	0.035	6.22	49016	31.0	876	88

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 $m = \% \text{ moisture} = 9.7\%$

$V_e$  = volume of extraction solvent used = 100 mL

$V_{std}$  = volume of fortification solution = 1 mL (for fortified samples only)

$V_i$  = HPLC injection volume = 0.1 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	6.25	2760	1.93
0.05	100	5	6.26	7966	5.20
0.1	100	10	6.26	15770	10.1
0.5	100	50	6.23	78638	49.6
1.0	100	100	6.23	159163	100

slope = 1591.75  
Y-intercept = -305.0952  
corr. coeff. = 0.99998

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

45

TABLE VI. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: C18 COLUMN

**GS-23199**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL11B	0 (blank)	20.00	95.5	4	0.477	-	0	< 2	< 10	-
VAL12B	0 (control)	20.01	88.5	4	0.434	-	0	< 2	< 10	-
VAL13B	10	20.00	90.0	4	0.437	9.37	6582	3.89	8.90	89
VAL14B	10	20.03	92.0	4	0.448	9.40	7472	4.46	9.96	100
VAL15B	10	20.05	93.5	4	0.455	9.29	7678	4.59	10.1	101
VAL16B	10	20.00	91.5	4	0.444	9.26	*	*	*	*
VAL17B	100	20.03	91.0	10	0.177	9.40	20660	12.8	72.4	72
VAL18B	100	20.03	91.5	10	0.178	9.39	21493	13.4	75.0	75
VAL19B	1000	20.05	91.0	50	0.035	9.38	47971	30.2	851	85
VAL20B	1000	20.00	91.0	50	0.035	9.37	49681	31.2	883	88

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 $m = \% \text{ moisture} = 9.7\%$

$V_e$  = volume of extraction solvent used = 100 mL

$V_{std}$  = volume of fortification solution = 1 mL (for fortified samples only)

$V_i$  = HPLC injection volume = 0.1 mL

\* Sample could not be accurately quantitated due to interference.

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	9.28	3716	2.07
0.05	100	5	9.29	8664	5.21
0.1	100	10	9.31	15845	9.77
0.5	100	50	9.31	79042	49.9
1.0	100	100	9.30	158190	100

slope = 1576.19  
Y-intercept = 446.599  
corr. coeff. = 0.99999

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

TABLE VI. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: C18 COLUMN

**CGA-294849**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL11B	0 (blank)	20.00	95.5	4	0.477	-	0	< 2	< 10	-
VAL12B	0 (control)	20.01	88.5	4	0.434	-	0	< 2	< 10	-
VAL13B	10	20.00	90.0	4	0.437	6.12	6611	3.90	8.93	89
VAL14B	10	20.03	92.0	4	0.448	6.10	6409	3.77	8.43	84
VAL15B	10	20.05	93.5	4	0.455	6.04	6734	3.98	8.74	87
VAL16B	10	20.00	91.5	4	0.444	6.01	6559	3.87	8.71	87
VAL17B	100	20.03	91.0	10	0.177	6.12	24690	15.4	86.8	87
VAL18B	100	20.03	91.5	10	0.178	6.11	23517	14.6	82.1	82
VAL19B	1000	20.05	91.0	50	0.035	6.13	47240	29.7	837	84
VAL20B	1000	20.00	91.0	50	0.035	6.12	48900	30.7	869	87

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 9.7%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.1 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	6.16	3703	2.06
0.05	100	5	6.15	8152	4.88
0.1	100	10	6.13	16608	10.2
0.5	100	50	6.15	78826	49.7
1.0	100	100	6.13	158422	100

slope = 1577.46  
 Y-intercept = 455.075  
 corr. coeff. = 0.99999

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

47



TABLE VII. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: CN COLUMN

**GS-23199**

Sample Code	Amount Added (ppbl)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppbl)	% Recovery
VAL11B	0 (blank)	20.00	95.5	4	0.239	4.49	393	< 1 (0.17)	< 10 (0.72)	-
VAL12B	0 (control)	20.01	88.5	4	0.217	4.56	350	< 1 (0.15)	< 10 (0.68)	-
VAL13B	10	20.00	90.0	4	0.219	4.36	3695	2.08	9.54	89
VAL14B	10	20.03	92.0	4	0.224	4.34	4066	2.30	10.3	96
VAL15B	10	20.05	93.5	4	0.228	4.36	3419	1.92	8.45	78
VAL16B	10	20.00	91.5	4	0.222	4.34	4196	2.37	10.7	100
VAL17B	100	20.03	91.0	10	0.089	4.35	11174	6.41	72.4	72
VAL18B	100	20.03	91.5	10	0.089	4.33	12221	7.02	78.9	78
VAL19B	1000	20.05	91.0	50	0.018	4.33	25748	14.9	838	84
VAL20B	1000	20.00	91.0	50	0.018	4.34	25499	14.7	832	83

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 9.7%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.05 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	50	1	4.47	1962	1.08
0.05	50	2.5	4.48	4504	2.55
0.1	50	5	4.36	8896	5.10
0.5	50	25	4.35	42542	24.6
1.0	50	50	4.34	86815	50.2

slope = 1727.48  
 Y-intercept = 95.0185  
 corr. coeff. = 0.99993

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

48

TABLE VII. CALIBRATION AND RECOVERY DATA FOR FORTIFIED GEORGIA SOIL: CN COLUMN

**CGA-294849**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL11B	0 (blank)	20.00	95.5	4	0.239	-	0	< 1	< 10	-
VAL12B	0 (control)	20.01	88.5	4	0.217	-	0	< 1	< 10	-
VAL13B	10	20.00	90.0	4	0.219	7.32	3207	1.85	8.44	84
VAL14B	10	20.03	92.0	4	0.224	7.30	2912	1.68	7.48	75
VAL15B	10	20.05	93.5	4	0.228	7.36	1828	1.05	4.61	46*
VAL16B	10	20.00	91.5	4	0.222	7.32	1940	1.11	5.01	50*
VAL17B	100	20.03	91.0	10	0.089	7.31	12485	7.20	81.4	81
VAL18B	100	20.03	91.5	10	0.089	7.27	10834	6.25	70.2	70
VAL19B	1000	20.05	91.0	50	0.018	7.28	22731	13.1	740.	74
VAL20B	1000	20.00	91.0	50	0.018	7.29	21621	12.5	706	71

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 9.7%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.05 mL

\* Samples went to dryness during final rotary evaporation. Losses of CGA-294849 seen when this happens.

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	50	1	7.48	1643	0.942
0.05	50	2.5	7.53	4481	2.58
0.1	50	5	7.30	8978	5.18
0.5	50	25	7.28	42684	24.6
1.0	50	50	7.28	86876	50.2

slope = 1731.79  
 Y-intercept = 11.6827  
 corr. coeff. = 0.99994

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

TABLE VIII. CALIBRATION AND RECOVERY DATA FOR FORTIFIED NEW YORK SOIL: C18 COLUMN

**GS-23199**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL21	0 (blank)	20.00	93.0	4	0.465	-	0	< 2	< 10	-
VAL22	0 (control)	20.02	86.0	4	0.421	9.27	1046	< 2 (0.66)	< 10 (1.56)	-
VAL23	10	20.00	89.5	4	0.433	8.96	6839	4.36	10.1	85
VAL24	10	20.00	89.5	4	0.433	9.02	5765	3.68	8.49	69
VAL25	10	20.00	88.5	4	0.428	8.89	8450	5.39	12.6	110
VAL26	10	20.02	89.0	4	0.431	8.92	9061	5.79	13.4	119
VAL27	100	20.00	90.0	10	0.174	8.98	23798	15.2	87.4	86
VAL28	100	20.00	88.0	10	0.170	8.93	24032	15.4	90.2	89
VAL29	1000	20.00	90.0	50	0.035	8.90	48976	31.3	899	90
VAL30	1000	20.00	90.0	50	0.035	8.89	44478	28.4	817	82

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 $m = \% \text{ moisture} = 11.7\%$

$V_e$  = volume of extraction solvent used = 100 mL

$V_{std}$  = volume of fortification solution = 1 mL (for fortified samples only)

$V_i$  = HPLC injection volume = 0.1 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	8.99	3143	2.00
0.05	100	5	8.98	8206	5.24
0.1	100	10	8.96	15567	9.95
0.5	100	50	8.94	77602	49.6
1.0	100	100	8.92	156556	100

slope = 1562.68  
Y-intercept = 21.1470  
corr. coeff. = 0.99998

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

TABLE VIII. CALIBRATION AND RECOVERY DATA FOR FORTIFIED NEW YORK SOIL: C18 COLUMN

**CGA-294849**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL21	0 (blank)	20.00	93.0	4	0.465	-	0	< 2	< 10	-
VAL22	0 (control)	20.02	86.0	4	0.421	5.83	1358	< 2 (0.82)	< 10 (1.94)	-
VAL23	10	20.00	89.5	4	0.433	5.90	6084	3.83	8.85	69
VAL24	10	20.00	89.5	4	0.433	5.94	6098	3.84	8.87	69
VAL25	10	20.00	88.5	4	0.428	5.85	6709	4.23	9.88	79
VAL26	10	20.02	89.0	4	0.431	5.86	6241	3.93	9.13	72
VAL27	100	20.00	90.0	10	0.174	5.88	23237	14.8	84.9	83
VAL28	100	20.00	88.0	10	0.170	5.86	20946	13.3	78.2	76
VAL29	1000	20.00	90.0	50	0.035	5.86	46991	30.0	860	86
VAL30	1000	20.00	90.0	50	0.035	5.88	47273	30.1	865	86

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 11.7%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.1 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	100	2	5.89	3174	1.98
0.05	100	5	5.88	7975	5.04
0.1	100	10	5.88	15769	10.0
0.5	100	50	5.88	78294	49.9
1.0	100	100	5.90	156740	100

slope = 1566.18  
 y-intercept = 80.2287  
 corr. coeff. = 1.00000

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

51

TABLE IX. CALIBRATION AND RECOVERY DATA FOR FORTIFIED NEW YORK SOIL: CN COLUMN

**GS-23199**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL21	0 (blank)	20.00	93.0	4	0.233	4.45	539	< 2 (0.09)	< 10 (0.37)	-
VAL22	0 (control)	20.02	86.0	4	0.210	4.34	1540	< 2 (0.69)	< 10 (3.27)	-
VAL23	10	20.00	89.5	4	0.217	4.35	7109	4.04	18.7	154
VAL24	10	20.00	89.5	4	0.217	4.36	3362	1.79	8.24	50
VAL25	10	20.00	88.5	4	0.214	4.34	5404	3.01	14.1	108
VAL26	10	20.02	89.0	4	0.216	4.36	5865	3.29	15.3	120
VAL27	100	20.00	90.0	10	0.087	4.36	13350	7.80	89.5	86
VAL28	100	20.00	88.0	10	0.085	4.35	12482	7.27	85.4	82
VAL29	1000	20.00	90.0	50	0.017	4.36	24779	14.7	842	83
VAL30	1000	20.00	90.0	50	0.017	4.36	23422	13.9	796	79

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 m = % moisture = 11.7%

V<sub>e</sub> = volume of extraction solvent used = 100 mL

V<sub>std</sub> = volume of fortification solution = 1 mL (for fortified samples only)

V<sub>i</sub> = HPLC injection volume = 0.05 mL

**Calibration standards**

Concentration (ng/mL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	50	1	4.36	1875	0.89
0.05	50	2.5	4.35	4532	2.49
0.1	50	5	4.36	9250	5.33
0.5	50	25	4.36	41351	24.6
1.0	50	50	4.36	83725	50.1

slope = 1661.70  
 Y-intercept = 396.410  
 corr. coeff. = 0.99992

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

52

TABLE IX. CALIBRATION AND RECOVERY DATA FOR FORTIFIED NEW YORK SOIL: CN COLUMN

**CGA-294849**

Sample Code	Amount Added (ppb)	Sample Weight Extr. (g)	(V <sub>a</sub> ) Aliquot Volume (mL)	(V <sub>f</sub> ) Final Volume (mL)	Sample Wt. Inj. (g)	Retention Time (min)	Peak Area	Analyte Found (ng)	Residue Found (ppb)	% Recovery
VAL21	0 (blank)	20.00	93.0	4	0.233	-	0	< 2	< 10	-
VAL22	0 (control)	20.02	86.0	4	0.210	-	0	< 2	< 10	-
VAL23	10	20.00	89.5	4	0.217	7.32	2477	1.58	7.32	73
VAL24	10	20.00	89.5	4	0.217	7.32	2298	1.48	6.82	68
VAL25	10	20.00	88.5	4	0.214	7.30	2383	1.53	7.14	71
VAL26	10	20.02	89.0	4	0.216	7.32	2350	1.51	7.00	70
VAL27	100	20.00	90.0	10	0.087	7.32	11285	6.90	79.2	79
VAL28	100	20.00	88.0	10	0.085	7.31	11727	7.16	84.1	84
VAL29	1000	20.00	90.0	50	0.017	7.32	25412	15.4	885	88
VAL30	1000	20.00	90.0	50	0.017	7.32	24263	14.7	845	85

These values are common to all samples and are used in the calculations detailed in Section II.K.3.0  
 $m = \% \text{ moisture} = 11.7\%$

$V_e$  = volume of extraction solvent used = 100 mL

$V_{std}$  = volume of fortification solution = 1 mL (for fortified samples only)

$V_i$  = HPLC injection volume = 0.05 mL

**Calibration Standards**

Concentration (ng/μL)	Injection Volume (μL)	Amount Inj. (ng)	Retention Time (min)	Peak Area	Amount Found (ng)
0.02	50	1	7.32	1536	1.02
0.05	50	2.5	7.29	4089	2.56
0.1	50	5	7.32	8608	5.28
0.5	50	25	7.31	40240	24.4
1.0	50	50	7.32	83268	50.3

slope = 1658.67  
Y-intercept = -151.4662  
corr. coeff. = 0.99983

\*Note: The displayed calculated values were taken from the Ciba Worksheet program. Attempts to duplicate these calculated values are subject to computer round-off error.

TABLE X. SUMMARY DATA FOR FORTIFIED SOIL: C18 COLUMN

Fortification Level (ppb)	% Recovery for GS-23199		
	California	Georgia	New York
10	78	73	85
10	90	89	69
10	87	120	110
10	90	89	119
10	79	89	
10	82	100	
10		101	
100	100	66	86
100	108	50*	89
100	105	72	
100	92	75	
1000	92	85	90
1000	90	87	82
1000	100	85	
1000	96	88	
Average	92	85 (87)*	91
Standard Deviation	9.1	16.4 (13.8)*	15.9
% Relative Std. Dev.	9.9	19.4 (15.9)*	17.5

\* Uncharacteristic low recovery. This sample was omitted from statistical calculations with an asterisk.

Pooled Recovery Data for all Soils by Fortification Level

	<u>10 ppb</u>	<u>100 ppb</u>	<u>1000 ppb</u>
Average	91	84 (88)*	90
Standard Deviation	14.7	18.5 (14.8)*	5.4
% Relative Std. Dev.	16.1	21.9 (16.8)*	6.1
Range	69 - 120	50 - 108	82 - 100
Number of Samples	17	10	10

TABLE X. SUMMARY DATA FOR FORTIFIED SOIL: C18 COLUMN  
(cont.)

Fortification Level (ppb)	<u>% Recovery for CGA-294849</u>		
	<u>California</u>	<u>Georgia</u>	<u>New York</u>
10	93	103	69
10	82	89	69
10	92	96	79
10	82	89	72
10	77	89	
10	75	84	
10	75	87	
10	89	87	
100	90	90	83
100	79	73	76
100	84	87	
100	84	82	
1000	93	90	86
1000	89	88	86
1000	96	84	
1000	92	87	
Average	86	88	78
Standard Deviation	6.9	6.3	7.1
% Relative Std. Dev.	8.1	7.2	9.2

Pooled Recovery Data for all Soils by Fortification Level

	<u>10 ppb</u>	<u>100 ppb</u>	<u>1000 ppb</u>
Average	84	83	89
Standard Deviation	9.2	5.6	3.7
% Relative Std. Dev.	11.0	6.8	4.1
Range	69 - 103	73 - 90	84 - 96
Number of Samples	20	10	10



TABLE XI. SUMMARY DATA FOR FORTIFIED SOIL: CN COLUMN

<u>Fortification Level (ppb)</u>	<u>% Recovery for GS-23199</u>		
	<u>California</u>	<u>Georgia</u>	<u>New York</u>
10	112	89	154*
10	99	96	50*
10	87	78	108
10		100	120
100	101	72	86
100	96	78	82
1000	117	84	83
1000	111	83	79
Average	103	85	95 (93)*
Standard Deviation	10.5	9.5	31.5 (16.9)*
% Relative Std. Dev.	10.2	11.2	33.0 (18.1)*

\* Statistics with these samples omitted. Uncharacteristic recovery values for unknown reason.

Pooled Recovery Data for all Soils by Fortification Level

	<u>10 ppb</u>	<u>100 ppb</u>	<u>1000 ppb</u>
Average	99 (99)*	86	93
Standard Deviation	26.1 (13.2)*	11.0	16.6
% Relative Std. Dev.	26.3 (13.3)*	12.8	17.8
Range	50 - 154	72 - 101	79 - 117
Number of Samples	11	6	6

TABLE XI. SUMMARY DATA FOR FORTIFIED SOIL: CN COLUMN


Fortification Level (ppb)	<u>% Recovery for CGA-294849</u>		
	<u>California</u>	<u>Georgia</u>	<u>New York</u>
10	62	84	73
10	70	75	68
10	61	46*	71
10		50*	70
100	79	81	79
100	93	70	84
1000	100	74	88
1000	97	71	85
Average	80	76	77
Standard Deviation	16.6	5.6	7.7
% Relative Std. Dev.	20.6	7.4	10.0

\* Samples went to dryness during rotary evaporation resulting in low recoveries. Samples will not be used for statistical calculations.

Pooled Recovery Data for all Soils by Fortification Level


	<u>10 ppb</u>	<u>100 ppb</u>	<u>1000 ppb</u>
Average	70	81	86
Standard Deviation	6.9	7.5	11.8
% Relative Std. Dev.	9.8	9.3	13.7
Range	61 - 84	70 - 93	71 - 100
Number of Samples	9	6	6

FIGURE 1. CHEMICAL NAMES AND STRUCTURES




GS-23199

CAS Name: 6-Methyl-1,2,4-triazine-3,5(2H,4H)-dione  
CAS No.: 932-53-6



CGA-294849

CAS Name: 4-Amino-6-methyl-1,2,4-triazine-3,5(2H,4H)-dione  
CAS No.: 16077-52-4



CGA-215944

CAS Name: (E)-4,5-Dihydro-6-methyl-4-[(3-pyridinyl  
methylene)amino]-1,2,4-triazin-3(2H)-one  
CAS No.: 123312-89-0

FIGURE 2. AG-666 FLOW DIAGRAM FOR SOIL

Weigh 20 gram soil sample. Fortify, if necessary.  
Add 100 mL of 20% water/methanol.  
Reflux extract for one hour.  
Centrifuge and filter sample.  
Measure and record volume of filtered extract.



Remove methanol from sample via rotary evaporation until  
approximately 20 mL remains.  
Basify sample with ammonium hydroxide.  
Pass through SAX SPE.  
Collect non-retained fraction in round bottom flask.



Remove methanol from sample via rotary evaporation until only  
water remains. Do not permit to go dry. Ensure 5-10 mL of  
water is present in sample.  
Acidify sample with formic acid.  
Pass through C18 SPE piggybacked to a SCX SPE.  
Collect non-retained fraction in round bottom flask.




Remove all solvent via rotary evaporation until only 3-5 mL of  
water remains. Do not permit to go dry.  
Add approximately 5 mL of methanol, swirl and transfer to a  
pre-calibrated concentration tube.  
Remove solvent via rotary evaporation until approximately  
2-3 mL of water remains.  
Dilute sample with water to 4-mL mark, or more for samples  
with high level residues.




Analyze by HPLC with UV detection.

FIGURE 3. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS:  
C18 ANALYTICAL COLUMN




Standard: 2 ng Injected




Standard: 5 ng Injected

60

FIGURE 3. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS:  
C18 ANALYTICAL COLUMN (cont.)




Standard: 10 ng Injected



Standard: 50 ng Injected


61

FIGURE 3. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS:  
C18 ANALYTICAL COLUMN (cont.)




Standard: 100 ng Injected

FIGURE 4. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: C18 ANALYTICAL COLUMN



Sample Code: VAL1C, Method Blank




Sample Code: VAL2C, Soil Control


63



FIGURE 4. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: C18 ANALYTICAL COLUMN (cont.)




Sample Code: VAL3C, Soil + 10 ppb




Sample Code: VAL7C, Soil + 100 ppb

64

FIGURE 4. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: C18 ANALYTICAL COLUMN (cont.)




Sample Code: VAL9C, Soil + 1000 ppb




Sample Code: VAL6, Soil + 10 ppb  
Example of late-eluting interference when column has gone bad  
or improperly cleaned between injections.

FIGURE 5. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: C18 ANALYTICAL COLUMN




Sample Code: VAL11B, Method Blank




Sample Code: VAL12B, Soil Control

66

FIGURE 5. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: C18 ANALYTICAL COLUMN (cont.)




Sample Code: VAL13B; Soil + 10 ppb



Sample Code: VAL17B, Soil + 100 ppb

67


FIGURE 5. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: C18 ANALYTICAL COLUMN (cont.)




Sample Code: VAL19B, Soil + 1000 ppb

68

FIGURE 6. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: C18 ANALYTICAL COLUMN




Sample Code: VAL21, Method Blank




Sample Code: VAL22, Soil Control

FIGURE 6. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: C18 ANALYTICAL COLUMN (cont.)




Sample Code: VAL23, Soil + 10 ppb



Sample Code: VAL27, Soil + 100 ppb


FIGURE 6. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: C18 ANALYTICAL COLUMN (cont.)




Sample Code: VAL29, Soil + 1000 ppb



FIGURE 7. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS:  
CN CONFIRMATION COLUMN




Standard: 1 ng Injected




Standard: 2.5 ng Injected

FIGURE 7. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS:  
CN CONFIRMATION COLUMN (cont.)




Standard: 5 ng Injected




Standard: 25 ng Injected

FIGURE 7. TYPICAL CHROMATOGRAMS FOR ANALYTICAL STANDARDS:  
CN CONFIRMATION COLUMN (cont.)




Standard: 50 ng Injected

FIGURE 8. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: CN CONFIRMATION COLUMN




Sample Code: VAL1C, Method Control




Sample Code: VAL2C, Soil Control

FIGURE 8. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: CN CONFIRMATION COLUMN (cont.)




Sample Code: VAL3C, Soil + 10 ppb




Sample Code: VAL7C, Soil + 100 ppb

FIGURE 8. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED CALIFORNIA SOIL: CN CONFIRMATION COLUMN (cont.)




Sample Code: VAL9C, Soil + 1000 ppb

FIGURE 9. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: CN CONFIRMATION COLUMN




Sample Code: VAL11B, Method Control




Sample Code: VAL12B, Soil Control

FIGURE 9. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: CN CONFIRMATION COLUMN (cont.)




Sample Code: VAL13B, Soil + 10 ppb



Sample Code: VAL17B, Soil + 100 ppb




FIGURE 9. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED GEORGIA SOIL: CN CONFIRMATION COLUMN (cont.)




Sample Code: VAL19B, Soil + 1000 ppb

FIGURE 10. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: CN CONFIRMATION COLUMN




Sample Code: VAL21, Method Control




Sample Code: VAL22, Soil Control

84

FIGURE 10. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: CN CONFIRMATION COLUMN (cont.)




Sample Code: VAL25, Soil + 10 ppb



Sample Code: VAL27, Soil + 100 ppb

82

FIGURE 10. TYPICAL CHROMATOGRAMS FOR CONTROL AND FORTIFIED NEW YORK SOIL: CN CONFIRMATION COLUMN (cont.)



Sample Code: VAL29, Soil + 1000 ppb

83

IX. REFERENCES

- (1) Vargo, J. D., Ciba Protocol 376-96, "Validation of "Draft" Analytical Method AG-666 for the Determination of GS-23199 and CGA-294849, Metabolites of CGA-215944, by High Performance Liquid Chromatography with UV Detection," including Protocol Amendment 1.
- (2) Vargo, J. D., Ciba Residue Test Report RI-MV-016-96, Report Number 1

84