

Cover Sheet for

ENVIRONMENTAL CHEMISTRY METHOD

Pesticide Name: Metolachlor

MRID #: 454996-01

Matrix: Water

Analysis: Immunoassay

This method is provided to you by the Environmental Protection Agency's (EPA) Environmental Chemistry Laboratory (ECL). This method *is not* an EPA method but one which was submitted to EPA by the pesticide manufacturer to support product registration. EPA recognizes that the methods may be of some utility to state, tribal, and local authorities, but makes no claim of validity by posting these methods. Although the Agency reviews *all* Environmental Chemistry Methods submitted in support of pesticide registration, the ECL evaluates only about 30% of the currently available methods. Most methods perform satisfactorily but some, particularly the older methods, have deficiencies. Moreover, the print quality of the methods varies considerably because the methods originate from different sources. Therefore, the methods offered represent the best available copies.

If you have difficulties in downloading the method, or further questions concerning the methods, you may contact Elizabeth Flynt at 228-688-2410 or via e-mail at flynt.elizabeth@epa.gov.

VOLUME 2 OF 14 OF SUBMISSION

CGA-77102: METHOD

TITLE

ANALYTICAL METHOD FOR THE SEMI-QUANTITATIVE DETERMINATION OF
CGA-77102 IN WATER BY ENZYME IMMUNOASSAY INCLUDING VALIDATION
DATA

DATA REQUIREMENT

Not Applicable

AUTHORS

Jody L. Needham, B.A.
James F. Brady, Ph.D.

COMPLETION DATE

June 23, 1999

PERFORMING LABORATORY

Analytical Resources Group
Development Resources/Chemical Support Department
Novartis Crop Protection, Inc.
Greensboro, NC 27419

LABORATORY STUDY IDENTIFICATION


Novartis Number 1004-98

SUBMITTER/SPONSOR

Novartis Crop Protection, Inc.
410 Swing Road
Post Office Box 18300
Greensboro, NC 27419

VOLUME 1 OF 1 OF STUDY

PAGE 1 OF 40



STATEMENT OF NO DATA CONFIDENTIALITY CLAIM

No claim of confidentiality is made for any information contained in this study on the basis of its falling within the scope of FIFRA Section 10 (d) (1) (A), (B) or (C).

Company: Novartis Crop Protection, Inc.


Company Representative: Karen Stumpf

Title: Senior Regulatory Manager

Signature: *Karen Stumpf* Date: 6/23/99

These data are the property of Novartis Crop Protection, Inc. and, as such, are considered to be confidential for all purposes other than compliance with the regulations implementing FIFRA Section 10.

Submission of these data in compliance with FIFRA does not constitute a waiver of any right to confidentiality which may exist under any other provision of common law or statute or in any other country.



GOOD LABORATORY PRACTICE COMPLIANCE STATEMENT

This study was performed in accordance with Good Laboratory Practice Standards as required by the EPA-FIFRA Good Laboratory Practice Standards, 40 CFR part 160 with the following exception(s):

1. Alachlor and Dimethanamid-ESA analytical standards did not comply with the above guidelines.

Jody L. Needham
Jody L. Needham, B.A.
Study Director
Development Resources/
Chemical Support Department

6/23/99
Date

m. w. Cheung
Max W. Cheung, Ph.D.
Manager, Development Resources/
Chemical Support Department
Representative of Submitter/Sponsor

06/23/99
Date

Submitter/Sponsor: Novartis Crop Protection, Inc.
410 Swing Road
Post Office Box 18300
Greensboro, NC 27419

REPORT APPROVAL

Jody L. Needham
Jody L. Needham, B.A.
Study Director
Development Resources/
Chemical Support Department

6/23/99
Date

m.w. Cheung
Max W. Cheung, Ph.D.
Manager, Development Resources/
Chemical Support Department
Representative of Submitter and Sponsor

06/23/99
Date

Submitter/Sponsor: Novartis Crop Protection, Inc.
410 Swing Road
Post Office Box 18300
Greensboro, NC 27419

CERTIFICATION

The analytical reports and experimental results included in Project H077102NAU950A, in support of analytical method 1004-98, are authentic accounts of the experiments.

Jody L. Needham
Jody L. Needham, B.A.
Study Director
Development Resources/
Chemical Support Department

6/23/99
Date

James F. Brady
James F. Brady, Ph.D.
Author
Development Resources/
Chemical Support Department

6-23-99
Date

Richard H. Ross, Jr.
Richard H. Ross, Jr., Ph.D.
Director, Development Resources/
Chemical Support Department
336-632-2384

6/23/99
Date

QUALITY ASSURANCE STATEMENT

Study Title: Analytical Method for the Semi-Quantitative Determination of CGA-77102 in Water by Enzyme Immunoassay Including Validation Data

Study Director: Needham, Jody L.

Method Number: 1004-98

Pursuant to Good Laboratory Practice Regulations, this statement verifies that the aforementioned study was inspected and/or audited and the findings reported to Management and to the Study Director by the Quality Assurance Unit on the dates listed below.

<u>Inspection/Audit Type</u>	<u>Inspection/Audit Dates</u>	<u>Reporting Date</u>
Audit Protocol	11/18/98 - 11/18/98	11/19/98
Inspect Analytical	05/05/99 - 05/05/99	05/06/99
Audit Final Report	05/18/99 - 05/20/99	05/24/99

Prepared by: Kathleen M. Campbell Date: 5-24-99

Kathleen M. Campbell
Staff Quality Assurance Auditor
Quality Assurance Unit
Novartis Crop Protection, Inc.

TABLE OF CONTENTS

	<u>Page No.</u>
TITLE PAGE	1
STATEMENT OF NO DATA CONFIDENTIALITY CLAIM	2
GOOD LABORATORY PRACTICE COMPLIANCE STATEMENT	3
REPORT APPROVAL	4
CERTIFICATION	5
QUALITY ASSURANCE STATEMENT	6
TABLE OF CONTENTS	7
GENERAL INFORMATION	10
I. SUMMARY AND INTRODUCTION	12
A. SCOPE	12
B. PRINCIPLE	12
II. MATERIALS AND METHODS	12
A. APPARATUS	12
B. REAGENTS	13
C. ANALYTICAL PROCEDURE	14
1.0 Sample Preparation	14
2.0 Enzyme Immunoassay	15
2.1 Inhibition of Enzyme Conjugate	15
2.2 Wash	16
2.3 Color Development	17

TABLE OF CONTENTS
(Continued)

	<u>Page No.</u>
D. INSTRUMENTATION	18
1.0 Description and Operating Conditions.....	18
2.0 Standardization	18
E. INTERFERENCES.....	18
F. CONFIRMATORY TECHNIQUES	19
G. TIME REQUIRED	19
H. MODIFICATIONS AND POTENTIAL PROBLEMS	19
I. PREPARATION OF STANDARD SOLUTIONS.....	20
J. METHODS OF CALCULATION	20
III. RESULTS AND DISCUSSION.....	20
A. VALIDATION.....	20
1.0. ACCURACY AND PRECISION.....	21
2.0. SPECIFICITY	21
3.0. LIMITS OF DETECTION AND QUANTITATION.....	21
B. LIMITATIONS	22
C. CIRCUMSTANCES AFFECTING THE STUDY.....	22
IV. CONCLUSION	22
V. TABLES AND FIGURES	
TABLE 1. RESULTS OF STANDARD ADDITION EXPERIMENTS.....	23

TABLE OF CONTENTS
(Continued)

	<u>Page No.</u>
TABLE 2. CROSS-REACTIVITY RESULTS FOR THE CGA-77102 ENZYME IMMUNOASSAY.....	24
FIGURE 1. FLOW DIAGRAM FOR METHOD 1004-98	25
FIGURE 2. PLATE LAYOUT CHART ON WHICH THE ANALYST WILL INDICATE THE POSITION OF SAMPLE AND STANDARD SOLUTIONS IN THE MICROWELL ASSAY PLATE	26
FIGURE 3. A TYPICAL STANDARD CURVE OF THE CGA-77102 MICROWELL ASSAY.....	27
FIGURE 4. STRUCTURES OF TEST SUBSTANCES EVALUATED FOR CROSS-REACTIVITY.....	28
VI. REFERENCE.....	32
APPENDIX 1 RESIDUE TEST REPORT NO. 1	33
APPENDIX 2. RESIDUE TEST REPORT NO. 2.....	37

GENERAL INFORMATION

Protocol Number: 1004-98 plus one amendment

Study Initiation Date: December 9, 1998

Study Personnel:

Management: Max W. Cheung, Ph.D.
Manager, Analytical Resources Group

Study Director and Analyst: Jody L. Needham, B.A.
Chemist I, Analytical Resources Group

Analyst: James F. Brady, Ph.D.
Scientist III, Analytical Resources Group

Reference Substance Information:

Novartis Code Number	CAS Number	CAS Name	Lot Number	Purity	Reassay Date
Metolachlor	51218-45-2	Acetamide, 2-chloro-N-(2-ethyl-6-methylphenyl)-N-(2-methoxy-1-methylethyl)-	S94-1719	99.3 %	2/2000
Alachlor	15972-60-8	Acetamide, 2-chloro-N-(2,6-diethylphenyl)-N-(methoxymethyl)-	S86-0926	98.4 %	NON-GLP
Alachlor-ESA (NOA-436644)	142363-53-9	Ethanesulfonic acid, 2-[(2,6-diethylphenyl)(methoxymethyl)amino]-2-oxo-	DAH-XXII-24	99.4 %	11/2000
Acetochlor	34256-82-1	Acetamide, 2-chloro-N-(ethoxymethyl)-N-(2-ethyl-6-methylphenyl)-	S97-2120	98 %	8/99
Acetochlor-ESA (NOA-436645)	187022-11-3	Ethanesulfonic acid, 2-[(ethoxymethyl)(2-ethyl-6-methylphenyl)amino]-2-oxo-	DAH-XXII-30	97.8 %	11/2000
Dimethanamid	87674-68-8	Acetamide, 2-chloro-N-(2,4-dimethyl-3-thienyl)-N-(2-methoxy-1-methylethyl)-	S97-2034	99.7 %	7/2000
Dimethanamid-ESA	205939-58-8	Ethanesulfonic acid, 2-[(2,4-dimethyl-3-thienyl)(2-methoxy-1-methylethyl)amino]-2-oxo-	BPM-XVIII-19	Approx 87.5 %	NON-GLP
Metalaxyl	57837-19-1	DL-Alanine, N-[2,6-dimethylphenyl)-N-(methoxyacetyl)-, methyl ester	S94-1763	97.1 %	10/2000

GENERAL INFORMATION
(continued)

Novartis Code Number	CAS Number	CAS Name	Lot Number	Purity	Reassay Date
CGA-354743	Not Assigned	Ethanesulfonic acid, 2-[(2-ethyl-6-methylphenyl)(2-methoxy-1-methylethyl)amino]-2-oxo-, sodium salt	DAH-XXVI-12	96.8 %	9/2000
CGA-50720	152019-74-4	Acetic acid, [(2-ethyl-6-methylphenyl)amino]oxo-	DPS-I-74-2	>99.9%	6/99
CGA-40172	131068-72-9	Acetamide, N-(2-ethyl-6-methylphenyl)-2-hydroxy-N-(2-methoxy-1-methylethyl)-	JAK-II-42(IV)	95.9 %	6/99
CGA-67869	66637-79-4	DL-Alanine, N-(2,6-dimethylphenyl)-N-(hydroxyacetyl)-, methyl ester	DAH-XV-35	>99.9%	6/99
CGA-67125	Not Assigned	Formamide, N-(2-ethyl-6-methylphenyl)-N-(2-methoxy-1-methylethyl)-	DPS-VII-41-1	99.9 %	12/99
CGA-37735	97055-05-5	Acetamide, N-(2-ethyl-6-methylphenyl)-2-hydroxy-	JAK-II-42(II)	99.2 %	3/2000
CGA-41638	65513-61-3	Acetamide, 2-chloro-N-(2-ethyl-6-methylphenyl)-N-(2-hydroxy-1-methylethyl)-	S94-1765	95 %	9/99
CGA-51202	152019-73-3	Acetic acid, [(2-ethyl-6-methylphenyl)(2-methoxy-1-methylethyl)amino]oxo-	JAK-I-63	98.8 %	3/2000
CGA-62826	87764-37-2	DL-Alanine, N-(2,6-dimethylphenyl)-N-(methoxyacetyl)-	BPM-I-48	98.0 %	11/2000
CGA-46576	Not Assigned	Not Assigned	AW-I-43	97.4 %	6/99
CGA-380168	Not Assigned	Ethanesulfonic acid, 2-[(2-ethyl-6-methylphenyl)(2-methoxy-1-methylethyl)amino]-2-oxo-, (S)-, sodium salt	DAH-XXIII-62	98.8 %	1/99
CGA-77102	87392-12-9	Acetamide, 2-chloro-N-(2-ethyl-6-methylphenyl)-N-(2-methoxy-1-methylethyl)-, (S)-	S95-1844	98.5 %	11/2000
CGA-37913	61520-53-4	1-Propanol, 2-[(2-ethyl-6-methylphenyl)amino]-	JAK-XVI-38	99.0 %	12/2000

I. SUMMARY AND INTRODUCTION

A. SCOPE

This method is to be used for the semiquantitative determination of CGA-77102 ([Acetamide, 2-chloro-*N*-(2-ethyl-6-methylphenyl)-*N*-(2-methoxy-1-methylethyl)-, (*S*)-], CAS # 87392-12-9) in water. The limit of detection (LOD, the smallest dose that yields a response that is statistically significant different than the response of the zero dose) is 0.05 ppb of CGA-77102. The limit of quantitation (LOQ, the lowest level of fortification with which an acceptable recovery can be obtained) for the method is 0.10 ppb of CGA-77102.

B. PRINCIPLE

A 100- μ l aliquot of a representative water sample is added to a microwell coated with CGA-77102 antibody. The assay is carried out by sequential addition of enzyme conjugate, wash solution, enzyme substrate, and stop solution. The reaction is terminated by acidification. Quantification is performed spectrophotometrically at 450 nm. A flow diagram for the method is presented in Figure 1.

II. MATERIALS AND METHODS

A. APPARATUS

- 1.0 CGA-77102 Plate kit, enzyme immunoassay microwell assay, Beacon Analytical Systems, Scarborough, ME, catalog #CPP-005.
- 2.0 Microtiter plate reader, Titertek Multiskan MCC/340 MK II, ICN catalog #78-626-00, or equivalent.
- 3.0 AUTOmate™ microplate software and manual, ICN cat. #78-599-02 and 78-599-98, respectively, or equivalent.
- 4.0 Microtiter plate shaker, Lab-Line Instruments catalog #4625 or equivalent.
- 5.0 ELP-40 Microplate Strip Washer, Bio-Tek Instruments, Bio-Tek Part No: 4091000, or equivalent (use of equipment is optional).
- 6.0 Eight-channel pipette capable of dispensing 100- μ l volumes, Elkay Labsystems catalog #4142-417 or equivalent.

- 7.0 Adjustable microliter pipette, capable of dispensing 20- to 200- μ l volumes, Rainin cat. #P-200 or equivalent.
- 8.0 Pipette tips for eight-channel pipette for volumes of 20- to 200- μ l, packaged in an eight by twelve array, Costar cat. #4865 or equivalent.
- 9.0 Pipette tips for single-channel pipette for volumes of 100- to 1000- μ l, packaged as 100 tips per rack, Costar cat. #4867 or equivalent.
- 10.0 Graduated pipette, 5- and 10-mL, Fisher cat. # 13-660E and 13-660F or equivalent.
- 11.0 Microtiter assay plate, 96 well, flat bottom, polystyrene, non-sterile, not treated for cell culture work, Flow Laboratories cat. #76-102-05 or equivalent.
- 12.0 Reagent reservoirs for multichannel pipettes, non-sterile, ICN Cat. #77-824-01 or equivalent.
- 13.0 pH meter, Corning M-90 pH Stick, VWR cat. #34102-081 or equivalent.
- 14.0 Magnetic stir bar, 0.5-1.0 inch in length, Fisher cat. #14-511-61 or 14-511-63 or equivalent.
- 15.0 Magnetic stirrer, Fisher cat. #11-496-21 or equivalent.
- 16.0 4-oz Amber Boston round bottles, Fisher cat. #03-320-4B or equivalent.

B. REAGENTS

- 1.0 Distilled, deionized water (H_2O).
- 2.0 Acetonitrile (ACN), HPLC Grade, Fisher cat. #A998SK-4 or equivalent.
- 3.0 Sodium Chloride (NaCl) crystal, certified A.C.S., CAS No. 7647-14-5, Fisher cat. # S271-500 or equivalent.
- 4.0 Polyoxyethylene-sorbitan monolaurate (TWEEN 20). Sigma cat. #P-1379 or equivalent.

- 5.0 50% H₂O/TWEEN 20 – Combined 1 mL of TWEEN 20 and 1 mL of H₂O.
- 6.0 H₂O/TWEEN 20 wash solution – Dissolve 9 g of NaCl and 1 mL of 50% H₂O/TWEEN 20 in 1 liter of H₂O.
- 7.0 Enzyme conjugate, included in the CGA-77102 plate kit, Beacon Analytical Systems.
- 8.0 Substrate solution, included in the CGA-77102 plate kit, Beacon Analytical Systems. Alternatively, a commercially available horseradish peroxidase tetramethylene benzidine substrate, such as “K-Blue” (Neogen Corporation cat. #300177) may be used.
- 9.0 Stop solution, included in the CGA-77102 plate kit, Beacon Analytical Systems. Alternatively, 1.0 N HCl, prepared by the analyst, may be used.
- 10.0 CGA-77102 analytical standard, Novartis Crop Protection, P.O. Box 18300, Greensboro, NC 27419-8300. Storage conditions: frozen.

C. ANALYTICAL PROCEDURE

The antibody-coated microwells, all reagents and sample and standard solutions must be warmed to room temperature, approximately 22°C, prior to use. The reaction kinetics are temperature dependent. Be certain all solutions are warmed to room temperature before running the assay.

1.0 Sample Preparation

1.1 Preparation of Water Samples

- 1.1.1 Water samples are received and stored refrigerated until use.
- 1.1.2 The pH of each sample should be measured. Set the sample container on a magnetic stirrer and add a small stir bar. Allow the sample to stir for a few seconds prior to inserting the electrode. Record the value obtained. Samples having a pH value between 5.0 and 9.0 may be analyzed by this method. Samples having

a pH outside this range are unsuitable for this assay.

2.0 Enzyme Immunoassay

2.1 Inhibition of Enzyme Conjugate

- 2.1.1 Approximately 200 μ l of each sample or standard solution is added individually to wells of the uncoated reservoir plate as indicated on the plate layout chart previously completed by the analyst. A sample chart is illustrated in Figure 2.
- 2.1.2 Decant the enzyme conjugate solution into a multichannel pipette reagent reservoir.
- 2.1.3 Adjust an eight-channel pipette to deliver 100 μ l (refer to the pipette operation manual or applicable Novartis Crop Protection SOP for operation of the multichannel pipette). Snugly fit eight tips to the pipette and transfer a 100- μ l aliquot from column 1, rows A through H, of the reservoir plate to the corresponding column on the antibody-coated assay plate. Eject the used tips to a waste receptacle. Attach eight new tips and proceed to transfer aliquots from columns 2 through 12 of the reservoir plate to the assay plate taking care to change tips after each transfer.

If an alternative, equivalent pipetting device is used in lieu of a hand-held multi-channel pipette, the 100- μ l aliquot from each well of the reservoir plate is transferred to the assay plate in a manner consistent with the operation of the device.

- 2.1.4 As soon as all sample and standard solutions have been transferred to the assay plate, successively pipette 100- μ l aliquots of the enzyme conjugate solution from the reagent reservoir to each column of the assay plate from left to right as

described above. The pipette tips must be changed after each transfer.

Note: The enzyme conjugate solution has a surfactant-like character and tends to bubble and foam when pipetted. Through preliminary trials the analyst should ensure that manual or automated pipetting devices reproducibly deliver 100- μ l aliquots without bubbling or foaming.

- 2.1.5 Place the plate on the plate shaker. Start the shaker with the power control set to the "constant" position. The speed setting should be adjusted to approximately ninety oscillations per minute or slightly less than "2" on the Lab Line shaker. Allow the plate to shake for one hour.

Note: The reaction kinetics of this assay are temperature dependent. To avoid exposing the plate to harsh fluctuations in temperature such as strong drafts, the plate shaker should be placed inside an opaque closed chamber. Placing a cardboard box over the shaker will suffice if a more elaborate apparatus, such as a darkened Plexiglas chamber, is not available.

2.2 Wash

- 2.2.1 Remove the plate from the shaker. Using a hand-held multichannel pipette, remove the reactants from each column left to right across the plate.
- 2.2.2 Starting with column one, add approximately 200 μ l of wash solution to the each well, gently agitate the plate and immediately remove the wash solution. Proceed across the plate from column two to twelve washing each column in a similar fashion. Wash the plate two additional times.

2.2.3 After the final wash, shake the plate vigorously to remove most of the remaining liquid. The plate may be inverted and blotted with dry paper towels. Blot only the exterior of the plate. Do not insert the blotting paper into the wells. Visually examine the plate to ensure little wash remains.

Should only part of the microtiter plate be used, the analyst should exercise care not to contaminate unused wells.

If an alternative washing device, such as the ELP-40 microplate strip washer, is used, be certain the reactants are removed in a manner consistent with their addition.

2.3 Color Development

2.3.1 Add 100 μ l of substrate to each well across the plate in the same fashion as the previous solutions were dispensed. Tips need not be changed between additions.

2.3.2 Place the plate on the plate shaker and shake for 0.5 hour with the controls set at the same settings as described in Section II.C.2.1.5. Individual wells will gradually turn varying shades of blue.

2.3.3 Terminate color development by adding 100 μ l of stop solution to each well in the same fashion as the color reagent was added. Mix the acidified solution well by repeated pipetting and dispensing of each well's contents, taking care not to form bubbles or to leave liquid behind in the pipette tips. Should bubbles form, they can be broken by manipulation with a pipette tip. Tapping the side of the strip holder gently may aid mixing of an acidified solution. A well-mixed solution will appear yellow to the eye with no traces of blue remaining. Be sure to change tips

between each addition to avoid contaminating a column of wells with residue from a previous column.

The absorbance of the final reaction product within each microwell should be analyzed within 30 minutes of addition of stop solution.

D. INSTRUMENTATION

1.0 Description and Operating Conditions

1.1 The Titertek Multiskan MCC/340 MK II eight-channel filter spectrophotometer is used to measure the absorbance of the final reaction solutions in each well. The operator should select filter 4 (450 nm) for this method. The instrument should be warmed up for ten minutes prior to use. Refer to the instrument operating instructions for further details on the operation of this instrument.

2.0 Standardization

- 2.1 Each assay consists of standards and samples run concurrently on the same plate. The absorbance values obtained from sample solutions may only be compared to the absorbance of standards run in the same analytical set. CGA-77102 standards range from 4.0 to 0.05 ng/mL in addition to a H₂O blank (zero dose).
- 2.2 Spectrophotometric analysis of the colored reaction products obtained from standard solutions will yield absorbance values measured at 450 nm (A_{450}). With a calculator or computer, use these data to generate a log/linear regression function. This curve is in the form of $y = m \log(x) + b$. The concentrations of the CGA-77102 standards are plotted on a logarithmic scale on the horizontal (x) axis and the absorbance values are plotted on a linear scale on the vertical (y) axis. An example of a typical standard curve is shown in Figure 3.

E. INTERFERENCES

- 1.0 The antibodies used in this assay bind primarily to CGA-77102 and metolachlor. Refer to Section III.A.2.0.

SPECIFICITY for a detailed discussion of the cross-reactivity parameters of this assay.

- 2.0 Slight inhibition was observed when glassware was not rinsed before use. Therefore, to obtain best results, all glassware must be pre-rinsed before use. This includes volumetric flasks, bottles, and scintillation vials used in the preparation of standards, beakers used to contain distilled water, and any other glassware that may be used. Pre-rinsing involves triple rinsing with Acetone and blowing dry with Nitrogen air. Vessels used to contain water should then be rinsed three times with de-ionized water before filling with distilled water.

F. CONFIRMATORY TECHNIQUES

- 1.0 The ability of the microwell assay to respond to CGA-77102 in water may be assessed by the analysis of samples fortified with a known amount of CGA-77102.

G. TIME REQUIRED

An analyst can analyze forty-two samples and six standards in duplicate in approximately three hours.

H. MODIFICATIONS AND POTENTIAL PROBLEMS

- 1.0 As previously stated in the introduction to the ANALYTICAL PROCEDURE, the microwell assay plates, reagents and sample solutions must be warmed to room temperature, approximately 22°C, before use. Placing the plate (in its plastic bag) and the reagents in a hood with the sash drawn low for approximately 30 minutes will bring the reagents to temperature efficiently.
- 2.0 The pH of the water sample, if extreme, may affect antibody binding. The pH of each sample should be examined. If a sample is found to have a pH outside the range of 5.0 - 9.0, it is unsuitable for analysis by this method.
- 3.0 The bottom surface of the microwell strips is the optical surface through which the absorbance of the final reaction product will be measured. The analyst should exercise care to prevent damage to this surface. The strip holder is designed to prevent the bottom of the strips from contacting flat surfaces, such as a bench top, upon which the strip

holder may be placed. Nevertheless, the analyst should maintain a clean work area as a preventative measure.

- 4.0 The analyst should take care to be certain all sample and standard solutions are positioned properly in the reservoir plate. The large sample load requires that special attention to detail be maintained throughout the analysis.
- 5.0 The analyst may observe absorbance values which are less than those of the highest standard, 4.0 ng/mL. In this event, the concentration of the corresponding sample cannot be calculated since its absorbance readings do not fall within the range of the standard curve. The sample should be diluted and re-assayed to obtain absorbance readings which lie within the bounds of the standard curve. The concentration of the undiluted sample can be calculated by multiplying the concentration of the diluted sample by the dilution factor.
- 6.0 The transfer of solutions by multichannel pipettes requires the analyst to constantly monitor his or her technique. Pipetting errors are the major source of error in immunoassay methodology.

I. PREPARATION OF STANDARD SOLUTIONS

Dissolve 5.0 mg of CGA-77102 in a minimal amount of ACN (approximately 15 mL), sonicate, and bring volume up to 100 mL with H₂O. This will make 100 mL of a 50 µg/mL solution. Serially dilute this solution with H₂O to make 5000-, 1000-, 10-, 4.0-, 1.0-, 0.3-, 0.1-, and 0.05-ng/mL standards. Prepare a blank consisting solely of H₂O. All standard solutions should be stored in amber Boston round bottles. Store these solutions at 4°C when not in use.

J. METHODS OF CALCULATION

CGA-77102 residues in sample solutions are determined by inserting the absorbance value of a given sample into the log/linear regression function generated by methods described in Section II.D.2.2. These calculations may be made on a computer or hand-held calculator.

III. RESULTS AND DISCUSSION

A. VALIDATION

The objective of this study was to validate "draft" analytical method 1004-98 by analysis of freshly fortified water samples.

1.0 ACCURACY AND PRECISION

The accuracy and precision of this method were assessed by conducting standard addition experiments using laboratory reagent water. Two groups of ten samples each were fortified with 0.10 and 1.0 ppb of CGA-77102, respectively, and analyzed pre- and post-fortification. The results of these experiments are shown in Table 1. The mean percent recovered (95% at 0.10 ppb and 123% at 1.0 ppb) indicate the method has acceptable accuracy. The standard deviations of these measurements (0.01 ppb at 0.10 ppb and 0.09 ppb at 1.0 ppb) indicate Method 1004-98 has acceptable precision.

2.0. SPECIFICITY

The specificity of Method 1004-98 was evaluated by determining the cross reactivity parameters of twenty-one structurally related test substances. Results of these experiments are shown in Table 2. These data indicate the antibodies used in Method 1004-98 are primarily reactive to CGA-77102. Metolachlor (CGA-24705) is strongly recognized (55%) because one-half of this material is, in fact, CGA-77102 (S-CGA-24705). Dimethanamid and CGA-380168, the ethanesulfonic acid metabolite of CGA-77102, are far less reactive (12.9 and 2.1%, respectively). The remaining test substances are less than 1% reactive or not reactive at all. Thus, this method has little or no reactivity with degradates of chloroacetanilide compounds. Refer to Figure 4 for the chemical structures of the test substances screened for cross reactivity.

3.0. LIMITS OF DETECTION AND QUANTITATION

The limits of detection for test substances screened were calculated using Brady's modification of Rodbard's method (1). This method was found to have an LOD of 0.05 ppb for CGA-77102, 0.05 ppb for Metolachlor, 0.20 ppb for Dimethanamid, and 0.23 ppb for CGA-380168. All other test substances were found to have less than one percent reactivity to the assay relative to CGA-77102 so additional limits of detection were not determined. An LOQ was determined only for CGA-77102 because this was the only test substance used to conduct standard addition experiments. Based on the data generated in these

experiments, CGA-77102 was determined to have an LOQ of 0.10 ppb.

B. LIMITATIONS

- 1.0 This method is intended to complement, not replace, chromatographic analyses. It can be used as an inexpensive screening technique to remove samples that yield responses below a pre-selected level of concern from further, more expensive analytical procedures.
- 2.0 As stated in Section II.H.2.0, the pH of a sample may affect antibody binding. Consequently, the pH of all samples must be checked prior to analysis. If the pH of a sample is outside the range of 5.0 – 9.0, that sample is unsuitable for analysis by Method 1004-98.
- 3.0 The cross reactivity data presented in Table 2 is not meant to be an exhaustive evaluation of the antibodies used in Method 1004-98. The twenty-one test substances screened represent a limited attempt to investigate the cross reactivity of the assay. Therefore, analysts should confirm residues of concern determined by this method by an additional analytical method.

C. CIRCUMSTANCES AFFECTING THE STUDY

There were no adverse circumstances affecting the quality or integrity of the data.

IV. CONCLUSION

This method has been found to be a rapid, reliable and accurate method for the semi-quantitative determination of CGA-77102 in water.

V. TABLES AND FIGURES

TABLE 1. RESULTS OF STANDARD ADDITION EXPERIMENTS

Distilled Water Samples	ppb CGA-77102 Equivalents Found				² Net	CGA-77102 Fortification Level (ppb)
	Pre-fortification		Post-fortification			
A	<0.05	<0.05	0.09	0.09	0.09	0.10
B	<0.05	<0.05	0.08	0.09	0.08	0.10
C	<0.05	<0.05	0.10	0.09	0.10	0.10
D	<0.05	<0.05	0.09	0.09	0.09	0.10
E	<0.05	<0.05	0.10	0.16	0.13	0.10
F	<0.05	<0.05	0.09	0.10	0.10	0.10
G	<0.05	<0.05	0.08	0.10	0.09	0.10
H	<0.05	<0.05	0.10	0.10	0.10	0.10
I	<0.05	<0.05	0.09	0.10	0.09	0.10
J	<0.05	<0.05	0.08	0.09	0.09	0.10

Mean ± SD = 0.10 ± 0.01

N = 10

% CV = 13.4

Distilled Water Samples	ppb CGA-77102 Equivalents Found				² Net	CGA-77102 Fortification Level (ppb)
	Pre-fortification		Post-fortification			
A	<0.05	<0.05	1.03	1.22	1.12	1.0
B	<0.05	<0.05	1.10	1.23	1.17	1.0
C	<0.05	<0.05	1.36	1.39	1.37	1.0
D	<0.05	<0.05	1.24	1.24	1.24	1.0
E	<0.05	<0.05	1.21	1.28	1.25	1.0
F	<0.05	<0.05	1.34	1.37	1.36	1.0
G	<0.05	<0.05	1.15	1.32	1.24	1.0
H	<0.05	<0.05	1.16	1.25	1.20	1.0
I	<0.05	<0.05	1.21	1.28	1.25	1.0
J	<0.05	<0.05	1.09	1.08	1.09	1.0

Mean ± SD = 1.23 ± 0.09

N = 10

% CV = 7.3

¹Refer to Appendix 2.

²Net values were calculated by subtracting the mean pre-fortification result from the mean post-fortification result.

TABLE 2. ¹CROSS-REACTIVITY RESULTS FOR THE CGA-77102 ENZYME IMMUNOASSAY

TEST SUBSTANCE	LOD	² LOQ	³ <i>I</i> ₅₀	PERCENT REACTIVITY RELATIVE TO CGA-77102	SET #	ANALYSIS DATES	NUMBER OF ANALYSES
⁴ CGA-77102	0.05	0.10	0.48	100	990104A 990105A 990107A 990128A	1/4/99 1/5/99 1/7/99 1/28/99	24
Metolachlor	0.05	--	1.0	55	990104A	1/4/99	6
Dimethanamid	0.20	--	4.2	12.9	990104A	1/4/99	6
CGA-380168	0.23	--	13.5	2.1	990105A	1/5/99	6
CGA-41638	--	--	46	<1.0	990105A	1/5/99	6
Acetochlor	--	--	68	<1.0	990104A	1/4/99	6
CGA-40172	--	--	84	<1.0	990107A	1/7/99	6
Alachlor	--	--	470	<1.0	990104A	1/4/99	6
Alachlor-ESA	--	--	⁵ NR	NR	990104A	1/4/99	6
Acetochlor-ESA	--	--	NR	NR	990104A	1/4/99	6
Metalaxyl	--	--	NR	NR	990107A	1/7/99	6
CGA-354743	--	--	NR	NR	990107A	1/7/99	6
CGA-50720	--	--	NR	NR	990107A	1/7/99	6
CGA-67869	--	--	NR	NR	990107A	1/7/99	6
CGA-67125	--	--	NR	NR	990107A	1/7/99	6
CGA-37735	--	--	NR	NR	990104A	1/4/99	6
CGA-51202	--	--	NR	NR	990105A	1/5/99	6
CGA-62826	--	--	NR	NR	990105A	1/5/99	6
CGA-46576	--	--	NR	NR	990105A	1/5/99	6
Dimethanamid-ESA	--	--	NR	NR	990128A	1/28/99	6
CGA-37913	--	--	NR	NR	990128A	1/28/99	6

¹ Refer to Appendix 1.

² Limit of Quantitation (LOQ), the lowest level of fortification (in units of ng/mL) with which an acceptable recovery can be obtained. An LOQ was determined for CGA-77102 only.

³ *I*₅₀, the concentration of test substance (in units of ng/mL) that yields half the response of the zero dose.

⁴ Parameters shown for CGA-77102 are the average of four determinations.

⁵ NR - not reactive.

FIGURE 1. FLOW DIAGRAM FOR METHOD 1004-98




FIGURE 2. PLATE LAYOUT CHART ON WHICH THE ANALYST WILL INDICATE THE POSITION OF SAMPLE AND STANDARD SOLUTIONS IN THE MICROWELL ASSAY PLATE

PLATE LAYOUT SHEET

	1	2	3	4	5	6	7	8	9	10	11	12	
A													A
B													B
C													C
D													D
E													E
F													F
G													G
H													H
	1	2	3	4	5	6	7	8	9	10	11	12	

Plate ID
 Date
 Analyst
 Notebook/Page

Comments

FIGURE 3. A TYPICAL STANDARD CURVE OF THE CGA-77102 MICROWELL ASSAY







FIGURE 4. STRUCTURES OF TEST SUBSTANCES EVALUATED FOR CROSS-REACTIVITY




CGA-77102
(S-CGA-24705)




Metolachlor
(CGA-24705)




Alachlor



Alachlor-ESA




Acetochlor




Acetochlor-ESA


FIGURE 4. STRUCTURES OF TEST SUBSTANCES EVALUATED FOR CROSS-REACTIVITY (continued)




Dimethanamid




Dimethanamid-ESA




Metalaxyl



CGA-354743




CGA-50720




CGA-40172


FIGURE 4. STRUCTURES OF TEST SUBSTANCES EVALUATED FOR CROSS-REACTIVITY (continued)




CGA-67869




CGA-67125




CGA-37735



CGA-41638




CGA-51202




CGA-62826


FIGURE 4. STRUCTURES OF TEST SUBSTANCES EVALUATED FOR CROSS-REACTIVITY (continued)



CGA-46576



CGA-380168
(S-CGA-354743)



CGA-37913

VI. REFERENCE

1. Brady, J. F. Interpretation of immunoassay data. In *Immunoanalysis of Agrochemicals: Emerging Technologies*. Karu, A., Nelson, J. Wong, R., Eds.; ACS Symposium Series 586; American Chemical Society: Washington, DC, 1995; pp 266-287.

APPENDIX 1

RESIDUE TEST REPORT NO. 1

Development Resources/Chemical Support Department
Novartis Crop Protection, Inc.
Greensboro, North Carolina

RESIDUE TEST REPORT
TEST NUMBER: RI-MV-002-99
INVENTORY NUMBER: 18425.1
REPORT NO.: 1

PROJECT NUMBER: H077102NAU950A

PROTOCOL NUMBER: 1004-98

TEST SUBSTANCES: CGA-77102, Metolachlor, Alachlor, Alachlor-ESA, Acetochlor, Acetochlor-ESA, Dimethanamid, Dimethanamid-ESA, Metalaxyl, CGA-354743, CGA-50720, CGA-40172, CGA-67869, CGA-67125, CGA-37735, CGA-41638, CGA-51202, CGA-62826, CGA-46576, CGA-380168, CGA-37913

TEST SYSTEM: Water

NO. OF ANALYSES: 144

LABORATORY: Development Resources/
Chemical Support Department

DESCRIPTION: Distilled water samples were fortified with a range of concentrations of each test substance listed above to assess the reactivity of the CGA-77102 enzyme immunoassay as part of the method validation of "Draft" Analytical Method 1004-98.

STUDY DIRECTOR: Jody L. Needham, B.A.

SIGNATURE: *Jody L. Needham*

DATE: *6/23/99*

MANAGEMENT REPRESENTATIVE: Max W. Cheung, Ph.D.

SIGNATURE: *M.W. Cheung*

DATE: *06/23/99*

SUBMITTED BY: Jody L. Needham and James F. Brady

QUALITY ASSURANCE AUDITOR:

SIGNATURE: *Veresa DCX for Kathleen Campbell*

DATE: *6/22/99*

DISTRIBUTION: J. Needham
J. Brady

M. W. Cheung
Main File

0029901 (S:77102\1004-98\RT\Doc) cb/sbh:6/21/99

RESIDUE TEST REPORT
 TEST NUMBER: RI-MV-002-99
 INVENTORY NO.: 18425.1
 REPORT NO.: 1

PROTOCOL NUMBER: 1004-98
 PROJECT NUMBER: H077102NAU950A
 TEST SUBSTANCES: See below
 TEST SYSTEM: Water

METHOD OF ANALYSIS: "Draft" 1004-98

TABLE 1. CROSS-REACTIVITY RESULTS FOR THE CGA-77102 ENZYME IMMUNOASSAY

TEST SUBSTANCE	¹ LOD	² LOQ	³ I ₅₀	PERCENT REACTIVITY RELATIVE TO CGA-77102	SET #	ANALYSIS DATE	NUMBER OF ANALYSES
⁴ CGA-77102	0.05	0.10	0.48	100	990104A	1/4/99	24
					990105A	1/5/99	
					990107A	1/7/99	
					990128A	1/28/99	
Metolachlor	0.05	--	1.0	55	990104A	1/4/99	6
Dimethanamid	0.20	--	4.2	12.9	990104A	1/4/99	6
CGA-380168	0.23	--	13.5	2.1	990105A	1/5/99	6
CGA-41638	--	--	46	<1.0	990105A	1/5/99	6
Acetochlor	--	--	68	<1.0	990104A	1/4/99	6
CGA-40172	--	--	84	<1.0	990107A	1/7/99	6
Alachlor	--	--	470	<1.0	990104A	1/4/99	6
CGA-354743	--	--	⁵ NR	NR	990107A	1/7/99	6
Metalaxyl	--	--	NR	NR	990107A	1/7/99	6
Acetochlor-ESA	--	--	NR	NR	990104A	1/4/99	6
Alachlor-ESA	--	--	NR	NR	990104A	1/4/99	6
Dimethanamid-ESA	--	--	NR	NR	990128A	1/28/99	6
CGA-50720	--	--	NR	NR	990107A	1/7/99	6
CGA-67869	--	--	NR	NR	990107A	1/7/99	6
CGA-67125	--	--	NR	NR	990107A	1/7/99	6
CGA-37735	--	--	NR	NR	990104A	1/4/99	6
CGA-51202	--	--	NR	NR	990105A	1/5/99	6
CGA-62826	--	--	NR	NR	990105A	1/5/99	6
CGA-46576	--	--	NR	NR	990105A	1/5/99	6
CGA-37913	--	--	NR	NR	990128A	1/28/99	6

¹ Limit of Detection (LOD), the smallest dose (in units of ng/ml) that yields a response that is significantly different from the response of the zero dose.

² Limit of Quantitation (LOQ), the lowest level of fortification (in units of ng/ml) that gives an acceptable recovery. An LOQ was determined for CGA-77102 only.

³ I₅₀, the concentration of test substance (in units of ng/ml) that yields half the response of the zero dose.

⁴ Parameters shown for CGA-77102 are the average of four determinations.

⁵ NR - not reactive.

0029901 (S:77102\1004-98\RT\Doc) sbh:6/21/99

RESIDUE TEST REPORT
 TEST NUMBER: RI-MV-002-99
 INVENTORY NO.: 18425.1
 REPORT NO.: 1

PROTOCOL NUMBER: 1004-98
 PROJECT NUMBER: H077102NAU950A
 TEST SUBSTANCES: See below
 TEST SYSTEM: Water

METHOD OF ANALYSIS: "Draft" 1004-98

TABLE 2. TEST SUBSTANCES USED IN PROTOCOL 1004-98

NAME	CODE NUMBER	% PURITY	EXPIRATION DATE
¹ CGA-77102	S95-1844	98.5	11/2000
Metolachlor	S94-1719	99.3	2/2000
² Alachlor	S86-0926	98.4	NON-GLP
² Alachlor-ESA	DAH-XXII-24	99.4	11/2000
² Acetochlor	S97-2120	98	8/99
² Acetochlor-ESA	DAH-XXII-30	97.8	11/2000
² Dimethanamid	S97-2034	99.7	7/2000
² Dimethanamid-ESA	BPM-XVIII-19	approx. 87.5	NON-GLP
Metalaxyl	S94-1763	97.1	10/2000
¹ CGA-354743	DAH-XXVI-12	96.8	9/2000
¹ CGA-50720	DPS-I-74-2	>99.9	6/99
¹ CGA-40172	JAK-II-42 (IV)	95.9	6/99
¹ CGA-67869	DAH-XV-35	>99.9	6/99
¹ CGA-67125	DPS-VII-41-1	99.9	12/99
¹ CGA-37735	JAK-II-42 (II)	99.2	3/2000
¹ CGA-41638	S94-1765	95	9/99
¹ CGA-51202	JAK-I-63	98.8	3/2000
¹ CGA-62826	BPM-I-48	98.0	11/2000
¹ CGA-46576	AW-I-43	97.4	6/99
¹ CGA-380168	DAH-XXIII-62	98.8	1/99
¹ CGA-37913	JAK-XVI-38	99.0	12/2000

¹These compounds lack common names.

²These compounds, or their metabolites, are not Novartis products.

APPENDIX 2

RESIDUE TEST REPORT NO. 2

Development Resources/Chemical Support Department
Novartis Crop Protection, Inc.
Greensboro, North Carolina

RESIDUE TEST REPORT

TEST NUMBER: RI-MV-002-99
INVENTORY NUMBER: 18425.1
REPORT NO.: 2

PROJECT NUMBER: H077102NAU950A
PROTOCOL NUMBER: 1004-98
TEST SUBSTANCE: CGA-77102
TEST SYSTEM: Water
NO. OF ANALYSES: 22 LABORATORY: Development Resources/
Chemical Support Department

DESCRIPTION: Distilled water samples were fortified with CGA-77102 and analyzed pre- and post-fortification as part of the method validation of "Draft" Analytical Method 1004-98.

STUDY DIRECTOR: Jody L. Needham, B.A.

SIGNATURE: *Jody L. Needham* DATE: *6/23/99*

MANAGEMENT REPRESENTATIVE: Max W. Cheung, Ph.D.

SIGNATURE: *m.w. Cheung* DATE: *06/23/99*

SUBMITTED BY: Jody L. Needham and James F. Brady

QUALITY ASSURANCE AUDITOR:

SIGNATURE: *Theresa Bliz for Kathleen Campbell* DATE: *6/22/99*

DISTRIBUTION: J. Needham
J. Brady
M. W. Cheung
Main File

0029902 (S:\77102\1004-98\RTR\Doc) sbh:6/21/99

Residue Test Report

Test Number: RI-MV-002-99
 Inventory No.: 18425.1
 Report No.: 2

Protocol Number: 1004-98
 Project Number: H077102NAU950A
 Method of Analysis: "Draft" 1004-98

Test Substance: CGA-77102
 Test System: Water

TABLE 1. RESULTS OF WATER STANDARD ADDITION EXPERIMENTS IN THE METHOD VALIDATION OF "DRAFT" 1004-98

Distilled Water Samples	ppb CGA-77102 Equivalents Found		Net	Fortification Level (ppb)	Analysis Set Name	Date	NB ref.
	Pre-fortification	Post-fortification					
A	<0.05	0.09	0.09	0.10	990126A	1/26/99	6593/64
B	<0.05	0.08	0.08	0.10	990126A	1/26/99	6593/64
C	<0.05	0.10	0.10	0.10	990126A	1/26/99	6593/64
D	<0.05	0.09	0.09	0.10	990126A	1/26/99	6593/64
E	<0.05	0.16	0.13	0.10	990126A	1/26/99	6593/64
F	<0.05	0.09	0.10	0.10	990126A	1/26/99	6593/64
G	<0.05	0.08	0.09	0.10	990126A	1/26/99	6593/64
H	<0.05	0.10	0.10	0.10	990126A	1/26/99	6593/64
I	<0.05	0.09	0.09	0.10	990126A	1/26/99	6593/64
J	<0.05	0.08	0.09	0.10	990126A	1/26/99	6593/64

Mean ± SD = 0.10 ± 0.01
 N = 10
 CV = 0.13

¹Net values were calculated by subtracting the mean pre-fortification result from the mean post-fortification result.

Residue Test Report

Test Number: RI-MV-002-99
 Inventory No.: 18425.1
 Report No.: 2

Protocol Number: 1004-98
 Project Number: H077102NAU950A
 Method of Analysis: "Draft" 1004-98

Test Substance: CGA-77102
 Test System: Water

TABLE 1. RESULTS OF WATER STANDARD ADDITION EXPERIMENTS IN THE METHOD VALIDATION OF "DRAFT" 1004-98 (continued)

Distilled Water Samples	ppb CGA-77102 Equivalents Found		Net ¹	CGA-77102 Fortification		Analysis Set Name	Date	NB ref.
	Pre-fortification	Post-fortification		Level (ppb)				
A	<0.05	1.03	1.12	1.0	990505A	5/5/99	6593/86	
B	<0.05	1.10	1.17	1.0	990505A	5/5/99	6593/86	
C	<0.05	1.36	1.37	1.0	990505A	5/5/99	6593/86	
D	<0.05	1.24	1.24	1.0	990505A	5/5/99	6593/86	
E	<0.05	1.21	1.25	1.0	990505A	5/5/99	6593/86	
F	<0.05	1.34	1.36	1.0	990505A	5/5/99	6593/86	
G	<0.05	1.15	1.24	1.0	990505A	5/5/99	6593/86	
H	<0.05	1.16	1.20	1.0	990505A	5/5/99	6593/86	
I	<0.05	1.21	1.25	1.0	990505A	5/5/99	6593/86	
J	<0.05	1.09	1.09	1.0	990505A	5/5/99	6593/86	

Mean ± SD = 1.23 ± 0.09
 N = 10
 CV = 0.07

¹Net values were calculated by subtracting the mean pre-fortification result from the mean post-fortification result.