

# Tools I Can't Live Without: Speed Up Method Development and Translation, Essentials for GC and GC/MS Users

Michelle Misselwitz, Julie Kowalski, Jack Cochran, Becca Stevens

[michelle.misselwitz@restek.com](mailto:michelle.misselwitz@restek.com)



*We need a new method developed...  
YESTERDAY!*



***You have:***  
Analyte list


***You need:***  
GC column and instrument  
conditions

# Speed Up and Simplify GC Method Development With Restek's EZGC® Online Suite



The image shows two overlapping software interfaces on a tablet. The top interface is the "EZGC Method Translator" and the bottom is the "EZGC Flow Calculator".

EZGC Method Translator		EZGC Flow Calculator	
Carrier Gas	Translation	Carrier Gas	
Helium	Helium	Helium	
30.00	11.00	Length	30.00 m
0.25	0.25	Inlet Diameter	0.25 mm
0.25	0.25	File Thickness	0.25 mm
200	200	Temperature	40.00 °C
1.40	1.40	Control Parameters	
42.74	40.44	Outlet Flow	0.40 mL/min
1.07	0.43	Average Velocity	42.74 cm/min
11.42	3.27	Hold Time	1.17 min
0.66	0.60	Inlet Pressure	11.42 PSI
		Outlet Pressure (set)	0.60 PSI
		Inlet	
		Temperature	200.00 °C
		Carrier Volume	0.00 mL
		Flow	1.40 mL/min
		Solvent Split Time	1.13055 min



New! EZGC® Method Translator and Flow Calculator



EZGC® Chromatogram Modeler

## EZGC<sup>®</sup> Chromatogram Modeler

[Send Us Your Feedback!](#)

Enter compound names, synonyms, or CAS #s.

**Examples:**

1,4-dioxane  
acetaldehyde, ethanol and isopropyl alcohol

[I need to paste a list. Show me a list box.](#)



Welcome to the Restek EZGC<sup>®</sup> Chromatogram Modeler, the easiest way to jump-start your method development! To get started, enter the compounds you wish to separate into the field above. To paste a list, first click the link "I need to paste a list" underneath the input field.

### Browser requirements

To get the best results from the EZGC<sup>®</sup> modeler, we recommend one of the following browsers:

- Firefox (desktop or Android tablet)
- Chrome (Windows or Mac desktop or Android tablet)
- IE 8 or IE 9 or above (Windows desktop)
- Safari (desktop or iPad tablet)
- Opera (desktop or mobile)

### Learn EZGC<sup>®</sup> in 5 Minutes!

Find out how easy it is to get started with the EZGC<sup>®</sup> app in this brief screencast.

[Watch the Video](#)

### It's Easy, but Powerful!

Restek Innovations Laboratory Manager Chris English goes into detail about how to get the most out of the EZGC modeler.

[Watch the Video](#)

## My EZGC<sup>®</sup> Models

Untitled

1,1,1,2-Tetrachloroethane, 1,1,1-Trichloroethane, 1,1,2-Trichloroethane, 1,1-Dichloroethene, 1,2,4-Trichlorobenzene...

Dec 14, 2012

[View/edit full list](#)



# EZGC® Chromatogram Modeler

Send Us Your Feedback!

Enter a compound name, synonym or CAS #, one per line.

Go

[Show me the standard entry field.](#)

2,4'-DDD[x], 2,4'-DDE[x], 2,4'-DDT[x], 4,4'-DDD[x], 4,4'-DDE[x], 4,4'-DDT[x], 4,4'-Dibromobiphenyl[x], Captan [x], Carbophenothion[x], Chlorobenzilate[x], Decachlorobiphenyl[x], Dibutyl chlorendate[x], Dicofol[x], Dieldrin[x], Endosulfan I[x], Endosulfan II[x], Endosulfan sulfate[x], Endrin[x], Endrin aldehyde[x], Endrin ketone[x], Heptachlor epoxide[x], Isodrin[x], Kepone[x], Methoxychlor[x], Mirex[x], [Show All...](#)

The following compounds were not found:

a-chlordane[x]  
acifluorfen me[x]  
γ chlordane[x]

[\[Clear All\]](#)



Need more help?

[Send this list](#) to Restek Technical Service with your questions.

Untitled. [Click here to edit.](#)



*Click title to edit*

Peaks	Ret. Time (tr)	Res. (R <sub>s</sub> )	Peak Width
1. Isodrin	30.21	8.31	0.10
2. 4,4'-Dibromobiphenyl	31.05	--	0.10
3. Heptachlor epoxide	31.05	--	0.10
4. 2,4'-DDE	32.14	3.58	0.10
5. Captan	32.49	3.58	0.10
6. Endosulfan I	33.00	1.02	0.10
7. cis-Nonachlor	33.10	1.02	0.10
8. 4,4'-DDE	34.08	2.51	0.10
9. Dieldrin	34.33	2.14	0.10
10. 2,4'-DDD	34.55	2.14	0.10
11. Perthane	35.01	3.97	0.10
12. Endrin	35.39	3.97	0.10
13. 2,4'-DDT	35.83	4.12	0.10
14. trans-Nonachlor	36.29	1.63	0.10
15. Kepone	36.46	0.13	0.11
16. Chlorobenzilate	36.47	0.13	0.10
17. 4,4'-DDD	36.76	2.94	0.10
18. Endosulfan II	37.13	1.98	0.10
19. Carbophenothion	37.34	1.98	0.10
20. Endrin aldehyde	37.70	3.26	0.10
21. 4,4'-DDT	38.04	3.26	0.10
22. Endosulfan sulfate	39.11	10.60	0.10
23. Methoxychlor	40.42	3.10	0.10
24. Dibutyl chlorendate	40.73	1.94	0.10
25. Endrin ketone	40.92	1.56	0.11
26. Dicofol	41.09	1.56	0.10
27. Mirex	42.69	13.10	0.11
28. cis-Permethrin	44.14	4.07	0.10
29. trans-Permethrin	44.54	4.07	0.10
30. Decachlorobiphenyl	47.47	29.80	0.10

Not found: a-chlordane, acifluorfen me, γ chlordane

## Product Suggestions

Rtx®-PCB Columns  
cat. # 13223  
\$537.0/ea.

Qty:

[Add To Cart](#)

## My EZGC® Models

Untitled


1,1,1,2-Tetrachloroethane, 1,1,1-Trichloroethane, 1,1,2-Trichloroethane, 1,1-Dichloroethane, 1,2,4-Trichlorobenzene...

Dec 14, 2012


[View/edit full list](#)

14. trans-Nonachlor	36.29	1
15. Kepone	36.46	0
16. Chlorobenzilate	36.47	0
17. 4,4'-DDD	36.76	2
18. Endosulfan II	37.13	1
19. Carbophenothion	37.34	1
20. Endrin aldehyde	37.70	3
21. 4,4'-DDT	38.04	3
22. Endosulfan sulfate	39.11	10
23. Methoxychlor	40.42	3
24. Dibutyl chlorendate	40.73	1
25. Endrin ketone	40.92	1
26. Dicofol	41.09	1
27. Mirex	42.69	13
28. cis-Permethrin	44.14	4
29. trans-Permethrin	44.54	4
30. Decachlorobiphenyl	47.47	29

Not found:  $\alpha$ -chlordane, acifluorfen me,  $\gamma$ -chlordane



**Column** Rtx®-PCB, 30 m, 0.25 mm ID, 0.25  $\mu$ m (cat.# 13223)  
**Carrier Gas:** Helium, constant flow @ 1.68 mL/min  
**Linear Velocity:** 40.0 cm/sec  
**Outlet Pressure (abs):** atmospheric pressure  
**Oven Temp:** 100 °C to 310 °C @ 4.0 °C/min



**Column** Rtx®-PCB, 30 m, 0.25 mm ID, 0.25  $\mu$ m (cat.# 13223)  
**Carrier Gas:** Helium, constant flow @ 1.68 mL/min  
**Linear Velocity:** 40.0 cm/sec  
**Outlet Pressure (abs):** atmospheric pressure  
**Oven Temp:** 100 °C to 310 °C @ 4.0 °C/min

*We need a new method developed...  
YESTERDAY!*



***You have:***  
Analyte list  
GC column

***You need:***  
GC instrument conditions

# We need to know...

Flow rate

Oven program

Accurate column length

Initial oven temperature

Splitless hold time



# Optimal Flow for GC

SOF EOF

# Optimal Flow for GC

## SOF

**S**peed

**O**ptimized

**F**low

## EOF

**E**fficiency

**O**ptimized

**F**low

# Theory of Fast Capillary Gas Chromatography – Part 3: Column Performance vs. Gas Flow Rate

Leonid M. Blumberg<sup>1)</sup>

Fast GC Consulting, PO Box 585, Hockessin, DE 19707, USA

J. High Resol. Chromatogr. **1999**, 22, (7) 403-413

Ms received: June 2, 1998; accepted: October 15, 1998

**Key Words:** *Constant length optimization; constant efficiency optimization; efficiency-optimized flow rate; fast GC; film inefficiency factor; high pressure drop; speed-optimized flow rate*

## Speed-Optimized Flow and Efficiency-Optimized Flow

### Summary

At the high pressure drop required for the fast analysis of complex mixtures, the equations for the column plate height,  $H$ , and plate duration,  $Q$ , as functions of the carrier gas velocity,  $\bar{u}$ , differ substantially from the equations for the same quantities expressed via the carrier gas flow rate,  $F$ . While  $\bar{u}$  as an independent pneumatic variable is more convenient for the theoretical studies,  $F$  is a more convenient as a control parameter in practical applications. Equations for  $H$  vs.  $\bar{u}$  and for  $Q$  vs.  $\bar{u}$  from Parts 1 and 2 are transformed here into expressions for  $H$  vs.  $F$  and  $Q$  vs.  $F$ . An efficiency-optimized flow rate (EOF) and a speed-optimized flow rate (SOF) are found. Expressions for these two quantities are considerably simpler than their velocity-based counterparts. In particular, SOF does not depend on column length, film thickness, and pressure drop.

### 1 Introduction

This report continues the study [1, 2] of the speed-separation performance of a capillary column having an arbitrary film thickness and operating at an arbitrary pressure drop.

At a low pressure drop, the structure of the equations for the plate height,  $h$ , and plate duration,  $q$ , does not depend on the type of the independent pneumatic variable. Thus,

$$h = \frac{b}{w} + cw, \quad q = \frac{b}{w^2} + c$$

are, respectively, Van Deemter and Purnell [5, 6] equations in which the independent variable  $w$  can be either  $\bar{u}$  or  $F$  (or the column pressure drop) while  $b$  and  $c$  are the coefficients whose values do not change with the change in the value of  $w$ .

At the high pressure drop, however, the structure of the expression for both the plate height,  $H$ , and the plate duration,  $Q$ , depends on the selection of the independent pneumatic variable. Thus, for example, as a function of  $\bar{u}$ ,  $H$  can be expressed as (see Part 1)

$$H = \frac{B}{\bar{u}^2} + C_1\bar{u}^2 + C_2\bar{u}$$



Contents lists available at SciVerse ScienceDirect

## Journal of Chromatography A

journal homepage: [www.elsevier.com/locate/chroma](http://www.elsevier.com/locate/chroma)



Discussion

### Plate height formula widely accepted in GC is not correct

Leonid M. Blumberg


Fast GC Consulting, PO Box 1243, Wilmington, DE 19801, USA

*Dear Editor,*

I would like to bring to attention of the readers of your Journal that


- widely accepted formula for a column plate height ( $H$ ) in GC is incorrect, and that
- carrier gas average velocity ( $\bar{u}$ ) is not the best variable for describing  $H$  as a function of a gas flow in GC columns.

# Van Deemter Plot





# Van Deemter Plot





Contents lists available at SciVerse ScienceDirect

## Journal of Chromatography A

journal homepage: [www.elsevier.com/locate/chroma](http://www.elsevier.com/locate/chroma)



Discussion

### Plate height formula widely accepted in GC is not correct

Leonid M. Blumberg

Fast GC Consulting, PO Box 1243, Wilmington, DE 19801, USA

Suggests optimal FLOW rate for GC,  
not optimal average velocity...

*Dear Editor,*

I would like to bring to attention of the readers of your Journal that

- widely accepted formula for a column plate height ( $H$ ) in GC is incorrect, and that
- carrier gas average velocity ( $\bar{u}$ ) is not the best variable for describing  $H$  as a function of a gas flow in GC columns.

# Optimal Flow for GC

## SOF

Speed **O**ptimized **F**low

Column ID × constant

Helium

8 x column ID

8 x 0.25 = 2.0 mL/min

Hydrogen

10 x column ID

10 x 0.25 = 2.5 mL/min

## EOF

Efficiency **O**ptimized **F**low

$\frac{\text{SOF}}{\sqrt{2}}$

Helium

2.0 /  $\sqrt{2}$  = 1.4 mL/min

Hydrogen

2.5 /  $\sqrt{2}$  = 1.8 mL/min

# Optimal Flow for GC

## SOF

Speed **O**ptimized **F**low

Helium, 0.25 mm ID

**2.0 mL/min**

## EOF

Efficiency **O**ptimized **F**low

Helium, 0.25 mm ID

**1.4 mL/min**



- ✓ Considerably simpler than velocity-based
- ✓ SOF does not depend on column length, film thickness, and pressure drop

# We need to know...

Flow rate



Oven program

Accurate column length

Initial oven temperature

Splitless hold time



# Optimal Heating Rate in Gas Chromatography

L. M. Blumberg,<sup>1</sup> M. S. Klee<sup>2</sup>

<sup>1</sup>*Fast GC Consulting, P.O. Box 585, Hockessin, DE 19707, USA*

<sup>2</sup>*Agilent Technologies, 2850 Centerville Road, Wilmington, DE 19808, USA*

Received 8 June 2000; accepted 15 September 2000

J. Microcolumn Separations, 12 (9) 508-514 (2000)

10 /  $t_M$

**Abstract:** Several optimization criteria for column heating rate in temperature programmed gas chromatography (GC), under different optimization constraints are identified. Applying these criteria to experimental data, it is shown that when column pressure drop is high, the optimum heating rate for *n*-alkanes and pesticides in a column with a silicone stationary phase of a typical thickness is about 10°C per void time. This heating rate is recommended as a default for all temperature programs in capillary GC. © 2000 John Wiley & Sons, Inc. J Micro Sep 12: 508–514, 2000

**Key words:** *gas chromatography; optimal heating rate; separation-speed tradeoff*


# Optimal Heating Rate for GCC


- Adequate separation of a required number of solutes in the shortest time
  - Maximizing peak capacity
  - Minimizing analysis time
- NOT aimed specifically at separating specific pairs of solutes

$$\frac{10^{\circ}\text{C}}{t_M}$$

Air injection  
100 split, 0.1  $\mu\text{L}$   
Scan 25-100

# Optimal Heating Rate for GCC



$$\frac{10^{\circ}\text{C}}{t_M} \rightarrow \frac{10^{\circ}\text{C}}{1.18 \text{ min}} = 8.5^{\circ}\text{C}/\text{min}$$


# We need to know...

Flow rate



Oven program



Accurate column length

Initial oven temperature

Splitless hold time

# EZGC™ Method Translator and Flow Calculator

- <http://www.restek.com/ezgc-mtfc>
- Available in a web-based version
- Or, download the app from the link above
- Windows 8 – 7 – Vista – XP

The screenshot displays two side-by-side software windows. The left window is titled 'EZGC Method Translator' and features a table for comparing 'Original' and 'Translation' values for various parameters. The right window is titled 'EZGC Flow Calculator' and shows a control panel with input fields and sliders for flow rate, velocity, and pressure, along with a 'Results' section at the bottom.

Carrier Gas	Original	Translation
Helium	Helium	Helium

Column	Original	Translation
Length	30.00 m	15.00 m
Inner Diameter	0.25 mm	0.25 mm
Film Thickness	0.25 µm	0.25 µm
Phase Ratio	230	230

Control Parameters	Original	Translation
Outlet Flow	1.40 mL/min	1.40 mL/min
Average Velocity	42.74 cm/sec	60.44 cm/sec
Holdup Time	1.17 min	0.41 min
Inlet Pressure	11.42 psi	3.77 psi
Outlet Pressure	0.00 psi	0.00 psi

Oven Program	Original	Translation
Number of Ramps (1-4)	40	40
Ramp Rate (°C/min)	8.5	330
Temp (°C)	1	24
Hold Time (min)	1	330
Ramp Rate (°C/min)	40	0.33

Control Method	Original	Translation
Constant Flow	Constant Flow	Constant Flow

Results	Original	Translation
Run Time	36.12 min	12.78 min
Speed	2.83 x	2.83 x

EZGC Flow Calculator	Original	Translation
Carrier Gas	Helium	Helium

Column	Original	Translation
Length	30.00 m	15.00 m
Inner Diameter	0.25 mm	0.25 mm
Film Thickness	0.25 µm	0.25 µm
Temperature	40.00 °C	40.00 °C

Control Parameters	Original	Translation
Outlet Flow	1.40 mL/min	1.40 mL/min
Average Velocity	42.74 cm/sec	60.44 cm/sec
Holdup Time	1.17 min	0.41 min
Inlet Pressure	11.42 psi	3.77 psi
Outlet Pressure	0.00 psi	0.00 psi

Inlet	Original	Translation
Temperature	250.00 °C	250.00 °C
Linear Volume	1.00 mL	1.00 mL
Flow	1.40 mL/min	1.40 mL/min
Splitless Valve Time	1.1 to 1.5 min	1.1 to 1.5 min



## EZGC™ Method Translator

### Carrier Gas

Helium  Helium

### Column

Length	30.00	15.00	m
Inner Diameter	0.25	0.25	mm
Film Thickness	0.25	0.25	µm
Phase Ratio	250	250	

### Control Parameters

Outlet Flow	<input type="text"/> 1.40	<input type="text"/> 1.40	mL/min
Average Velocity	42.74	60.44	cm/sec
Holdup Time	1.17	0.41	min
Inlet Pressure	psi <input type="text"/> 11.42	3.77	psi
Outlet Pressure (abs)	<input type="text"/> 0.00	<input type="text"/> 0.00	psi
	<input type="button" value="Atm"/> <input type="button" value="Vacuum"/>	<input type="button" value="Atm"/> <input type="button" value="Vacuum"/>	

### Oven Program

- Isothermal  
 Ramps

Number of Ramps (1-4)	Ramp 1			Ramp 2		
	Ramp Rate (°C/min)	Temp (°C)	Hold Time (min)	Ramp Rate (°C/min)	Temp (°C)	Hold Time (min)
1	8.5	330	1	24	330	0.35

### Control Method

Constant Flow

### Results

Solve for  Efficiency  Speed  Translate  Custom

Run Time	36.12	12.78	min
Speed		2.83	x

## EZGC™ Flow Calculator

### Carrier Gas

Helium

### Column

Length	30.00	m
Inner Diameter	0.25	mm
Film Thickness	0.25	µm
Temperature	40.00	°C

### Control Parameters

Outlet Flow	<input type="text"/> 1.40	mL/min
Average Velocity	42.74	cm/sec
Holdup Time	1.17	min
Inlet Pressure	psi <input type="text"/> 11.42	psi
Outlet Pressure (abs)	<input type="text"/> 0.00	psi
	<input type="button" value="Atm"/> <input type="button" value="Vacuum"/>	

### Inlet

Temperature	250.00	°C
Liner Volume	1.00	mL
Flow	1.40	mL/min
Splitless Valve Time	1.1 to 1.5	min

[Download EZGC™ Method Translator and Flow Calculator](#)

For Windows 8/7/Vista/XP

# GC Column Length Determination

- Using Holdup Time and the Flow Calculator
- Upon initial GC column installation
- Enables accurate Electronic Pneumatic Control of GC carrier gas flow
  
- Also used after column trimming

# EZGC™ Flow Calculator

## Carrier Gas

Helium

## Column

Length	30.00	m
Inner Diameter	0.25	mm
Film Thickness	0.25	μm
Temperature	90	°C

## Control Parameters

Outlet Flow	Optimum Range 1.4 to 2.0 mL/min	→	1.40	mL/min
Average Velocity			43.73	cm/sec
Holdup Time			1.14	min
Inlet Pressure (gauge)			14.91	psi
Outlet Pressure (abs)			0.00	psi

Atm


Vacuum

## Inlet

Temperature	250	°C
Liner Volume	0.99	mL
Flow	1.23	mL/min
Splitless Valve Time	1.2 to 1.7	min

Use MT Original Values

Use MT Translation Values



# EZGC™ Flow Calculator

## Carrier Gas

Helium

## Column

Length 30.00 m  
Inner Diameter 0.25 mm  
Film Thickness 0.25 μm  
Temperature 90 °C

## Control Parameters

Outlet Flow Optimum Range 1.4 to 2.0 mL/min 1.40 mL/min  
Average Velocity 43.73 cm/sec  
Holdup Time 1.14 min  
Inlet Pressure (gauge) 14.91 psi  
Outlet Pressure (abs) 0.00 psi

## Inlet

Temperature °C  
Liner Volume 0.99 mL  
Flow 1.23 mL/min  
Splitless Valve Time 1.2 to 1.7 min

Use MT Original Values

Use MT Translation Values

Double click here to "lock" Inlet Pressure

# EZGC™ Flow Calculator

## Carrier Gas

Click the "spinner" to increase Length

## Column

Length 30.10 m  
Inner Diameter 0.25 mm  
Film Thickness 0.25 μm  
Temperature 90 °C

## Control Parameters

Outlet Flow Optimum Range 1.4 to 2.0 mL/min 1.40 mL/min  
Average Velocity 43.58 cm/sec  
Holdup Time 1.15 min  
Inlet Pressure (gauge) 14.91 psi  
Outlet Pressure (abs) 0.00 psi

## Inlet

Temperature 250 °C  
Liner Volume 0.99 mL  
Flow 1.23 mL/min  
Splitless Valve Time 1.2 to 1.7 min

Use MT Original Values

Use MT Translation Values

The "spinner" advances the Holdup Time

# EZGC™ Flow Calculator

## Carrier Gas

Helium ▾

## Column

Length	30.70 m
Inner Diameter	0.25 mm
Film Thickness	0.25 μm
Temperature	90 °C

## Control Parameters


Outlet Flow	Optimum Range 1.4 to 2.0 mL/min	1.37 mL/min
Average Velocity		42.74 cm/sec
Holdup Time		1.20 min
Inlet Pressure (gauge)	→	14.91 psi ▾
Outlet Pressure (abs)		0.00 psi
	Atm Vacuum	

## Inlet

Temperature	250 °C
Liner Volume	0.99 mL
Flow	1.20 mL/min
Splitless Valve Time	1.2 to 1.7 min

Use MT Original Values

Use MT Translation Values





# EZGC™ Flow Calculator

## Carrier Gas

Helium

## Column

Length	30.70 m
Inner Diameter	0.25 mm
Film Thickness	0.25 μm
Temperature	90 °C


## Control Parameters

Outlet Flow	Optimum Range 1.4 to 2.0 mL/min	1.37 mL/min
Average Velocity		42.74 cm/sec

## Control Parameters

Outlet Flow	Optimum Range 1.4 to 2.0 mL/min	1.37 mL/min
Average Velocity		42.74 cm/sec
Holdup Time		1.20 min
Inlet Pressure (gauge)	→	14.91 psi
Outlet Pressure (abs)		0.00 psi

Atm Vacuum



# Combining Optimized Flow and Optimal Heating Rate

GC Column	Detector	He EOF mL/min	Avg Vel cm/sec	Holdup time (min)	OHR °C/min	Anal time min
60m x 0.25mm x 0.25µm	MS	1.4	31	3.23	3.1	74.2
30m x 0.25mm x 0.25µm	MS	1.4	44	1.14	8.8	26.1
15m x 0.25mm x 0.25µm	MS	1.4	62	0.40	25.0	9.2
20m x 0.18mm x 0.18µm	MS	1.0	39	0.74	13.5	17.0
60m x 0.25mm x 0.25µm	ECD	1.4	27	3.71	2.7	85.1
30m x 0.25mm x 0.25µm	ECD	1.4	35	1.45	6.9	33.3
15m x 0.25mm x 0.25µm	ECD	1.4	42	0.60	16.7	13.8

Mass spectrometer is a vacuum-outlet detector.

Holdup time is at 90°C.

Analysis time is based on 90 to 320°C oven program.

Methiocarb

30.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

1.4 mL/min

8.5 $^{\circ}$ C/min

Fenthion

Chlorpyrifos

Dichlofluanid

Pirimiphos  
methyl

Malathion

Time (s)

895

900

905

910

915

920

925

930

— 168


— 290

— 167

— 173

— 278


— 197



30.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms


1.0 mL/min; 7.1 $^{\circ}$ C/min

18.41min




1.4 mL/min; 8.5 $^{\circ}$ C/min

15.45min



2.0 mL/min; 10.1 $^{\circ}$ C/min


13.00min



30.7m x 0.25mm x 0.25µm Rxi-5ms

1.0 mL/min; 7.1°C/min


18.41min



Efficiency-optimized flow

1.4 mL/min; 8.5°C/min


15.45min



Speed-optimized flow

2.0 mL/min; 10.1°C/min

13.00min




30.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

1.4 mL/min; 8.5 $^{\circ}$ C/min

Efficiency-optimized flow


15.45min



2.0 mL/min; 10.1 $^{\circ}$ C/min


Speed-optimized flow

13.00min



4.0 mL/min; 14.3 $^{\circ}$ C/min

9.23min



# We need to know...

Flow rate



Oven program



Accurate column length




Initial oven temperature

Splitless hold time




# Initial Oven Temperature



# Initial Oven Temperature

## SOLVENT FOCUSING



Boiling point of 1<sup>st</sup> analyte


Boiling point of solvent

**20-40°C lower**

Initial Oven Temperature

# Initial Oven Temperature

## ANALYTE FOCUSING



Boiling point of 1<sup>st</sup> analyte

20°C lower



60-80°C

20°C higher

Boiling point of solvent


# Solvent Focusing: Solvent Polarity Mismatch

## Mismatched Polarity



o-Phenylphenol







## Polarity Match



o-Phenylphenol

**If solvent boiling point is NOT below that of lowest analyte boiling point focusing of both solvent and analytes**

# Initial Oven Temperature



**Solvent Focused:**  
Acetonitrile enters column as liquid and beads because of solvent/stationary phase mismatch

**Analyte Focused:**  
Acetonitrile is flash vaporized and not cold trapped (condenses) on the front of the column

**No Cold Trapping**  
Both the acetonitrile and analyte (o-Phenylphenol) are NOT cold trapped effectively on the front of the column

o-Phenylphenol, XIC at m/z 170

# We need to know...

Flow rate



Oven program



Accurate column length



Initial oven temperature




Splitless hold time


# Splitless Valve Time

Calculating the Splitless Valve Time with  
the Flow Calculator



# Calculating the GC Inlet Liner Volume for Splitless Valve Time



- Volume of a cylinder
  - $V = \pi r^2 h$
- Inlet liner for Agilent GC
  - 3.1416 ( $\pi$ ) 
  - 4 mm diameter, 2 mm radius ( $r$ )
  - 78.5 mm height ( $h$ )
- Liner  $V = 986 \mu\text{L}$  (**0.99 mL**)
  - Ignore wool and taper

# Calculating the GC Inlet Liner Volume for Splitless Valve Time


Liner volume 0.99 mL  
Want 1.5 to 2x sweep...

Inlet	
Temperature	250 °C
Liner Volume	0.99 mL
Flow	1.22 mL/min
Splitless Valve Time	1.2 to 1.7 min

EZGC™ Flow Calculator  
Splitless Valve Time  
range is 1.2 to 1.7 min.

EZGC™ Flow Calculator	
Carrier Gas	
	Helium
Column	
Length	30.70 m
Inner Diameter	0.25 mm
Film Thickness	0.25 µm
Temperature	90 °C
Control Parameters	
Outlet Flow	Optimum Range 1.4 to 2.0 mL/min → 1.40 mL/min
Average Velocity	43.23 cm/sec
Holdup Time	1.18 min
Inlet Pressure (gauge)	15.25 PSI
Outlet Pressure (abs)	0.00 PSI
	<input type="button" value="Atm"/> <input type="button" value="Vacuum"/>
Inlet	
Temperature	250 °C
Liner Volume	0.99 mL
Flow	1.22 mL/min
Splitless Valve Time	1.2 to 1.7 min





# We need to know...

Flow rate



Oven program



Accurate column length



Initial oven temperature



Splitless hold time



# EZGC™ Method Translator / Flow Calculator



*Make our methods  
better, faster, cheaper!•*

- Increasing speed of analysis
  - Decreasing column L and/or ID
  - Switching to a faster carrier gas (e.g., He to H<sub>2</sub>)
- Updating oven temperature program after column trimming for maintenance
- Improving Original methods in separation and/or speed of analysis
- Translating methods from GC-FID (or other atmospheric outlet detector) to GC-MS (vacuum outlet) or vice versa

# EZGC™ Method Translator

**Carrier Gas**  
 Original: Helium Translation: Helium

**Column**  
 Length:   m  
 Inner Diameter:   mm  
 Film Thickness:   μm  
 Phase Ratio:

**Control Parameters**  
 Outlet Flow:   mL/min  
 Average Velocity:   cm/sec  
 Holdup Time:   min  
 Inlet Pressure: psi   psi  
 Outlet Pressure (abs):   psi  
Atm Vacuum Atm Vacuum

**Oven Program**  
 Isothermal  
 Ramps  

Number of Ramps (1-4)	Ramp Rate (°C/min)	Temp (°C)	Hold Time (min)	Ramp Rate (°C/min)	Temp (°C)	Hold Time (min)
1	8.5	330	1	24	330	0.35

**Control Method**  
Constant Flow

**Results**  
 Solve for:  Efficiency  Speed  Translate  Custom  
 Run Time:   min  
 Speed:  x

Use Flow Calculator Values

**Carrier Gas**  
Helium

**Column**  
 Length:  m  
 Inner Diameter:  mm  
 Film Thickness:  μm  
 Temperature:  °C

**Control Parameters**  
 Outlet Flow: Optimum Range 1.4 to 2.0 mL/min  mL/min  
 Average Velocity:  cm/sec  
 Holdup Time:  min  
 Inlet Pressure: psi  psi  
 Outlet Pressure (abs):  psi  
Atm Vacuum

**Inlet**  
 Temperature:  °C  
 Liner Volume:  mL  
 Flow:  mL/min  
 Splitless Valve Time:  min

Use MT Original Values Use MT Translation Values

Download  
**EZGC™ Method Translator and Flow Calculator**  
 For Windows 8/7/Vista/XP





## **EZGC™ Method Translator and Flow Calculator Glossary**

The EZGC™ Method Translator is a tool built for gas chromatography (GC) method development. Generally, the goal of Method Translation is to allow alteration of GC column format, carrier gas, flow, etc., while keeping peak elution order—NOT retention times—the same. (Note that Method Translation assumes that the GC stationary phase type remains the same between Original and Translation methods.)

Some of the most practical uses for Restek's EZGC™ Method Translator are listed below:

- Increasing speed of analysis through decreasing column length and/or decreasing inner diameter and/or switching to a faster carrier gas (e.g., going from helium to hydrogen).
- Updating the oven temperature program through Translation after column trimming for maintenance so peak elution orders do not change.
- Improving Original methods in separation and/or speed of analysis by solving for Efficiency or Speed in Translation.
- Translating methods from GC-FID (or other atmospheric outlet detector) to GC-MS (vacuum outlet) or vice versa.

### **Basic Navigation in the EZGC™ Method Translator and Flow Calculator**

"White" cells are user-entry cells. "Blue" cells are locked cells that contain calculated values. In the Method Translator, the Translation's Control Parameters can be unlocked by selecting the Custom translation method in the Results section.


Highlighting numerical values using the mouse allows easy user entry of new values. A double mouse click in any user-entry cell highlights the value automatically for user entry. Hitting the Tab key while in a cell updates the cell with the user entered value and moves to the next cell for additional user entry, if necessary.

In the Control Parameters section for both the Method Translator and Flow Calculator, a double mouse click in the Outlet Flow, Average Velocity, Holdup Time, or Inlet Pressure cell will make that cell the "set point" around which the other control parameters are calculated. Column dimensions (and Temperature, in the Flow Calculator) can then be changed, and the set point value will remain fixed. A blue arrow denotes the "set point" cell.




# Column Trimming for Maintenance

Before



After injecting that DIRTY extract!!



# No Method Translation

- Trimming column for maintenance
- But, not updating column length for flow control
- And, updating column length, but not translating oven temperature program



I just want to trim my column  
and change **NOTHING!**

# What if you don't translate?




I just want to trim my column  
and change **NOTHING!**



I will input my new column  
length but **NOT** my oven  
program rate!


30.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

Original Method



29.6m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

No column length  
update or method  
translation



27.6m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

No column length  
update or method  
translation

Malathion

Aldrin


Fenthion

Chlorpyrifos

Parathion

23.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

No column length  
update or method  
translation







23.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms



I will change my column length but that is IT!!




30.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms




23.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

Method Translated




30.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms



23.7m x 0.25mm x 0.25 $\mu$ m Rxi-5ms

NOT Translated



# EZGC™ Method Translator

Carrier Gas	Original	Translation				
	Helium	Helium				
Column						
Length	30.70	23.70 m				
Inner Diameter	0.25	0.25 mm				
Film Thickness	0.25	0.25 μm				
Phase Ratio	250	250				
Control Parameters						
Outlet Flow	→ 1.40	→ 1.40 mL/min				
Average Velocity	43.23	49.20 cm/sec				
Holdup Time	1.18	0.80 min				
Inlet Pressure (gauge)	15.25	11.62 PSI				
Outlet Pressure (abs)	0.00	0.00 PSI				
	Atm Vacuum	Atm Vacuum				
Oven Program						
<input type="radio"/> Isothermal	Ramp (°C/min)	Temp (°C)	Hold (min)	Ramp (°C/min)	Temp (°C)	Hold (min)
<input checked="" type="radio"/> Ramps		90	0.1		90	0.05
Number of Ramps						
1 (1-4)	8.5	330	1	12.5	330	0.7
Control Method						
	Constant Flow					
Results						
Solve for	<input type="radio"/> Efficiency	<input type="radio"/> Speed	<input checked="" type="radio"/> Translate	<input type="radio"/> Custom		
Run Time	29.34	19.95	min			
Speed		1.47	x			
Use FC Values for Original		Use FC Values for Translation				

- After column trimming, use “Speed” to predict new retention times
- Actual retention times previous method divided by “Speed” factor = predicted retention times for translated method

# EZGC™ Method Translator

Carrier Gas      Original      Translation

Helium      Helium

## Column

Length	30.70	23.70	m
Inner Diameter	0.25	0.25	mm
Film Thickness	0.25	0.25	µm
Phase Ratio	250	250	

## Control Parameters

Outlet Flow	→ 1.40	→ 1.40	mL/min
Average Velocity	43.23	49.20	cm/sec
Holdup Time	1.18	0.80	min
Inlet Pressure (gauge)	15.25	11.62	PSI ▾
Outlet Pressure (abs)	0.00	0.00	PSI
	Atm Vacuum	Atm Vacuum	

## Oven Program

<input type="radio"/> Isothermal	Ramp (°C/min)	Temp (°C)	Hold (min)	Ramp (°C/min)	Temp (°C)	Hold (min)
<input checked="" type="radio"/> Ramps		90	0.1		90	0.05
Number of Ramps						
1 (1-4)	8.5	330	1	12.5	330	0.7

## Control Method

Results      Solve for       Run Time       Speed

Results      Solve for       Efficiency       Speed       Translate       Custom

Run Time      29.34      19.95 min

Speed      1.47 x

Use FC Values for Original

Use FC Values for Original

Use FC Values for Translation


- After column trimming, use “Speed” to predict new retention times
- Actual retention times previous method divided by “Speed” factor = predicted retention times for translated method

# EZGC™ Method Translator / Flow Calculator

Increasing speed of analysis

Decreasing column L and/or ID

Switching to a faster carrier gas (e.g., He to H<sub>2</sub>)



# EZGC™ Method Translator

## Carrier Gas

Original: Helium    Translation: Helium

## Column

Length	30.70	20.10	m
Inner Diameter	0.25	0.18	mm
Film Thickness	0.25	0.18	µm
Phase Ratio	250	250	

## Control Parameters

Outlet Flow	→ 1.40	→ 1.01	mL/min
Average Velocity	43.23	45.33	cm/sec
Holdup Time	1.18	0.74	min
Inlet Pressure (gauge)	15.25	24.97	psi ▾
Outlet Pressure (abs)	0.00	0.00	psi

## Oven Program

<input type="radio"/> Isothermal	Ramp (°C/min)	Temp (°C)	Hold (min)	Ramp (°C/min)	Temp (°C)	Hold (min)	
<input checked="" type="radio"/> Ramps		90	0.1		90	0.05	
Number of Ramps	1 (1-4)	8.5	330	1	13.6	330	0.65

## Control Method

Constant Flow ▾

## Results

Solve for  Efficiency  Speed  Translate  Custom

Run Time	29.34	18.35	min
Speed		1.60	x

Method Translator sets flow and oven rate for Translation column.

→ 1.40	→ 1.01	mL/min
43.23	45.33	cm/sec
1.18	0.74	min
15.25	24.97	psi ▾
0.00	0.00	psi



    

Ramp (°C/min)	Temp (°C)	Hold (min)	Ramp (°C/min)	Temp (°C)	Hold (min)
	90	0.1		90	0.05
8.5	330	1	13.6	330	0.65


## Constant Flow ▾

Efficiency  Speed  Translate  Custom

29.34	18.35	min
	1.60	x








**Results**   Solve for    Efficiency    Speed    Translate    Custom



Run Time	29.34	18.35 min
Speed	<b>15.45 / 9.59 = 1.61</b>	1.60 x

**5ms**



# Helium to Hydrogen Carrier Gas

- 20m x 0.18mm x 0.18 $\mu$ m



20.1m x 0.18mm x 0.18 $\mu$ m Rxi-5ms

Methiocarb

Helium at 1.0 mL/min  
Oven at 13.6 $^{\circ}$ C/min

Fenthion

Chlorpyrifos

Pirimiphos  
methyl

Dichlofluanid

Malathion

9.59min

Time (s)

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278

197

555

168

290

167

173

278


197

555

168



290

Efficiency



20.1m x 0.18mm x 0.18µm Rxi-5ms

Hydrogen 1.3 mL/min  
Oven at 22.8°C/min



5.75min

Results

Solve for

Efficiency

Speed

Translate

Custom

Run Time

11.63

9.74

min


Speed

1.19

x

Speed

Hydrogen 1.8 mL/min  
Oven at 27.2°C/min



4.85min

Efficiency

Methiocarb

20.1m x 0.18mm x 0.18µm Rxi-5ms

Helium at 1.0 mL/min  
Oven at 13.6°C/min

Fenthion

Pirimiphos  
methyl

Dichlofluanid

Malathion

Chlorpyrifos

9.59min

Time (s)

555

168

290

167

173

570

278

575

197

Speed

Hydrogen 1.8 mL/min  
Oven at 27.2°C/min

4.85min

Time (s)

280

282

168

290

167

286

173

288

278

290

197

292



# Simplify Method Development



- ✓ SOF, EOF and OHR
- ✓ *EZGC* Chromatogram Modeler
- ✓ *EZGC* Flow Calculator
- ✓ *EZGC* Method Translator

[www.restek.com/Landing-Pages/EZGC-R-Method-Development-Tool-Suite](http://www.restek.com/Landing-Pages/EZGC-R-Method-Development-Tool-Suite)

THANK YOU!!



[www.restek.com/Landing-Pages/EZGC-R-Method-Development-Tool-Suite](http://www.restek.com/Landing-Pages/EZGC-R-Method-Development-Tool-Suite)

**RESTEK**

Innovative Chromatography Solutions

[www.restek.com](http://www.restek.com)